ISMATEC®

IDEX
HEALTH & SCIENCE

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Operating Manual

Mode d'emploi

BVP Standard

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ISM 444

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Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
QP Q0 / Q1 / Q2 / Q3
Easy-Load / Standard / Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Deutsch

BVP Standard

with variable speed drive

ISM 444

For pump-heads:

Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
QP Q0 / Q1 / Q2 / Q3
Easy-Load / Standard / Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

English

BVP Standard

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ISM 444

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Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
QP Q0 / Q1 / Q2 / Q3
Easy-Load / Standard / Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Français

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Contents

Safety precautions	4–6
Warranty terms	7
Product	8
Operating panel	9
Rear panel	9
Mains voltage	10–11
Voltage setting	11
Changing the fuses	11
Starting the pump	12–13
Pumping against back-pressure	14
Overload protection	15
When the pump is not in use	16
Running-in period for tubing	17
Tubing life	17

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Contents

Analog interface	18–20
Accessories	
Foot switch	21
Spare cassettes	22
Click'n go cassettes	22
Cassettes with pressure lever	23
Pump-heads	24–25
Flow-rates per channel for multi channel pump-heads	26–28
Flow-rates per channel for single-channel pump-heads	29–31
Maintenance	32
Repairs	32
Spare-parts	32
Disposal	33
Technical specifications	34

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Caution: rotating pumphead creates a pinch and crush hazard

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Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pump-head). Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC® from date of delivery: 2 years

All other parts, excluding consumables: 1 year

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in case of:

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third-party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfac-ing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing

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Product

Package contents

- BVP drive
- Pump-head as ordered
- Power cord
- Operating instructions

➔ For operating the selected pump-head, please refer to the instruction manual separately supplied with the corresponding pump-head.

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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

Operating panel

- 1 Mains switch
- 2 Control light
- 3 Speed selector
999 = 240 rpm (max. speed)
adjustable in 0.1% steps
- 4 RUN/STOP
- 5 Rotation direction
- 6 MAX/SLOW key for drive speed
(ideal for fast filling or emptying the
system)

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⚠ Mains voltage

Mains voltage	Voltage setting	Fuse rating
220–240 V _{AC}	230 V 50/60Hz	2 x T1.25A/250V
110–120 V _{AC}	115 V 50/60Hz	2 x T2.50A/250V

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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
⚠ Mains voltage

Voltage setting 115V/230V and changing the fuses

➔ Switch the pump off,
pull out the mains plug.

1 Pull out the fuse-holder by opening it
with a small screw-driver (size 0).

2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse-
holder so that the required voltage
rating is facing the window of the
fuse-holder.

3 Insert 2 new fuses
➔ 230 V: 2 x T1.25A/250V
➔ 115 V: 2 x T2.50A/250V

⚠ Use always 2 slow-blow fuses of
the same type complying with the local
mains voltage. In case of a defective
fuse, always replace both.

4 Shut the fuse-holder. The voltage rating
is visible in the window.

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
Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel).
If necessary, change the rating as indicated on page 10–11.
- Mount the pump-head according to the instruction manual supplied with the pump-head.
- Insert the tubing into the peristaltic pump-head.
- If using an FMI piston pump-head, adjust the angle for the piston stroke volume.
- Connect the pump tubing to the system.
- Connect the pump to the mains.

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Starting the pump

- 1 Press in the mains switch, the control light is on.
- 2 Set the required speed on the 3-digit speed selector.
(999 = 240 rpm)
➔ The rotation speed can be adjusted while the pump is running.
- 3 Choose the rotation direction
➔ Never change the rotation direction while the pump is running.
- 4 The MAX/SLOW- button allows rapid filling (or emptying) of the tubes while the pump is running.
- 5 Start the pump by switching the RUN/STOP button.

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Pumping against pressure

The BVP *Standard* can be used for continuous duty at a differential pressure of max. 1.5 bar. For short-time operation higher differential pressures may be managed depending on the mounted pump-head and tubing i.d. (smaller i.d. is preferable).

Depending on the application conditions the FMI pump-heads allow operations at differential pressures of up to 7 bar.

In case of any doubts please contact your ISMATEC® agent.

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Overcurrent protector

The tubing pump BVP *Standard* features an overload protector.

When the overcurrent protector is activated, the pump is stopped and the control light flashes (see also the description of pin 12 of the analog interface on page 19).

Whenever this situation occurs, **the pump must be switched off immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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

When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

- a Tubing used with cassettes do not need to be removed. Releasing the cassette on one side is sufficient.
- b The single-channel pump-heads 360, 380 and Easy Load allow opening the tube-bed in order to release the tubing pressure.

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⚠ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.
Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Analog interface

Pin 1, (GND) Ground

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

Pin 5, speed IN

External speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Input impedance and input range can be selected via a dip-switch inside the pump (see page 20).

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
Analog interface

Pin 7, +36 V_{dc}

About +36 V_{dc} (unstabilised) are available the pump motor (max. load 1 A).

Pin 9, speed OUT

The default setting is 0-10 V_{dc}, proportionally to the motor speed 2.4-240 rpm. Alternatively a frequency range from 0 to 12 kHz is available. The input range can be selected via a dip-switch inside the pump (see page 20).

Pin 10, +5V_{dc}

This pin is connected to the stabilised +5V_{dc} voltage supply of the pump.

Pin 12, idle, overload

This pin is at +5V_{dc} while the pump runs. When idle or in overload, the voltage changes to 0V_{dc}.
➔ Switch off the mains voltage and let the device cool down for 2 minutes.

Pin 13, internal speed

If pin 13 is connected with pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at pin 5.

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
Pin 9 Ausgang/Output/Sortie
 A = 0–10 V_{DC} (Standard)
 B = 0–12 kHz

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Pins		Imp.	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0–5V	100kΩ	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0–10V	20kΩ	OFF	ON	OFF	OFF	OFF	ON
	0–20mA	250Ω	OFF	OFF	ON	OFF	ON	ON
	4–20mA	250Ω	OFF	OFF	OFF	ON	ON	OFF

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➔ **Make sure that the pump is disconnected from the mains supply.**

⚠ The instrument should only be opened by a qualified technician!
 Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Foot switch

Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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MS/CA Click'n'go Cassette
MS/CA Click'n'go cassette

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CA Click'n'go
Order No. IS 3710

CA pressure lever*
Order No. IS 0122

Adaptor for type CA IS 0123
(order 2 adaptors per cassette)

Click'n'go-Cassettes

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

Accessoires

Cassettes de rechange en POM-C

MS/CA Click'n'go
No de commande IS 3510

MS/CA levier de pression*
No de commande IS 0649

CA Click'n'go cassette
No de commande IS 3710

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CA pressure lever*

Order No. IS 3820

* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

Materials POM-C

Polyoxymethylene-Copolymer

- good chemical resistance to many organic solvents and strong alkaline chemicals
- is affected by strong acids and oxidizing substances
- UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

Materials PVDF

Polyvinylidene fluoride

- very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
- not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
- stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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Mounting

Each pump-head is supplied with an instruction manual.

For selecting the correct tubing, please refer to our detailed tubing documentation.

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


CA 12 / ISM 733

Interchangeable Pump-heads

These interchangeable pump-heads can also be ordered separately and mounted on the BVP drive.

➔ for the **piston pump-heads** the suffix „-W“ resp. „-WT“ means: with isolation (W) or temperature gland (WT) for crystallizing media.

0.0033–100 ml/min
MS/CA8-6
ISM 724
Erweiterungsblock:
Extension block:
Bloc d'extension:
ISM 185

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
⚠ Flow rates

In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®

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Maintenance

Provided the BVP drive is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective BVP drive to your ISMATEC® representative.

Please give information on defect and date of purchase.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Temperature +41 to +104°F

Rel. humidity max. 80%

– not condensing, at normal laboratory conditions

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EN 61326-1/EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products

619 Oak Street
Oak Harbor, WA 98277
United States

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Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH