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
Operating Manual

Mode d'emploi

Ecoline
VC-MS/CA8-6
VC-MS/CA4-12

Ecoline VC-280
Ecoline VC-380

Ecoline VC-360

Ecoline
VC-Easy-Load

ISMATEC®

IDEX
HEALTH & SCIENCE

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Deutsch

Tubing pumps with variable speed drive

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Contents

Safety precautions	4
Warranty terms	7
Product	8
Ordering Information	9
Mounting the pump-head	9
Rear panel	10
Mains voltage	10
Voltage setting	12
Changing the fuses	12
Operating panel	13
Starting the pump	14
Pumping against pressure	15
When the pump is not in use	15
Overload	16
Running-in period for tubing	17
Tubing life	17
Inserting the tubing	
Models VC-280/1 and VC-380/1	18
Changing the rotor	20
Model VC-360	22
Model VC-Easy-Load	23
Models VC-MS/CA8-6 and 4-12	24

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Contents

Flow rates		
Models	VC-280/1, VC-380/1 VC-360 VC-Easy-Load	25
Models	VC-MS/CA8-6 VC-MS/CA4-12	26
Analog interface		27
Settings of switch S1		29
Operation via foot switch		29
Accessories		30
Maintenance		32
Repairs		32
Spare-parts		32
Disposal		32
Technical specifications		
Models	VC-280/1, VC-380/1	33
Model	VC-360	34
Model	VC-Easy-Load	34
Model	VC-MS/CA8-6	35
Model	VC-MS/CA4-12	35

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when cassettes or tubing are inserted or changed.

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Caution: rotating pumphead creates a pinch and crush hazard

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Safety precautions

- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pump-head). Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed is completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 11 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

1 year from date of delivery

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

Our warranty becomes invalid in case of:

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third-party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing.

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
Product

Package contents

- Pump (model as ordered)
- Power cord
- Operating instructions

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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


Rear panel

- 1 Mains socket
- 2 Mains switch
on/off (start/stop)
- 3 Opening for screw-driver
- 4 Fuse-holder containing
2 fuses (change always both fuses)
- 5 Voltage selector (see page 12)
- 6 Analog input: Speed control
0–5 or 0–10 V / 0–20 V or
4–20mA)
Digital input (TTL level):
Start/Stop, rotation direction

⚠ Ventilation slits (a)

For sufficient ventilation keep the slits at the back and the base slits clear.

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Mains voltage

Mains voltage	Voltage setting	Fuse rating (slow-blow)
220–240 V _{AC}	230V 50/60Hz	2 x T1A/250V
110–120 V _{AC}	115V 50/60Hz	2 x T2A/250V

Socket/Power cord

Use exclusively the originally supplied power cord.

The socket must be earthed (protective conductor contact).

Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by opening it
with a small screw-driver (size 0).
 - 2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse-hol-
der so that the required voltage rating is
facing the window of the fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x T1A/250V
115 V_{AC}: 2 x T2A/250V
⚠ Use always 2 slow-blow fuses of
the same type complying with the local
mains voltage.
 - 4 Shut the fuse-holder. The voltage rating
is visible in the window.

⚠ **Before starting-up**
Check if the voltage setting visible in the
window of the fuse-holder complies with
your local mains voltage.

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Operating panel

- 1 Operating LED
on/off / Overload (blinking)
- 2 Rotation direction
left / stand-by / right
- 3 Speed setting
2-digit potentiometer
1–99%, 1% steps

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Ecoline VC-280/1
VC-380/1
VC-360
VC-Easy-Load
VC-MS/CA8-6
VC-MS/CA4-12

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Starting the pump

- Check if the voltage indicated on the identification label of the pump complies with your local mains voltage. Change if necessary as indicated on page 12.
- Insert the type of tubing required for the specific application. (For the tubing handling see pages 18–24).
- Connect the pump tubing to the system.
- Set the rotation direction.
- Connect the pump to the mains and switch it on.

➔ Rotor starts to revolve immediately unless the potentiometer is set to «00» or the switch for the rotation direction is in the stand-by position.

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Pumping against pressure

The Ecoline pumps and the tubing are destined for applications at differential pressure conditions not exceeding 1.5 bar.

When the pump is not in use

During downtimes we recommend to protect the tubing from unnecessary pressure and tension.

VC-MS/CA8-6 and VC-MS/CA4-12
Cassettes can easily be released by pushing up the fixing-tongue
VC-360 Open the tube-bed at the notch lever without removing the tubing
VC-280/1 and VC-380/1 Remove the tubing from the tube-bed
Easy-Load Open the tube bed

Caution

Please be aware that siphoning may occur at the moment the cassettes are released or the tube-bed is opened!

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Typical causes are:

excess backpressure, highly viscous medium, plugging or mechanical blockage, wrong tubing wall thickness, or low mains voltage.

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Overload condition

If the pump can't reach its set speed, the 'POWER' LED starts flashing. The pump continues to deliver but not at the desired flowrate. It can be operated for some time under a small overload. However, a reduced lifetime of motor and gearbox must be taken into account. As soon as the load decreases and the pump reaches its set speed, it stops flashing.

Should the load increase further instead, a motor current delimitter takes action and stops the device to prevent damage. The pump stands still but its 'POWER' LED keeps flashing. To reset, turn the device off. Make sure to remove the cause for the overload before operating the pump again.

Overload situations occur mostly in the combined presence of large tubing inner diameters and high speed.

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⚠ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Models VC-280/1, VC-380/1

- 1 ➔ Switch the pump off.

Loosen the screws of the protection plate and remove the plate

- 2 Move the upper tubing clamp to the middle position.

- 3 Insert the tubing and secure it with the upper tubing clamp.

- 4 Move the lower tubing clamp to the middle position.

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


Inserting the tubing

Models VC-280/1, VC-380/1

- 5 Turn the rotor counter-clockwise and simultaneously insert the tubing along the tube-bed until the tubing is squeezed by two rollers (tubing must not be twisted).
- 6 Secure the tubing with the lower tubing clamp.
- 7 Put back the protection cover and fasten it with the two screws.
- 8 Switch the pump on.
 - ➔ Run the pump for a few seconds; depending on the rotation direction open the upper or lower tubing clamp and restretch the tubing carefully.

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Interchangeable rotor

for VC-280/1 and VC-380/1

2 rollers for higher flow rates
3 rollers for lower pulsation

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

Changing the rotor


- 1 ➔ Switch the pump off.
- 2 Loosen the screw on the rotor using an Allen key (size 2 mm).
- 3 Remove the rotor by pulling it towards you.
➔ Be careful that the wedge and spacer remain on the shaft.
- 4 Insert the new rotor
➔ The screw on the rotor must be aligned with the wedge on the shaft.
- 5 Fasten the rotor screw tight using the Allen key.

⚠ For manipulating a pump-head also observe the safety precautions listed in the operating manual.

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Inserting the tubing

Model VC-360

- 1 ➔ Switch the pump off.
- 2 Press the notch lever upwards.
- 3 Open the tube-bed.
- 4 Push the upper tubing clamp upwards; insert the tubing.
- 5 Loop the tubing round the rotor; press the lower tubing clamp downwards; insert the tubing.
- 6 Close the tube-bed and press it to the right so that the notch lever locks into place.
- 7 ➔ Switch the pump on.
- 8 Run the pump for a few seconds; depending on the rotation direction open the upper or lower tubing clamp; restretch the tubing carefully.

Important

When the pump is idle open the tube-bed in order to protect the tubing from unnecessary load

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
Inserting the tubing

Model Easy-Load®

- 1 ➔ Switch the pump off.
- 2 Open the tube bed.
➔ Rotate the lever to the left.
- 3 Load the tubing. Center the tubing between the retainers.
- 4 Close the tube bed.
➔ Rotate the lever to the right.
- 5 To prevent the tubing from moving, fix it with the tubing retainers on the left and right side of the pump-head.
➔ Press the retainer slightly towards the head and simultaneously downwards.
The retainer locks in the position according to the tube size.
➔ For optimum tubing life, keep tubing straight where it enters and exits the pump-head.
- 6 ➔ Switch the pump on.

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Inserting the tubing

Model VC-MS/CA8-6 & 4-12

- ➔ Switch the pump off.
- Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously.
- Insert the tubing with one collar into the cassette.
- Let the tubing hang down (prevent it from being twisted).
- Insert the tubing with the second collar at the other end of the cassette.
- Reinsert the cassette into the roller-head.
- ➔ Switch the pump on.

Important

When the pump is idle, release all cassettes at the fixing-tongue.

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


⚠ Flow rates

In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water at 22°C and TygonST® tubing, without differential pressure

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

➔ Please refer to our detailed tubing documentation.

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Analog interface

➔ When using the analog interface, the switch for the rotation direction (see page 13, paragraph 2) must be in the neutral middle position.

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, shield (GND)

This pin couples the shielding of the cable the ground.

Pin 3, stop/run

The pump starts when connected to Pin 1 (GND).

Pin 4, direction

In open position the pump operates clockwise; when connected to pin 1 (GND) it pumps counter-clockwise.

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
Pin 2, shield
 Pin 3, stop/run
 Pin 4, direction
 Pin 6, speed intern

Analog-Eingang
Analog input/ Entrée analogique
 Pin 5, speed IN
 0–5 V_{DC} / 0–10 V_{DC} / 0–20 mA / 4–20 mA

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Pin 8

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Analog interface

➔ Switch for the rotation direction must be in the neutral middle position.

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA und 4–20mA).

Input impedance and input range can be selected via a DIP-switch inside the pump (see page 29).

Pin 6, speed intern

When connected to the ground, the drive speed can be controlled by the potentiometer integrated in the control panel.

Pin 7, +5V

This pin provides the internal supply voltage of 5 Volts for the supply of a potentiometer (RI = 330Ω).

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
Settings of switch S1

Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	100kΩ	OFF*	OFF*	OFF*	OFF*	OFF*	OFF*
	0-10V	20kΩ	OFF	ON	OFF	OFF	OFF	OFF
	0-20mA	240Ω	ON	OFF	OFF	OFF	OFF	OFF
	4-20mA	240Ω	ON	OFF	ON	OFF	OFF	OFF

* Default setting

➔ **Check that the pump is disconnected from the mains.**

For reaching the DIP-switch you need to unscrew the 4 screws on both sides of the pump. Then, lift the housing carefully (mind the cables) and remove it.

Operation via foot-switch

With DIP switch 4 the user can switch between two possibilities:

Dip-Switch 4 off = On/Off

Dip-Switch 4 on = on, as long as pressure is maintained on the foot-switch.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	100kΩ	OFF*	OFF*	OFF*	OFF*	OFF*	OFF*
	0-10V	20kΩ	OFF	ON	OFF	OFF	OFF	OFF
	0-20mA	240Ω	ON	OFF	OFF	OFF	OFF	OFF
	4-20mA	240Ω	ON	OFF	ON	OFF	OFF	OFF

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Accessories

Foot-switch

Order No. IS 3572
This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

Cleaning agent

Agent for cleaning and protecting the stainless steel pump housing, 1 PE bottle with screw cap, containing 2.5 dl
Order No. IS 3380

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MS/CA Click'n'go

MS/CA Anpresshebel
MS/CA pressure lever
MS/CA levier de pression

Zubehör

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MS/CA Click'n'go	
Bestell-Nr.	IS 3510
MS/CA Anpresshebel*	
Bestell-Nr.	IS 0649

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Accessories

Spare-cassettes in POM-C

MS/CA Click'n'go	
Order No.	IS 3510
MS/CA pressure lever*	
Order No.	IS 0649

Spare-cassettes in PVDF

MS/CA pressure lever*	
Order No.	IS 3629

* The cassettes with pressure lever are available on request. This type of cassette provides better results at elevated differential pressure conditions (>1 bar).

⚠ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

Accessoires

Cassettes de rechange en POM-C

MS/CA Click'n'go	
No de commande	IS 3510
MS/CA levier de pression*	
No de commande	IS 0649

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

Provided the Ecoline pump is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective Ecoline pump to your ISMATEC® representative.

Please give information on defect and date of purchase.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Differential pressure

¹⁾ Possible differential pressure with appropriate tubing material, tubing with small i.d. 's may enable higher pressures.

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Ecoline VC-280/1, VC-380/1

DC Motor

microprocessor controlled	
Speed range	3.5–350 rpm
Channels	1
Pump rollers	2 or 3
Differential pressure	max. 21.7 psi ¹⁾
Remote control	(see page 27–29)

Mains connection/Fuse rating

230 V (50 Hz)	2 x T1A/250V
115 V (60 Hz)	2 x T2A/250V
Power consumption	100 Watts
Protection rating	IP 30

Operating conditions

Temp.	+41° to 104°F
Rel. humidity	max. 80%
not condensing, at normal laboratory conditions	

Dimensions/Weight

DxWxH (in.)	10 1/8" x 6 3/4" x 5 3/8"
Weight (lb)	12 / 12

CE-compatibility proved according to:
EN 61326-1 / EN 61010-1

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Differential pressure

¹⁾ Possible differential pressure with appropriate tubing material, tubing with small i.d.'s may enable higher pressures.

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Ecoline VC-360, VC-Easy-Load

DC Motor

microprocessor controlled	
Speed range	3.5–350 rpm
Channels	1
Pump rollers	3
Differential pressure	max. 21.7 psi ¹⁾
Remote control	(see page 27–29)

Mains connection/Fuse rating

115 V (60 Hz)	2 x T2A/250V
230 V (50 Hz)	2 x T1A/250V
Power consumption	100 Watts
Protection rating	IP 30

Operating conditions

Temp.	+41° to 104°F
Rel. humidity	max. 80%
not condensing, at normal laboratory conditions	

Dimensions/Weight

DxWxH (in.)	9 3/8" x 6 6/8" x 5 1/2"
	11 3/16" x 6 6/8" x 5 1/2"
Weight (lb)	11 / 11.5

CE-compatibility proved according to:
EN 61326-1 / EN 61010-1

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Ecoline VC-MS/CA8-6 & 4-12

DC Motor

microprocessor controlled	
Speed range	3.5–350 rpm
Channels	8 4
Pump rollers	6 12
Differential pressure	max. 21.7 psi
Remote control	(see page 27–29)

Mains connection/Fuse rating

230 V (50 Hz)	2 x T1A/250V
115 V (60 Hz)	2 x T2A/250V
Power consumption	100 Watts
Protection rating	IP 30

Operating conditions

Temp.	+41° to 104°F
Rel. humidity	max. 80%
not condensing, at normal laboratory conditions	

Dimensions/Weight

DxWxH (in.)	12 5/16" x 6 3/4" x 5 1/2"
	11 1/16" x 6 3/4" x 5 1/2"
Weight (lb)	12 / 11.9

CE-compatibility proved according to:
EN 61326-1 / EN 61010-1

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IDEX Health & Science
Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

IDEX Health & Science
Ismatec Laboratoriumstechnik GmbH
Futtererstrasse 16
97877 Wertheim
Germany
Tel: +49 9377 9203-0
Fax: +49 9377 1388
fluidics.europe@idexcorp.com

www.idex-hs.com

Verlangen Sie den
ISMATEC®-Katalog.

Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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