Flowmaster FMT300

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3 Rollen

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Deutsch

**Tubing
pump**
Microprocessor
controlled

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⚠ Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump **must not** be operated outside the destined operating and environmental conditions.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes
- The pump must be switched off when pump-heads or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- The Pump-head consists of revolving parts. Hence, the pump may only be operated when the tube-bed is closed.

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Transporting the pump

The pump weighs 26 kg. (58 lb).

On principal, it should only be transported by using a transporting aid, e.g. a cart.

For lifting the pump ask for assistance.

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Caution: rotating pumphead creates a pinch and crush hazard.

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Safety precautions

- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend placing the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 32 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

2 years from date of delivery.

Tubes are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packing or in a packet of equal quality. The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- 1 Flowmaster drive
Order No. ISM 1020
incl. integrated power cord,
2 m long, with IEC 320 plug
(male)
- 1 power cord with an IEC 320 con-
nector (female plug) and a coun-
try specific mains plug
- 1 gasket (in plastic bag) for
analog interface (AD0069)
- 1 Operating manual

If ordered:
other accessories

Please check the package
and its contents for transport
damage. If you find any signs of
damage, please contact your local
ISMATEC® representative
immediately.

➔ Complaints can only be
accepted within 8 days from
receipt of the goods.

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230 V_{AC}: mind. 10A
115 V_{AC}: mind. 20A

Rear panel

1 Cover plate for the fuses and
the voltage selector
115/230 V

2 Power cord

3 Analog interface

Input for:

- speed control
0–5 V or 0–10 V, and
0–20 mA or 4–20 mA
- rotation direction
- Run/Stop
- Foot switch

Output for:

- speed 0–10 V_{DC} or 0–5 V_{DC}
- potential-free status out-
puts for busy, ok and error
- auxiliary supply
24 V_{DC}/0.5 V

⚠ For using the analog-interface under IP 65 conditions please refer to page 32.

⚠ Mains voltage

230 V_{AC} / 50 Hz or
115 V_{AC} / 60 Hz

Power consumption

Max. 500 W

Fuse rating of pump

230 V_{AC}: 2 x T4A/250V
1 x T1.6A/250V

115 V_{AC}: 2 x T8A/250V
1 x T1.6A/250V

⚠ **Socket/Power cord**

The **socket must be earthed**.
(protective conductor contact)

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
⚠ Before starting-up

Check if the voltage setting on the rear panel complies with your local mains voltage.

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Starting the pump

1. Open the tube-bed by lifting the lever and insert the tubing into the pump-head.

2. Close the pump-head and connect the pump tubing to the system.

3. Connect the pump to the mains and switch it on using the power supply switch. The following values are displayed: LED test and the firmware version.

⚠ Fuse protection of mains supply:
230 V_{AC}: min. 10A
115 V_{AC}: min. 20A

⚠ Inserting the tube

Please ensure the tubing is straight and kink free for a distance of at least 50 cm on either side of the pump-head inlet and outlet.

➔ A kinked tube can chafe at the tube clamps, which may lead to leakage. When reusing the same tube, ensure it is used in the same orientation as the last time it was used. It is preferable to leave the tube in the pump-head during idle times.

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Pumping with timer function

1 When the pump is idle, you can switch between pumping and timer function by using the MAX key.

➔ In the settings menu, the parameter MODE must be set to TOGL (see page 18).

2 Use the ▲▼ keys for entering the required time.

The time is indicated in the following time units:

- seconds 0.2" to 899"
- minutes 15.0' to 899'
- hours 15.0 h to 959 h
- days 40.0 d to 980 d

3 Start: press RUN/STOP key
The pump runs during the time indicated in the display.

4 Interrupt: press RUN/STOP key again (remaining time blinks in the display).

5 Continue: press RUN/STOP key again.

6 Reset: press escape key (remaining time is deleted).

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When the pump is idle

When the pump is idle, we recommend releasing the tubing from pressure. Opening the pump head is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

⚠ When reusing the same tube, ensure it is used in the same orientation as the last time it was used. It is preferable to leave the tube in the pump-head during idle times.
Opening the pump head is sufficient. This helps to protect the tubing.

Pumping against pressure

The Flowmaster pump can be used for continuous duty at a differential pressure of max. 2.0 bar. In addition the max. differential pressure is limited by the tubing. See the table below (page 13).

⚠ IMPORTANT

Please do not exceed the max. differential pressure (page 13). If the differential pressure is too high, the tube can chafe and be cut at the tube clamps.

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Max. operating pressure

Information supplied by the tubing manufacturers:

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Overload protector

The Flowmaster pump features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' (for 'Overload') and the pump is stopped.

OL I = Current

OL S = Speed

OL t = Temperature

Whenever this situation occurs, the pump must be switched off immediately.

Before the pump is restarted, it is most important to check the reason for the overload (e.g. too high differential pressure etc.).

➔ Only when the operating conditions have been changed accordingly may the pump be started again.

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
Additional functions for basic settings

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2 Digital LED display

3 ▲ Key

- Increase speed
- Increase dispense time²⁾

4 MAX key

- max. rotation speed
- Keep key pressed while the pump is running (ideal for filling or emptying the tubing)
- Toggle between Pump and Timer mode.²⁾

default¹⁾

Resets to its default value

5 RUN/STOP

- starts the pump at the preset speed
- stops the pump
- Pauses dispensing {remaining Time flashes}.²⁾

ok¹⁾

Confirms a data entry

6 Rotation direction

- Changes the rotation direction
- Changes the rotation direction when pump is at rest.²⁾
- Cancels a paused dispensing.²⁾

escape¹⁾

- Cancels an entry
- Quits the settings menu

7 ▼ Key

- Reduce speed
 - Reduce dispense time²⁾
- settings¹⁾**

Access to the settings menu

8 Indication of rotation direction

9 Operating mode

- setting: basic settings mode¹⁾
- pump is controlled via analog interface, see page 22–23

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Tubing life expectancy

If this mode is active, the pump warns you when the life expectancy of the tubing has elapsed. *)

➔ The display alternates between tUbE and RPM (Pump mode) or tUbE and dispensing duration (Timer mode) and all LED's are blinking.

*) The expectet tubing life time can be entered by the user, see OPT-LIFE on page 20.

Resetting the tubing life expectancy

Open the tube bed. The display shows OPEN and the LED's are blinking.

- Push the default key
- Confirm with the ok key
- Close the tube bed

➔ If the run/stop function is controlled via the analog interface, the RUN/STOP key functions as an emergency stop switch. After an emergency stop has occured, the display shows StOP and all LED's blink. To return to the normal operating mode, push the RUN/STOP key again.

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

Analog interface

Pin 1, +24 V_{DC}, 0.5 A

Auxiliary voltage for peripherals (max. current 0.5A)

Pin 2, status error

This pin is connected to 'status 0' (pin 9) if the pump is switched off or in the case of an error.

max.: AC 70 V / 1 A

DC 48 V / 1 A

Pin 3, output 0

Reference point for the 'busy' output, potential-free.

Pin 4, busy

See page 23.

Pin 5, direction

If this input is set to high (5, 12 or 24 V) relative to 'GND' (Pin 8) the pump rotates in the counter-clockwise direction. Configuration see page 17.

Pin 7, speed IN

For external speed control 0-5 V (1MΩ), 0-10*V (20kΩ) 0-20 mA (250Ω), 4-20 mA (250Ω) Configuration see page 17.

Pin 8, GND (ground)

Reference potential for all analog and digital inputs and the analog output. GND lies on earth potential.

Pin 9, status 0

Reference potential for both status outputs, potential-free.

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
Analog interface

Pin 10, status ok

If the pump is switched on and provided there is no error, this output is connected to 'status 0' (Pin 9).

Pin 11/12 (4) busy

This outputs are connected to 'output 0' (Pin 3) if the pump is running.

Pin 13, start

Two Run/Stop modes:

- Toggle mode: The pump starts or stops when the input changes from low to high. (Default value)
- Direct mode: The pump runs as long as this input is pulled up to high potential (5/12/24*V) relative to GND (Pin 8). Configuration see page 17.

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
Please also refer to our detailed tubing documentation.

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Please note

A kinked tube can chafe at the tube clamps, which may lead to leakage (see page 9).

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Tubing life

The tubing-life depends on the application and the tubing material used.

Example:

Chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

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We recommend using only the following tubing materials:

- Tygon Standard
- Tygon LFL
- PharMed
- Tygon Silicone

We will be pleased to help the user select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®.

The values indicated in the table above are only approximate values and determined as follows:

Litre/min, with water and PharMed® tubing, without differential pressure

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Foot switch

Order No. IS10279

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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Maintenance

Provided the pump is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective pump with detailed description of the failure to your ISMATEC® representative.

Please use the original ISMATEC® packing or a packet of equal quality.

Please give information on defect, date of purchase, serial-no., and model.

Spare-parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® representative:

- replacement parts
- parts lists
- wiring diagrams

Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Opening the casing

⚠ In order to maintain the best possible seal of the casing, the pump should only be opened for replacing the fuses or changing the DIP-switches. We recommend to proceed strictly as follows:

Removing the casing hood

1. Remove the 5 screws with a Phillips screw driver (size 2) on each side of the casing hood.
2. Spread **slightly** the casing hood at the bottom and remove it by carefully pulling it upwards.

Mounting the casing hood

1. Before replacing the casing hood we recommend you lubricate the sealing on the contact surface. In case that the two sealing ends have come apart, they must be stuck together again with Silicone sealing compound.
2. Spread the casing hood slightly at the bottom and push it carefully over the pump. Be careful that no cables are trapped.

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Opening the casing

3. Press the casing hood downwards and fasten the screws evenly from top to bottom.

When tightening the screws, be careful that the profile sealing is not trapped. If the sealing should be damaged, it must be replaced, otherwise the seal may no longer be tight.

Sealing the casing hood

When the casing hood fits tight, each single screw must be removed again and stopped with silicone sealing compound.

Then, the screws must be refastened again evenly from top to bottom.

➔ Remove the screw - fill the hole with silicone - fasten the screw again - remove the last bit of silicone with a rag. If necessary clean up with acetone.

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

1. Open the pump-head with the lever.

Remove the tubing.

2. Loosen the two screw plugs of the protection plate and pull off the plate.

- 3 Unscrew the roller head (screw driver size 6) and remove it from the motor shaft.

4. Remove the tube retainer spring by pressing it down on both sides and pulling it out.

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Disassembling the pump-head

5. Pull the button at the hinge. Remove the tube bed by pulling it upwards.

6. Take off the lever.

7. Detach the torsion spring.

➔ Cleaning the pump-head:

Without problems it can be rinsed it with warm tap water.

However, it is not designed for being cleaned neither in an autoclave nor in any kind of washing machine.

(The pump drive is IP 65-rated and as such dust-tight and washdown-proof).

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
Assembling the pump-head

1. Slip the torsion spring (b) over the ball bearings (a). Position the endings of the spring under the lateral pins (c).

Slip the sleeve of the lever over the ball bearings (a) and push it to the stop.

2. Insert the tube spring between the two distance bolts and the spring pin.

3. Push the roller head onto the motor shaft. The wedge (e) on the motor shaft must fit into the groove (d) in the center of the pump-head. Tighten the roller head with a screw driver (size 6).

4. Pull the button (f) at the hinge. Insert the tube bed from above. The three T-grooves (g) must slip over the three pins (c).

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

b c

Assembling the pump-head

5. Put the covering plate in place and fix it with the two screw plugs.

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Greasing the pump-head

When the pump-head is disassembled for cleaning, the following parts should be greased with an all-purpose grease:

a the shaft (inside and outside)

b the sleeve of the lever (inside and outside)

c the inner side of the lever contacting the rollers.

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

Plug gaskets

If the interface port is used under IP 65 conditions, the originally supplied gasket must be inserted into the male plug of the corresponding connecting cables.

- a Gasket for plug (15 pin, male) for analog interface (female) Spare part No. AD0069

⚠ Choosing the correct interface cable (sealed plug, e.g. moulded plug) according to the IP-rating used, is entirely the user's responsibility.

Changing the fuses

⚠ Switch the pump off, pull out the mains plug.

1. Remove the covering plate.

2. Remove the quarter-turn type cap of the fuse-holder. The fuse is located in the cap.

Fuses

230 V_{AC}: 2 x T4A/250V
1 x T1.6A/250V
115 V_{AC}: 2 x T8A/250V
1 x T1.6A/250V

➔ Follow the specifications on the pump.

3. Insert and fix the cap with the fuse. Refasten the covering plate with its four screws.

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

Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Motor type DC motor

Speed

5.0 – 500.0 rpm
digitally adjustable in steps
of 0.1 rpm

Pump rollers

3 rollers, 18/8 stainless steel

Differential pressure

Continuous max. 2.0 bar (30 psi)
Table with max. operating
pressures of tubing on page 12

Remote control

via analog interface
(see pages 22–23)

Mains connection

230 V_{AC} / 50 Hz or
115 V_{AC} / 60 Hz

Fuse rating control board

für 230 V: 2 x T4A/250V
1 x T1.6A/250V
für 115 V: 2 x T8A/250V
1 x T1.6A/250V

Power consumption

Max. 500 W

Protection rating

IP 65

Operating conditions

Temperature +5 to +40°C,
(41 to 104°F)
Rel. humidity max. 80%
– not condensing, at normal
environmental conditions

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Appendix Activating a foot switch

Two foot switch modes are available:

- FSt0 (Toggle mode)
The pump is started by pressing the foot switch. For stopping the pump, the foot switch must be pressed again.
- FSdr (direct mode)
The pump runs as long as pressure is maintained on the foot switch.

Entering the settings menu

- a) Keep the settings key pressed while switching the pump on **or**;
- b) press the settings key when the pump is already switched on and the tube-bed is open.

Activating the foot switch

- 1 - Select AnA with the ▲▼ keys
- ok (confirm your selection)
- 2 - Select PSEt with the ▲▼ keys
- ok (confirm your selection)
- 3 - Select r S. with the ▲▼ keys
- ok (confirm your selection)
- 4 - Select rEM with the ▲▼ keys
- ok (confirm your selection)

Now pin 13 of the analog interface is defined as RUN/STOP input.

- 5 Return to AnA by pressing the escape key several times.
- ok (confirm your selection)

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

Appendix Activating a foot switch

- 6 – Select LOGI with the ▲▼ keys
– ok (confirm your selection)
- 7 – Select 24 U with the ▲▼ keys
– ok (confirm your selection)
- 8 – Return to AnA by pressing several times the escape key
– ok (confirm your selection)
- 9 – Select FS with the ▲▼ keys
– ok (confirm your selection)
- 10 – Select FStO (Toggle mode) or FSdr (direct mode) with the ▲▼ keys (page 10)
– ok (confirm your selection)
- 11 Quit the settings menu by pressing the escape key several times.

Deactivating the foot switch

- Repeat position a) or b)
- Repeat steps 1 to 3
- ok (confirm your selection)
- Select LOCL with the ▲▼ keys
– ok (confirm your selection)

Quit the settings menu by pressing the escape key several times.

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IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products

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Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH