IPC / IPC-N

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High precision
Tubing Pump
with planetary drive

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IPC
4 channels ISM 930
8 channels ISM 931
12 channels ISM 932
16 channels ISM 933
24 channels ISM 934

IPC
4 Canaux ISM 930
8 Canaux ISM 931
12 Canaux ISM 932
16 Canaux ISM 933
24 Canaux ISM 934

IPC-N
4 Kanäle ISM 935
8 Kanäle ISM 936
12 Kanäle ISM 937
16 Kanäle ISM 938
24 Kanäle ISM 939

IPC-N
4 channels ISM 935
8 channels ISM 936
12 channels ISM 937
16 channels ISM 938
24 channels ISM 939

IPC-N
4 Canaux ISM 935
8 Canaux ISM 936
12 Canaux ISM 937
16 Canaux ISM 938
24 Canaux ISM 939

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® pumps are designed for pumping applications in laboratories and industry.

As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.

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Caution: rotating pumphead creates a pinch and crush hazard.

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Safety precautions

- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before all cassettes are inserted and fully snapped-in.
- Do not manipulate the pump-head or the cassettes before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 9 in this manual, must be used.
- The fuse-holder must not be short-circuited.

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Safety precautions

- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

2 years from date of delivery.

Tubes are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your Ismatec® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- IPC or IPC-N pump (type as ordered)
- inkl. Click'n'go cassettes (POM-C)
- 1 power cord with an IEC 320 connector (female plug) and a country specific mains plug
- Operating manual

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➔ Complaints can only be accepted within 8 days from receipt of the goods.

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

Window for voltage setting
If necessary, the voltage setting must be changed.

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Rear panel

- 1 RS232 IN
- 2 Analog interface
 - input for:
 - speed control
0–5 V or 0–10 V, and
0–20 mA or 4–20 mA
 - Run/Stop
 - Rotation direction
 - output for:
 - speed 0–10 V_{DC}
or 0–10.9 kHz
- 3 Mains switch (on/off)
- 4 Fuse-holder with voltage selector 115/230 V
- 5 Mains socket

⚠ **For further details on the analog interface please refer to page 33.**

⚠ Mains voltage

Mains voltage	Fuse rating (slow-blow)
220–240 V _{AC} 50/60 Hz	2 x T500mAL/250V
110–120 V _{AC} 50/60 Hz	2 x T500mAL/250V

⚠ **Before starting-up**
Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.
If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ **Socket/Power cord**
Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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Voltage setting 115V / 230V Changing the fuses

➔ Switch the pump off, pull out the mains plug.

- 1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).
- 2 Take out the voltage selector plate. Turn it and reinsert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
- 3 Insert 2 new fuses
➔ 230 V_{AC}: 2 x T500mAL/250V
➔ 115 V_{AC}: 2 x T500mAL/250V

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

- 4 Shut the fuse-holder. The voltage rating is visible in the window.

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
z.B. / e.g. / p.ex.

Starting the pump

- Check if **the voltage setting** visible in the window of the fuseholder complies with your local mains voltage (page 10).
- Prepare the cassettes (page 12)
- Insert the tubing (page 12)
- Insert the cassettes (page 12)
- Connect the pump tubing
- Connect the pump to the mains and switch it on with the power supply switch.

Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Pre-set tubing i.d.

➔ The pump now changes to the operating mode last used.

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


➔ Pump with planetary drive
The geared rollers rotate simultaneously with the shaft of the roller-head.

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Preparing the cassettes

- ➔ Switch the pump off.
- ➔ Remove the cassette by slightly pressing the release lever on both sides of the pump.
- ➔ The cassettes have slits on both sides. Insert a cassette adaptor into the slit on each side of the cassette.
Press the adaptor firmly against the bottom of the slit.

➔ The better you insert the adaptor the easier you can insert the cassettes loaded with the tubing.

Inserting the tubing

- ➔ Insert the tubing (2-stop type) with one stop into an adaptor. Stretch the tubing slightly so that it can be inserted with the second stop into the other adaptor.
➔ Ensure that the tube is not twisted.
- ➔ Reinsert all cassettes into the pump. Ensure that they snap in on both sides.
➔ The direction of the arrow engraved on the cassette indicates the flow direction providing the least pulsation.
- ➔ Switch the pump on and start it.

➔ When the pump is idle, we recommend you release the tubes from pressure by disengaging all cassettes on one side.

➔ Cassettes with pressure lever see page 44

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

* = Grundeinstellungen S. 15–17

Operating panel


- 1 Digital LED display
- 2 Increase value
- 3 – Reduce value
– settings*
- 4 **RUN/STOP**
 - Starts and stops the pump (the rotation direction is briefly indicated with a running dash)
 - Interrupts and continues a dispensing cycle (p. 22/25)
 - ➔ The function Start/Stop can also be triggered by an external device such as a foot switch or hand dispenser (accessories on page 42).
- 4a Combination with Pause•Time and RUN/STOP = stand-by function (pump runs at 1% of the max. flow rate)
- 5 **MAX/CAL**
 - Max. speed in the modes PUMP•rpm and PUMP%•Flow rate (ideal for filling or emptying the tubes; keep the MAX key pressed down while the pump is running)
 - Calibrating key for the flow rate or dispensing volume (pages 19 and 23)
 - ok*
- 6 **Rotation direction**
 - Changes the rotation direction (the rotation direction is briefly indicated in the display with a running dash)
 - Stops an interrupted dispensing cycle for good (pages 22 and 25)
 - escape*

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Operating panel

7 MODE

- Changes between the operating modes
- default*

8 LEDs for active operating mode

- **PUMP% Flow rate**
Pumping by flow rate in $\mu\text{l}/\text{min}$ and ml/min , or
– in % of max. speed
- **DISP Volume**
Dispensing by volume
in μl , ml or l
- **DISP Time**
Dispensing by time
0.1s to 999h
- **PAUSE Time**
– Pause time for
intermittent dispensing
0.1s to 999 h
– Stand-by mode

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Selecting the basic settings


- 1 Switch the pump off.
- 2 Keep the **settings** key pressed and switch the pump on.
- 3 Select the basic settings by using the **▲▼** keys
 - Tube
 - Init (resets all basic settings to the default values)
 - Address (RS232, pump address)
 - Digits (display)
 - Titration
 - FS (foot switch)
 - Back-steps (roller back-steps)
 - Cycles (Number of dispensing steps)
- 4 **ok**
= confirm your selection
- 5 Select the value by using the **▲▼** key, or
 - press **escape** for cancelling wrong input
 - certain values can be reset to the **default** value by using the default key (for setting values see the pages hereafter)
- 6 **ok**
= confirm entered value
- 7 **escape**
Quit the menu basic settings

➔ For selecting the basic settings refer to the following pages.

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Basic settings

Tube

- For selecting the tubing i.d. when pumping by flow rate (1.02 mm = default value)
- Percent (% of max. rotation speed)

Initializing

Pressing the ok key resets the following basic settings in the currently used program to the default values:

Mode	PUMP% Flow rate
Speed	30 % of max. speed
Direction	clockwise
Dispensing time	4.5 sec.
Pause time	2.0 sec.
Dispensing by volume	
IPC	27 roller steps
IPC-N	6 roller steps
Dispensing cycles	0
Roller back-steps	0
Foot switch	toggle
Digits	3
Titration	off (pump)

➔ Initializing does not affect the pre-set tubing i.d.

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Foot switch

or hand dispenser

- Toggle (default value)
 - Press on foot switch
= pump starts running
 - Press foot switch again
= pump stops
- Direct
the pump runs as long as the foot switch is pressed

Roller back-steps

Number of roller back-steps (for drop-free dispensing in mode DISP Volume)

1..100 (0 = default value)

Cycles

Number of dispensing steps by time or volume (pages 22/25), 0..9999 (0 = default value, infinity)

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⚠ First enter the tube i.d. (pages 15, 16)

Default setting = 1.02 mm

- 1 Change mode to PUMP% Flow rate
- 2 Enter the required flow rate with ▲▼ keys (displayed in µl/min and ml/min)
Examples of displayed values:
90.0µ = 90.0 µl/min
0.100 = 0.1 ml/min
9.00 = 9.0 ml/min
40.0 = 40.0 ml/min
- 3 Start pump with RUN/STOP
➔ Flow rate can also be adjusted while the pump is running.
➔ For an accurate flow rate we recommend you calibrate the pump (see p. 19).
- 4 Maintain pressure on MAX-key while the pump is running.
➔ Pump is running at maximum speed (ideal for fast filling or emptying the system).

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
Calibrating the flow rate

- 1 Change mode to PUMP% Flow rate (see pumping by flow rate, page 18)
Enter the required flow rate with the ▲▼ keys
- 2 Change mode to DISP Time
– Enter 60 seconds by using the ▲▼ keys
– Start pump with RUN/STOP
➔ Pump stops automatically after 60 seconds
➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:
- 3 Return to mode PUMP% Flow rate
- 4 Press the MAX/CAL button (displayed value blinks)
Enter the weighed or measured value with the ▲▼ keys
- 5 Confirm with the MAX/CAL-key (the flow rate setting returns automatically to the initially entered set point)
- 6 Start with the RUN/STOP-key
➔ Depending on the application conditions repeated calibration may be necessary.

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Pumping by drive speed


⚠ First select % speed (pages 15 and 16, TUBE)

- 1 Change mode to PUMP% Flow rate
- 2 Enter the required speed with the ▲▼ keys
 - 1.0 – 100.0% of max. speed (IPC = 45 rpm, IPC-N = 11.25 rpm)
 - adjustable in steps of 0.1%
- 3 Start pump with RUN/STOP
 - running dash indicates the rotation direction
 - ➔ The rotation speed can also be adjusted while the pump is running.
- 4 Maintain pressure on the MAX/CAL key while the pump is running.
 - ➔ Pump runs at maximum speed (ideal for fast filling or emptying the system).

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Stand-by mode

The purpose of this mode is to prevent the tubing system from **drying out** (e.g. over night).

The pump maintains in operation and runs at a max. flow rate of 1%.

- 1 Change the mode to PAUSE•Time
- 2 Press the RUN/STOP key
 - running dash indicates the rotation direction
 - running dot indicates the stand-by mode
 - The pump runs at a max. rotation speed of 1%
- 3 Press the RUN/STOP key again. The pump returns from the stand-by mode to Pause.

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Dispensing by volume

⚠ First enter the tube i.d. (pages 15, 16)

- Change the mode to DISP Volume
- Use the ▲▼ keys for entering the required dispensing volume. Three display readings are available:
 Microlitres: e.g. 17.0 µl = **17.0µ**
 Millilitres: e.g. 17.0 ml = **17.0**
 Litres: e.g. 17.0 l = **17.0L**
 (Display in litres from 999 or 9999 ml, resp.)
 ➔ For dispensing accurately we recommend you calibrate the pump (see p. 23).
- Start pump with RUN/STOP The pump stops after dispensing the set volume.
 ➔ The dispensing process can be interrupted by pressing RUN/STOP.
 a) pressing RUN/STOP again continues the dispensing process, or
 b) pressing the escape key aborts the interrupted dispensing process
 ➔ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process.
 ➔ Changing to the MODE PUMP% Flow rate during the dispensing process allows the user to retrieve and adjust the following values with the ▲▼ keys:
 – flow rate (µl or ml/min) or
 – rotation speed (in %)

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

Calibrating the volume

- 1 Change the MODE to DISP Volume
 ➔ Enter the dispensing volume with the ▲▼ keys
- 2 Start with RUN/STOP
 ➔ Pump stops automatically
 ➔ Measure the dispensed liquid by volume or weight
- 3 Stay in the same mode (DISP Volume) and press the MAX/CAL key (display blinks)
- 4 Enter the effectively dispensed volume with the ▲▼ keys
- 5 Save setting with the MAX/CAL key
 ➔ Pump is calibrated and returns automatically to the initially entered set point
- 6 Start pump with RUN/STOP
 ➔ Depending on the application conditions repeated calibration may be necessary.

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➔ At the same time, the volume is also reset to the default value.

Default calibration (volume)

- 1 Change the mode to DISP Volume.
- 2 Press the MAX/CAL key (displayed value blinks).
- 3 Press the default key (default value blinks).
- 4 Confirm with the MAX/CAL key.

➔ At the same time, the flow rate is also reset to the default value.

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Dispensing by time

The dispensing time can be entered from 0.1s to 999h.

- Change the mode to DISP Time
- Enter the required dispensing time with the ▲▼ keys
 ➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.
 Seconds: 0.1–899.9 (in steps of 0.1s)
 Minutes: 15'–899' (in 1min steps)
 Hours: 15h–999h (in 1h steps)
- Start pump with RUN/STOP. The pump stops after the dispensing time has elapsed.
 ➔ The dispensing process can be interrupted by pressing RUN/STOP
 a) pressing RUN/STOP again continues the dispensing process, or
 b) pressing the escape key aborts the interrupted dispensing process

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process.

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




Dispensing by volume within a pre-set time

- Change the mode to DISP Volume
Enter the required volume with the ▲▼ keys
 - Change mode to DISP Time
Enter the required time with the ▲▼ keys (0.1s – 999h)
 - Push the MODE key
➔ Both LEDs light up simultaneously
 - Start pump with RUN/STOP
➔ The pump dispenses the required volume within the pre-set time.
- a **Changes / Corrections**
Press the escape key
➔ For changing the dispensing time or dispensing volume, MODE returns directly to DISP Volume
➔ Repeat steps 1 to 3

b **Error message**
LEDs are blinking alternately

9999 = volume too large
time too short
1111 = volume too small
time too long

Despite volume and time adjustments, the LEDs keep on blinking:

➔ Change the tubing i.d. in order to obtain other flow rate.

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
Please note

For entering the number of dispensing
cycles see page 29.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pause

⚠ First enter the tube i.d. (pages 15, 16)

- 1 Change to PAUSE Time
Enter the pause time
(0.1s – 999h)
- 2 Press the MODE key once
➔ Both LEDs for PAUSE and
DISP Volume light up
simultaneously.
Enter the volume with the key.
- 3 Start the pump with the
RUN/STOP key.
➔ Pump dispenses at intervals
➔ Active mode blinks
➔ The dispensing or pause
cycle can be interrupted by pres-
sing RUN/STOP
a) pressing RUN/STOP again
continues the dispensing or pau-
se cycle, or
b) pressing the escape key
aborts the interrupted
dispensing process

➔ The dispensing volume or pause
time can be changed even
during the dispensing process. The
newly entered set point
appears shortly on the display.

- ➔ Changing to the MODE PUMP%
Flow rate during the
dispensing process allows the user
to retrieve and adjust the
following values with the ▲▼ keys:
- flow rate (µl/min or ml/min) or
 - rotation speed (in %)

If necessary, calibrate the volume as
indicated on page 23.

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

- 1 Change to PAUSE Time
Enter the pause time (0.1s – 999h) with ▲▼ keys
- 2 Press the MODE key twice
➔ Both LEDs for PAUSE and DISP Time light up simultaneously
Enter the dispensing time with the ▲▼ keys
- 3 Start pump with RUN/STOP
➔ Pump dispenses at intervals
➔ Active mode blinks
➔ The dispensing or pause cycle can be interrupted by pressing RUN/STOP
 - a) pressing RUN/STOP again continues the dispensing or pause cycle, or
 - b) pressing the escape key aborts the interrupted dispensing process

➔ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

➔ Changing to the MODE PUMP% Flow rate during the dispensing process allows the user to retrieve and adjust the following values with the ▲▼ keys:

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
Please note

For entering the number of dispensing cycles see page 29.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume). See basic settings on page 17.

➔ During the pauses the remaining number of dispensing cycles is indicated on the display.

➔ During long pauses the remaining pause time can be retrieved by pressing the MAX/CAL key.

➔ The dispensing or pause cycle can be interrupted by pressing RUN/STOP.

a) pressing RUN/STOP again continues the process, or

b) pressing the escape key aborts the interrupted dispensing process

➔ The number of dispensing cycles remains stored until changes are made in the basic settings (page 17).

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Pumping against pressure

The IPC pump can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures)

In case of any doubts please contact your ISMATEC® agent.

When the pump is idle

When the pump is idle, we recommend you release the tubing from pressure. Releasing the cassette on one side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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
Please note

Please also refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period.

If constant and repeatable flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated our detailed tubing documentation (without obligation).

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The image shows a black rectangular box with the white text 'OL' inside, representing an error code on a device's display.

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Overcurrent protector

The IPC tubing pump features an overload protector. When an overload condition occurs, it is indicated on the display by the letters 'OL' and the pump is stopped.

Whenever this situation occurs, the pump must be switched off immediately.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again)

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure, etc.).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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Analog interface

Pin 1, GND (ground)

Connected to the protective earth conductor. Reference potential for all other pins.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

- In normal operation (pin 2 open) pin 3 operates the foot-switch (see basic settings, page 17)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Alternatives by means of DIP switches inside the pump (see page 35)

Pin 7, +30V_{DC}

About +30 V_{DC} are available (max. current 1A).

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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- Pin 2, remote
- Pin 3, start
- Pin 4, direction
- Pin 13, internal speed

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
Analog interface

Pin 9, speed OUT

The default setting is 0–10 V_{DC}, proportionally to the motor speed of 45.0 rpm (IPC) or 11.25 rpm (IPC-N). Alternatively a frequency range from 0–10.9 kHz is available. Alternatives with sliding switch S2 inside the pump (see page 35).

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.1 A)

Pin 11, roller count

One impulse is available for each roller step.

Pin 13, internal speed

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated (Normal operation, i.e. pin 2 is open)
Pin 13 serves as **auto-start function**.
If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated (Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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Serial interface

RS232 IN (female)
A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1–8). It is factory-set at 1 and can be changed with @ (also refer to the basic settings described on page 16). The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterisk *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Programming Software

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.
<http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website:
<http://www.ismatec.com>

Connecting cable

For connecting an IPC pump to a PC via the 9-pin RS232 connector.
Order No. AG0013

Foot-switch

Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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
Please refer to our detailed tubing documentation or our website.

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Spare-cassettes

Spare-cassettes in POM-C

CA Click'n'go IS 3710
CA pressure lever ¹⁾ IS 0122
Adaptor POM-C IS 0123
(order 2 adaptors per cassette)

Spare-cassettes in PVDF

CA pressure lever ¹⁾ IS 3820
Adaptor PVDF IS 3861
(order 2 adaptors per cassette)

¹⁾ The cassettes with pressure lever are available on request. This type of cassette provide better results at elevated differential pressure conditions (>1 bar). Preparing these cassettes see page 44.

⚠ Click'n'go-Cassettes

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

Materials POM-C

Polyoxymethylene-Copolymer
– good chemical resistance to many organic solvents and strong alkaline chemicals
– is affected by strong acids and oxidizing substances
– UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

Materials PVDF

Polyvinylidene fluoride
– very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
– not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
– stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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

Preparing the cassettes with pressure lever

- 1 ➔ Switch the pump off
- 2 Push the pressure lever down (it must be relaxed)
Insert the tubing with 2 colour-coded collars
- 3 Reinsert the cassette(s) onto the rollers. Ensure that it snaps in on both sides.
- 4 ➔ Switch the pump on and start it.
- 5 Lift the pressure lever step by step until the medium is drawn in and a continuous flow is achieved (no further step than necessary!).

Note

- Under certain circumstances (e.g. new tubing, elevated differential pressure, etc.) we recommend the user to lift the lever by just one more step.
- Never over-tighten the pressure lever.
- The setting of the pressure lever determines the occlusion.
- The optimum occlusion depends on tubing material, inner diameter, wall thickness, age of tubing, and application conditions.
- As a side-effect the occlusion setting also affects the flow rate this, however, should not encourage the user to use the pressure lever for flow rate adjustments because long-term behaviour is not predictable.
- If, however, the occlusion setting is used to adjust the flow rate, the channel-to-channel flow rate conformity may lead to good results in the beginning of the operation, yet the individual flow rates may increasingly diverge at a later stage.
- For an optimum long-term flow rate stability and channel-to-channel conformity we suggest using the Click'n'go cassettes the flow rate.

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Flow rates

In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure.

We will be pleased to help the user select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from Ismatec®

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

Provided the IPC is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective IPC pump to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., and model.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

Technical specifications

Drive

Motor type DC motor

Speed

0.4 – 45.0 rpm (IPC)
0.11 – 11.25 rpm (IPC-N)
adjustable in steps of 0.1% or as flow rate

Pump rollers

8 rollers, 18/8 stainless steel
actively driven

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Please note

We also recommend you observe our **Warranty Terms as well as our Terms and Conditions of Sale.**

In case of any queries, please contact your local ISMATEC® representative.

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Technical specifications

Differential pressure

max. 1.0 bar (100 psi) (For short term use, higher pressures are possible provided that tubing with a small i.d. and/or cassettes with pressure levers are used)

Remote control

via RS232 and analog interface (pages 33–37)

Mains connection

115/230 V_{AC} / 50/60 Hz

Fuse rating

2 x T500mAL/250V (230 V_{AC})
2 x T500mAL/250V (115 V_{AC})

Power consumption

Max. 30 W

Protection rating

IP 30

Operating conditions

Temperature +41 to 104°F
Rel. humidity max. 80%
– not condensing, at normal laboratory conditions

Dimensions/Weight

Channels	DxWxH (inch)	Weight
4	7.1x5.7/6.9x5.1	10.1 lb
8	8.7x5.7/6.9x5.1	11.2 lb
12	10.2x5.7/6.9x5.1	12.8 lb
16	11.8x5.7/6.9x5.1	14.3 lb
24	15x5.7/6.9x5.1	17.4 lb

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conformément à:
EN 61326-1 / EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

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Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH