## MCP-CPF *Process*

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QP.Q0  
QP.Q1  
QP.Q2  
QP.Q3

Deutsch

**Valveless  
Rotary Piston  
Dispensing Pump**  
Microprocessor  
controlled

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## Piston pump- heads require thorough flushing after use

We recommend flushing the pump-head after **every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

**Do not pump fluids containing suspended solids that are abrasive.**

Please observe our note regarding filtration in the chapters Safety Precautions and Maintenance.

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## Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® ceramic piston pumps are designed for pumping and dispensing applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes.
- The circuit between mains supply and pump has to be earthed.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- The pump must be switched off before taking off the pump-head.
- Do not open or remove the housing while the pump is operating.

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### Media containing particles

**Do not pump media containing particles that exceed a diameter of 0.8 mm.**

**Piston stroke or stroke volume must be set large enough with respect to the particle size!**

**Piston fractures are not covered by our warranty.**

### Corrosive liquids

Corrosive liquids may eventually produce leak paths around the sealing surfaces of the pump-head.

The pump should be inspected for leaks on a regular basis.

### Dry running

Before starting to run the pump we advise you fill the pump-head with liquid.

Ensure that for priming the pump does not run dry for more than 15 sec.

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## Safety precautions

- Only new fuses, according to the specifications stated on page 56 in this manual, must be used.
- The fuse holder must not be short-circuited.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Pump drive: 2 years

### Pump-head:

We advise you read the enclosed operating manual of the pump-head manufacturer.

**Piston fractures are not covered by our warranty (see our instructions for flushing the pump-head on page 3).**

In case of any queries, please contact your nearest ISMATEC® representative.

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## Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

### Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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## Product

### Contents of the package

- MCP-CPF *Process* drive  
Magnet included  
Order No. ISM 919  
incl. integrated power cord,  
2 m long, with IEC 320 plug  
(male)
- 1 power cord with an  
IEC 320 connector (female) and  
a country specific mains plug
- 2 gaskets (in a plastic bag)  
for the analog (AD0069) and  
RS232-IN (AD0071) interfaces
- Operating manual

If ordered:

- Pump-head
- Software
- other accessories

**Please check the package**  
and its contents for transport  
damage. If you find any signs of  
damage, please contact your local  
ISIMATEC® representative  
immediately.

- ➔ Complaints can only be accepted  
within 8 days from receipt of the  
goods.

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
### ⚠ Please note

For using the analog or RS232 interface under IP-65 conditions please refer to page 56.

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## Rear panel

- 1 Mains switch (on/off)
- 2 Power cord
- 3 RS232 IN (female)
- 4 RS232 OUT (male)
- 5 Analog interface
  - inputs for:
    - speed control  
0 – 5 V or 0 – 10 V, and  
0 – 20 mA or 4 – 20mA
    - Rotation direction
    - Run/Stop
    - Foot switch
    - 2 digital inputs (TTL-level)
    - outputs for:
      - speed 0 – 10 V<sub>DC</sub>  
or 0 – 7.2 kHz
      - 2 universal outputs

### ⚠ Mains voltage

100 - 240 V<sub>AC</sub>  
50 – 60 Hz

no adjustments necessary

### Power consumption

Max. 100 W

### Fuse rating on control board\*

F1.6A/250V  
T4A/250V

### Fuse rating on power supply\*

T3.15A/250V

\* See also page 57

### ⚠ Socket/Power cord

The socket must be earthed (protective conductor contact).

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
## Starting the pump

- Connect the pump to the mains and switch it on with the power supply switch.
- Enter the ID Codes of the mounted pump-heads in the program memories currently used (see pages 15 – 16).
- Adjust the angle for the piston stroke volume (see page 20).
- Connect the pump tubing to the system.

## Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Identification Code (ID-Code) of the pump-head entered in the program last used (e.g. rH00)
- 4 Reverse piston strokes  
➔ are only displayed if at least 1 or more back-steps are entered
- 5 Program last used
- 6 Displays the operating mode last used

**⚠ Before starting the pump for the first time or after connecting a new pump-head,** the ID-Code must be entered in each program used with this pump-head (see basic settings, page 15 – 16).

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3 Increase value

4 **RUN/STOP**

- Starts and stops the pump
- Interrupts or continues a dispensing cycle (page 27).
- ➔ The Start/Stop function can also be triggered by an external device such as a foot switch.

4a Combination with Pause•Time and RUN/STOP = stand-by function (pump runs at min. speed 10 rpm), ideal for avoiding the drying out of the pump-head and the tubing system (page 26).

5 **CAL**

- Calibrating key for the flow rate or dispensing volume (pages 24/28) \* ok

6 **MAX/MIN/RESET**

- Multi-function key for:
  - Max. speed<sup>1)</sup> (maintain pressure on this key while pump is running) ideal for filling or emptying the tubing system
  - min. speed (when pump is idle) turns slowly
  - Stops an interrupted dispensing cycle for good
  - Resets the accumulated volume in mode TOTAL to zero.
  - Resets a flow rate or volume calibration to the default value (page 29)  
\* default

<sup>1)</sup> functions in the operating modes PUMP rpm, PUMP Flow rate and TOTAL

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



## Operating panel

- 7 **Rotation direction**  
Changes the rotation direction  
→ Each time the rotation direction is changed, TOTAL is reset to »0«.  
\* escape
- 8 **MODE**  
Changes between the operating modes.
- 9 **LEDs for active operating mode**
- **STROKE div (divisions)**  
Setting the stroke volume  
1.0 – 450.0 RH-type pump-heads  
0.10 – 10.00 Q-type pump-heads
  - **PROGRAM**  
Program selection 1 – 4
  - **PUMP rpm**  
Pumping by speed  
10.0 – 1800 rpm
  - **PUMP Flow rate**  
Pumping by flow rate in µl/min, ml/min or l/min
  - **DISP Time**  
Dispensing by time  
0.10 s to 999 h
  - **DISP Volume**  
Dispensing by volume in µl, ml or litres
  - **PAUSE Time**  
for intermittent dispensing from 0.10 s up to 999 h
  - **TOTAL**  
Read-out of totally delivered volume in µl, ml or litres

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
## Selecting the basic settings

- 1 Select the program (p. 19)
- 2 Switch the pump off
- 3 Keep the settings key pressed and switch the pump on
- 4 Select the basic settings by using the ▲▼ keys

### Can be set separately for each program

- Cycles (Number of dispensing steps)
- Backstep (Number of piston stroke back-steps)
- Head = ID code (pump-head)
- FS (foot switch)
- Digits (display)
- Init (resets all basic settings to the default values)

### Valid for all 4 programs

- Init4: resets in all 4 programs the above mentioned basic settings to their default values
  - Address (RS232, pump address) and
  - Baud rate (RS232)
- 5 ok = confirm your selection
  - 6 Select the value by using the ▲▼ keys, or
    - press **escape** for cancelling wrong input
    - certain values can be reset to the **default** value by using the default key (for values, refer to the following pages)
  - 7 ok = confirm entered value
  - 8 escape = quit the menu basic settings

➔ For selecting the basic settings refer to the following pages.

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## Program specific basic settings

→ Individually adjustable for each  
of the 4 programs.

### Cycles

Number of dispensing steps by time  
or volume 0..9999  
(0 = Default value, infinity)  
(See intermittent dispensing  
pages 32 and 33)

### Backsteps

Number of piston stroke back-steps  
(for drip-free dispensing in the  
mode DISP Volume)  
from 1 .. 100 (0 = Default value)

### Head

Identification code (ID code) of the  
pump-head:  
RH head: rH0 or rH00, rH1  
Q head: q 0 or q 1, q 2, q 3  
→ In order to obtain correct  
pumping and dispensing values, the  
individual ID code of the mounted  
pump-head must be entered;  
especially when a pump is used  
for the first time and each time the  
pump-head is changed.

→ In each of the 4 programs a dif-  
ferent pump-head can be stored!

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
## Program specific basic settings

➔ **Individually** adjustable for each of the 4 programs.

### Foot switch

- Toggle (default value)
  - Step on foot switch = pump starts running
  - Press foot switch again = pump stops
- Direct  
the pump runs as long as the foot switch is pressed

### Digits (display)

Number of significant digits, 3 or 4 (3 = default value)

### Initializing

Pressing the ok key resets the following basic settings **in the currently used** program to the default values:

Mode:	PUMP rpm
Speed:	1000 rpm
Dispensing time:	4.50 sec.
Dispensing volume and flow rate:	depends on the pump-head
Pause time:	2.00 sec.
Number of dispensing cycles	0
piston stroke back-steps	0

➔ **The pump-head ID code entered in the currently used program remains stored.**

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## Selecting the program

When switching the pump on, it always selects the previously used program.

- 1 Change mode to PROGRAM by using the MODE key.
- 2 Press the CAL button. The currently set program starts to blink. Change to the required program with the  $\blacktriangle$  $\blacktriangledown$  keys.
- 3 Confirm with the CAL key.
- 4 The pump automatically returns to the last entered operating parameters of the selected program.

$\triangle$  **From now on, any changes carried out in the operating modes are automatically stored in the currently selected program as soon as the RUN/STOP key is pressed.**

## TOTAL volume

Change the mode to TOTAL. The totally delivered volume is displayed.

Four display readings are available:

Microlitre: e.g. 17.0  $\mu$ l = **17.0  $\mu$**   
Millilitre: e.g. 17 ml = **17.0**  
Litre: e.g. 17.0 l = **17.0L**  
Litre: e.g. 1700 l = **1700.**  
(Display in litres from 9999 ml)

For zero-setting, either

- press the reset key twice, or
- switch the pump off.

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STROKE.	o div	o PROGRAM
PUMP	o rpm	o Flow rate
DISP	o Time	o Volume
PAUSE	o Time	*TOTAL

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
**Setting the stroke volume**

Enter the stroke volume in the settings menu (see basic settings, page 15).

The flow rate is determined by the stroke volume and the rotation speed.

➔ Flow rate =  
 Stroke volume x speed

However, there is no recommendation of how to set the stroke volume and rotation speed in an „optimum“ way.

Basically the following applies:

- **Small stroke volume/ high speed:**
  - low pulsation
  - better resolution of the flow rate and dispensing volume
  - tendency to cavitation (insufficient flow-through of the liquid in the pump-head)
- **Large stroke volume/ low speed:**
  - for elevated viscosity
  - for differential pressures
  - for reduced resolution of the flow rate and dispensing volume
  - for easily degassing media, or such with a low steam pressure

Important: The pump always dispenses integer multiples of the stroke volume. If the stroke volume is adjusted to 10 µl and the targeted volume is 75 µl, the pump permits only 70 or 80 µl when selecting the dispensing volume with the ▲▼ keys.

➔ Adjust stroke volume mechanically to 7.5 µl, i.e. for a 25 µl pump-head from 180 units to 135 units.

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➔ Enter the required speed in the mode PUMP rpm (page 25).

- 1 Change the mode to STROKE div.
- 2 Press the CAL key. A blinking stroke volume appears in the LED display.
- 3 Enter the required stroke volume with the ▲▼ keys.
- 4 Confirm with the CAL key.
- 5 Adjustments on the head:  
**RH pump-heads:** Turn the knurled adjustment ring until the correct value on the fixed and revolving flow control ring equals the digital LED read-out. (The fixing screw must not be loosened!) (p. 22).

**Q pump-heads:** Adjust the required value on the revolving flow control knob.

### For pump-heads RH and Q:

➔ For accurate dispensing you should not set the stroke volume below 4 % (Q pump-heads) and 10% (RH pump-heads) of the max. value.

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## RH pump-heads Zero point calibration

Before using the pump the first time, or whenever necessary, we recommend you calibrate the 0-point of the flow control ring on the RH pump-head.

- 1 Start the pump in the mode PUMP Flow rate with the RUN/STOP key and pump the liquid to be delivered.
- 2 While the pump is in operation, turn the knurled adjustment ring to the right until the pump stops delivering the liquid (the motor is still running).

Loosen the fixing screw of the revolving flow control ring (black) and turn it to the value '0'.

Tighten the fixing screw again.

- 3 Set the required stroke volume and configure the basic setting stroke volume accordingly (page 21).

➔ This calibration is only valid for the RH pump-head. The Q pump-head does not require a zero point calibration.

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

### Please note

For accurate dispensing we recommend you calibrate the pump (page 24).

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## Pumping by flow rate

- Change the mode to PUMP Flow rate.
- Enter the required flow rate with the ▲▼ keys. Depending on the pump-head the flow rate can be displayed as follows:
 

Micro litres:	e.g. <b>1.70μ</b> = 1.70 μl
Micro litres:	e.g. <b>17.0μ</b> = 17.0 μl
Micro litres:	e.g. <b>0.170</b> = 170 μl
Millilitres:	e.g. <b>1.70</b> = 1.7 ml
Millilitres:	e.g. <b>17.0</b> = 17 ml
Millilitres:	e.g. <b>170</b> = 170 ml
Millilitres:	e.g. <b>999</b> = 999 ml
Litres:	e.g. <b>1.70L</b> = 1.70 l
- Start the pump with RUN/STOP.
  - ➔ The flow rate can be adjusted while the pump is running.
  - ➔ During the pumping process the following values can be retrieved with the MODE key:
    - Rotation speed (PUMP rpm)
    - Totally delivered volume (TOTAL)
- Maintain pressure on the MAX/MIN key while the pump is **running**.
  - ➔ Pump is running at max. speed (ideal for fast filling or emptying the system).
- Maintain pressure on MAX/MIN-key when pump is **idle** (power switched on!).
  - ➔ pump is turning at 10 rpm (for mounting the pump-head).

⚠ Entering the the preset stroke volume and the ID code of the currently mounted pump-head allows the user to set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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## Calibrating the flow rate

➔ Check and, if necessary, enter the ID code of the mounted pump-head (page 16).

1 Change the mode to PUMP Flow rate. Enter the required flow rate with the ▲▼ keys.

2 Change the mode to DISP Time. Enter 60 seconds by using the ▲▼ keys. Start the pump with RUN/STOP.  
➔ The pump stops automatically after 60 sec.  
➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:

3 Return to the mode PUMP Flow rate.

4 Press the CAL button (displayed value blinks).  
➔ The resolution of the LED display changes automatically to 4 digits. Enter the weighed or measured value with the ▲▼ keys.

5 Confirm with the CAL key (the flow rate setting returns automatically to the initially entered set point).

6 Start with the RUN/STOP key.

➔ Depending on the application conditions repeated calibration may be necessary.

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
## Pumping by speed

- Change mode to PUMP rpm. The speed is digitally adjustable in steps of  
 $0.1 \text{ rpm}$ :  $10.0 - 999.9 \text{ rpm}$   
 $1 \text{ rpm}$ :  $1000 - 1800 \text{ rpm}$
- Enter the required speed with the ▲▼ keys. Start the pump with RUN/STOP.
  - ➔ The speed can also be adjusted while the pump is running.
  - ➔ During the pumping process the following values can be retrieved with the MODE-key:
    - Flow rate (ml/min)
    - Totally delivered volume (TOTAL)
- Maintain pressure on the MAX/MIN key while the pump is **running**.
  - ➔ Pump runs at maximum speed (ideal for fast filling or emptying the system).
- Hold the MAX/MIN-key when pump is **idle** (power switched on!).
  - ➔ Pump is turning at  $10 \text{ rpm}$  (for mounting the pump-head).

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## Stand-by mode

The purpose of this mode is to prevent the pump-head and the tubing system from **drying out** (e.g. over night).

The pump remains in operation and runs at a min. speed of 10 rpm.

- 1 Change the mode to PAUSE Time
- 2 Press the RUN/STOP key
  - running dash indicates the rotation direction
  - running dot indicates the stand-by mode
  - The pump runs at a speed of 10 rpm.
- 3 Press the RUN/STOP key again. The pump returns from the stand-by mode to Pause.

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- 1 Change the mode to DISP Volume
- 2 Use the ▲▼ keys for entering the required dispensing volume. The following display readings are available:

Microlitres:	e.g.	1.70 µ =	<b>1.70 µl</b>
Microlitres:	e.g.	17.0 µ =	<b>17.0 µl</b>
Microlitres:	e.g.	0.170 =	<b>170 µl</b>
Millilitres:	e.g.	1.70 =	<b>1.7 ml</b>
Millilitres:	e.g.	17.0 =	<b>17 ml</b>
Millilitres:	e.g.	999 =	<b>999 ml</b>
Millilitres:	e.g.	1.70L =	<b>1.7 l</b>
Litre:	e.g.	170L =	<b>170 l</b>
Litre:	e.g.	1700. =	<b>1700 l</b>

➔ For accurate dispensing we recommend you calibrate the pump (page 28).

- 3 Start the pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process. The newly entered set point appears shortly in the display.

➔ During the dispensing process the following values can be retrieved with the MODE key:

- Rotation speed (PUMP rpm)
- Flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the rotation speed or flow rate can be changed even during the dispensing process.

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1 Change the mode to DISP Volume.  
➔ Use the ▲▼ keys for entering the required dispensing volume.

2 Start the pump with RUN/STOP.  
➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.

➔ Measure the dispensed liquid by volume or weight.

3 Stay in the same mode (DISP Volume) and press the CAL-key (display blinks).  
➔ The resolution of the LED-Display changes automatically to 4 digits.

4 Enter the effectively dispensed volume with the ▲▼ keys.

5 Save the entered value with the CAL key.  
➔ The pump always dispenses integer multiples of the stroke volume. If the stroke volume is adjusted e.g. to 10 µl and the targeted volume is 75 µl, the pump permits only 70 or 80 µl when selecting the dispensing volume with the ▲▼ keys.

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➔ At the same time, the volume is also reset to the default value.

### Default calibration of volume

- 1 Change the mode to DISP Volume.
- 2 Press the CAL button (displayed value blinks).
- 3 Press the RESET button (displayed value blinks).
- 4 Confirm with the CAL button.

➔ At the same time, the flow rate is also reset to the default value.

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



## Dispensing by time

The dispensing time can be entered from 0.10 s – 999 h.

- 1 Change mode to DISP Time
- 2 Enter the required dispensing time with the ▲▼ keys  
 ➔ The display accelerates when holding ▲ or ▼ key.

Seconds: 0.10" – 99.9" (in 0.01 s steps)  
 0.1" – 899.9" (in 0.1 s steps)  
 Minutes: 15' – 899' (in 1 min steps)  
 Hours: 15 h – 999 h (in 1 h steps)

- 3 Start pump with RUN/STOP

- ➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect from the subsequent dispensing cycle.
- ➔ During the dispensing process the following values can be retrieved by pressing the MODE-key:
  - speed (PUMP rpm)
  - flow rate (PUMP Flow rate)
  - totally dispensed volume (TOTAL)
- ➔ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.
- ➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible.

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

## Dispensing by volume within a pre-set time

- 1 Change the mode to DISP Time. Enter the required time with the ▲▼ keys (0.10 s – 999 h).
  - 2 Change the mode to DISP Volume. Enter the required volume with the ▲▼ keys.
  - 3 Push the MODE button.  
➔ Both LEDs light up simultaneously.
  - 4 Start the pump with RUN/STOP.  
➔ The pump dispenses the required volume within the pre-set time.
- a Changes / Corrections**  
Press the Reset key  
➔ For changing the dispensing time or dispensing volume, the pump returns directly to the mode DISP Time  
➔ Repeat steps 1 to 3
- b Error message**  
LEDs are blinking alternately  
9999 = volume too large  
time too short  
1111 = volume too small  
time too long
- Despite volume and time adjustments, the LEDs keep on blinking:  
➔ Adjust the stroke volume or change the pump-head model in order to obtain another flow rate range.  
➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not repeatable.

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### Please note

For entering the number of dispensing  
cycles see page 34.

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## Intermittent dispensing (by volume)

### Intermittent dispensing by volume with a pause

1 Change to PAUSE Time  
Enter the pause time  
(0.10 s – 999 h)

2 Press the MODE key twice  
➔ Both LEDs for PAUSE and  
DISP Volume light up  
simultaneously. Enter the  
volume with the ▲▼ keys.

3 Start the pump with the  
RUN/STOP key  
➔ Pump dispenses at intervals  
➔ Active mode blinks  
➔ The dispensing or pause cycle  
can be interrupted by pressing  
RUN/STOP  
a) pressing RUN/STOP again  
continues the dispensing or  
pause cycle, or  
b) pressing the RESET key  
aborts the interrupted dispensing  
process.

➔ The dispensing volume or  
pause time can be changed even  
during the dispensing process. The  
newly entered set point appears  
shortly on the display.

➔ During the dispensing process the  
following values can be  
retrieved with the MODE-key:  
– speed (PUMP rpm)\*  
– flow rate (PUMP Flow rate)\*  
– totally dispensed volume  
(TOTAL)

\*The rotation speed and flow rate  
can be changed with the  
▲▼ keys.

If necessary, calibrate the volume as  
indicated on page 28.

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## Intermittent dispensing (by time)

### Intermittent dispensing by time with a pre-set pause time

- 1 Change to PAUSE Time  
Enter the pause time (0.10 s – 999 h) with ▲▼ keys
- 2 Press the MODE key  
→ Both LEDs for PAUSE and DISP Time light up simultaneously.  
Enter the dispensing time with the ▲▼ keys.
- 3 Start pump with RUN/STOP  
→ Pump dispenses at intervals  
→ Active mode blinks  
→ The dispensing or pause cycle can be interrupted by pressing RUN/STOP  
a) pressing RUN/STOP again continues the dispensing or pause cycle, or  
b) pressing the RESET key aborts the interrupted dispensing process.

→ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

→ During the dispensing process the following values can be retrieved with the MODE-key:

- rotation speed (PUMP rpm)\*
- flow rate (PUMP Flow rate)\*
- totally dispensed volume (TOTAL).

\*The rotation speed and flow rate can be changed with the ▲▼ keys.

If necessary, calibrate the volume as indicated on page 28.

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
settings

STROKE	o div	o PROGRAM
PUMP	o rpm	o Flow rate
DISP	*Time	*Volume
PAUSE	*Time	o TOTAL

STROKE	o div	o PROGRAM
PUMP	o rpm	o Flow rate
DISP	o Time	*Volume
PAUSE	*Time	o TOTAL

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## Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

After entering the dispensing cycles (basic settings page 16) the pump can be started with the Run/Stop key.

According to the pre-set pause and dispensing time or dispensing volume, the pump now dispenses at intervals and stops when the number of dispensing cycles has elapsed.

➔ During the pauses the remaining number of dispensing cycles is indicated on the display.

- During long pauses the remaining pause time can be retrieved by pressing the CAL key.
- The dispensing or pause cycle can be interrupted by pressing RUN/STOP.
- Pressing RUN/STOP again continues the process.
- Pressing the RESET key aborts the interrupted dispensing process.

➔ The number of dispensing cycles remains stored until changes are made in the basic settings (page 16).

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## Drip-free dispensing

Reverse piston strokes (back-steps) allow drip-free dispensing (1 – 100 strokes). We specially recommend this operating mode if corrosive media is pumped or in case that only few sample material is available.

- 1 Switch the pump off
- 2 Hold the settings key and switch the pump on again
- 3 Select bStP by using the ▲▼ keys (basic settings, page 16)  
ok = confirm selection
- 4 Enter the number of reverse piston strokes with the ▲▼ keys  
➔ in case of wrong input return to bStP with the escape key  
ok = confirm selection
- 5 escape = quit menu basic settings

### Reverse piston strokes are active in the following modes:

- a Dispensing by volume (see page 27)
- b Dispensing a volume within a pre-set time (see page 31)
- c Intermittent dispensing by volume (see page 32)

⚠ Reject the first dispensing step. Afterwards, the number of reverse piston strokes are compensated.

### Deleting the reverse piston strokes

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b 

STROKE	o div	o PROGRAM
PUMP	o rpm	o Flow rate
DISP	* time	* Volume
PAUSE	o time	o TOTAL

c 

STROKE	o div	o PROGRAM
PUMP	o rpm	o Flow rate
DISP	* time	o Volume
PAUSE	* time	o TOTAL

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### Pumping against pressure

Depending on the pump head the MCP-CPF *Process* can be used for continuous duty at a differential pressure of max. 6.9 bar (100 psi).

### When the pump is idle

We recommend flushing the pump-head **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

### Overcurrent protector

The drive MCP-CPF *Process* features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL S = Speed set-point cannot be reached

OL I = Current Overload

Whenever this situation occurs, the pump must be switched off immediately.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again)  
Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
## Analog interface

### Pin 1, GND (ground)

Connected to the protective earth conductor. Reference potential for all other inputs.

### Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

### Pin 3, start

- In normal operation (pin 2 open) pin 3 operates the foot-switch (see basic settings page 17)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

### Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

### Pin 5, speed IN

External speed control (0 – 5 V, 0 – 10 V, 0 – 20 mA, 4 – 20 mA) Input impedance and input range can be selected via a DIP switch inside the pump (see page 39)

### Pin 6, input 1

### Pin 8, input 2

Digital inputs (TTL-level). Active when connected to pin 1. Are only retrievable via RS232 interface; e.g. with a software like the ProgEdit.

### Pin 7, +26V<sub>DC</sub>

About +26 V<sub>DC</sub> are available (max. current 1A).

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## ⚠ Please note

Depending on the operating temperature, analog inputs and outputs may have a slight drift. If highest performance is required with respect to speed setting and speed constancy, we recommend to control the pump via the RS232 interface. For monitoring the speed via the output "speed OUT" we suggest using the "frequency mode".

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## Analog interface

### Pin 9, speed OUT

The default setting is 0 – 10 V<sub>DC</sub>, proportionally to the motor speed.

Alternatively a frequency range from 0 – 7.2 kHz is available. Alternatives with sliding switch S2 inside the pump (see page 39)

### Pin 10, +5 V<sub>DC</sub>

About +5 V<sub>DC</sub> are available. (max. current 0.5 A)

### Pin 13, speed intern

Depending on how the pump is operated, pin 13 has different functions:

### Analog interface not activated

(Normal operation, i.e. pin 2 is open) Pin 13 serves as **auto-start function**.

If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

### Analog interface activated

- (Pin 2 on GND)
- Pin 13 open:  
The rotation speed is adjusted via pin 5 (speed IN).
  - Pin 13 on GND:  
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

### Pin 14, output 1 (PWM)

### Pin 15, output 2 (PWM)

Digital outputs in connection with pin 7 (+26 V<sub>DC</sub>), implemented as **pulse-width-modulation** (open collector). Can only be activated via RS232 interface; e.g. with a software like the ProgEdit.

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This switch affects Pin 9, motor speed

Position A: 0 – 10 V<sub>DC</sub> (Standard)  
Position B: 0 – 7.2 kHz

➔ In order to maintain the optimal casing seal, the pump should not be opened without a particular reason. This is, however, necessary in order to gain access to the switches. For this purpose we recommend that you proceed exactly as indicated on pages 55 – 56.


➔ **Make sure that the pump is disconnected from the mains supply.**

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### Please note

**The instrument should only be opened by a qualified technician!**  
(see pages 55 – 56)

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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## Please note

When switching on the MCP-CPF *Process* drive, the serial interface will only respond after a delay of 3 seconds.

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## Serial interface

### RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

### Setting the pump address

Each command string begins with the pump address (1 – 8). It is factory-set at 1 and can be changed with @ (also refer to the basic settings described on page 18). The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one RS232 interface.

### Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterix \*. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

**Incorrect** command strings are answered by #.

**Numerical values** are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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## Cascading several pumps (daisy-chain)

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address (also see basic settings on page 18).**

➔ Alternatively, the addressing can also be accomplished via the serial interface.

For that, each pump must firstly be connected individually to a PC and the address must be allocated with @. Only then may the pumps be daisy-chained. Each pump requires a definite address in order to avoid interference of their responses.

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
**Each pump requires an individual address**  
(see page 40).

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
Bestell-Nr. AG0013

## Programming Software

### LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.

<http://www.ni.com/>

### LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website:

<http://www.ismatec.com>

### Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments

Order No. IM0015C168

### Adaptor 25 F / 9M

25 pins, female / 9 pins, male for IM0015C168 (1 adapter for each controlled instrument)

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## Programming Software

### Software ProgEdit

- Create your individual application program.
- Download the program via RS232 interface to your MCP-CPF *Process* pump drive.
- Now run your application PC-independently wherever needed.

#### System requirements:

3 1/2" floppy drive, Win95 or higher, Pentium or higher.

- ➔ **Free full version**  
[www.ismatec.com](http://www.ismatec.com)
- ➔ **Downloads**

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\* For entering the pump-head ID code via basic settings (page 15)

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\* Pumpenkopf

⚠ In order to reach the max. flow rate, the stroke volume must be set to the max. value (depending on the pump-head model).

➔ If, depending on the application, this stroke volume must be reduced, we recommend you **not to go below 10% (RH\*) or 4% (Q\*) of the maximum value.**

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\* tête de pompe

Pumpenkopf Pump-head Tête de pompe	Bestell-Nr Order No. No comm.	Kolben Piston	Zylindergehäuse Cylinder case Boîtier du cylindre	Zylinder Cylinder Cylindre	Mediubenetzte Teile Wetted parts Parties mouillées	Temperatur max.	Differenzdruck diff. pressure Pression différentielle max.	Ein- / Auslass Main flow ports Raccords d'entrée / sortie	Spülanschluss Gland flow ports Raccord pour rinçage
RH00.CKC-LF	FMI 009	Ceramic	Kynar Fluorocarbon (PVDF)	Ceramic	Ceramic PVDF	100° C/212° F	6.9 bar / 100 psi	Kynar UNF 1/4"-28 (female)	--
RH00.SKY-LF	FMI 010	SS316	Kynar Fluorocarbon (PVDF)	Carbon	SS316 PVDF / Carbon	60° C/140° F	6.9 bar / 100 psi	Kynar UNF 1/4"-28 (female)	--
RH0.CKC	FMI 005	Ceramic	Kynar Fluorocarbon (PVDF)	Ceramic	Ceramic PVDF	100° C/212° F	6.9 bar / 100 psi	for tubing with 6.4 mm (1/4") o.d.	--
RH1.CKC	FMI 007	Ceramic	Kynar Fluorocarbon (PVDF)	Ceramic	Ceramic PVDF	100° C/212° F	6.9 bar / 100 psi	for tubing with 6.4 mm (1/4") o.d.	--
QP.Q0.SSY	FMI 202	SS316	SS316	Carbon	SS316 / Carbon	60° C/140° F	6.9 bar / 100 psi	1/4 NPT (female)	--
QP.Q0.SKY	FMI 316	SS316	Kynar Fluorocarbon (PVDF)	Carbon	SS316 PVDF / Carbon	60° C/140° F	4.1 bar / 60 psi	for tubing up to 12.7 mm (1/2") i.d.	--
QP.Q1.CSC	FMI 205	Ceramic	SS316	Ceramic	SS316 PVDF / Carbon	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	--
QP.Q1.CSC-W	FMI 320	Ceramic	SS316	Ceramic	Ceramic / SS316	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	10-32 (female)
QP.Q1.CSC-WT	FMI 219	Ceramic	SS316	Ceramic	Ceramic / SS316	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	1/8 NPT (female)
QP.Q2.CSC	FMI 212	Ceramic	SS316	Ceramic	Ceramic / SS316	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	--
QP.Q2.CSC-W	FMI 321	Ceramic	SS316	Ceramic	Ceramic / SS316	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	10-32 (female)
QP.Q2.CSC-WT	FMI 218	Ceramic	SS316	Ceramic	Ceramic / SS316	177° C/350° F	6.9 bar / 100 psi	1/4 NPT (female)	1/8 NPT (female)
QP.Q3.CKC	FMI 217	Ceramic	Kynar Fluorocarbon (PVDF)	Ceramic	Ceramic PVDF	100° C/212° F	1.7 bar / 25 psi	for tubing up to 12.7 mm (1/2") i.d.	--

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= With isolation gland for crystallizing media

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= With isolation gland for crystallizing media and threads for the use of 2 cartridge heaters and a thermocouple

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**Low Flow Kit for Q-type pump-heads with stainless steel casing**

is recommended for flow rates below 50 ml/min (e.g. for chromatography)

R479 FMI 056

**Port/tubing adaptors**

Allow tubing arrangements other than 6.4 mm (1/4 inch) o.d.

1	R412-0K	FMI 050	for tubing with 3.2 mm (1/8 inch) i.d.
2	R412-1K	FMI 051	for tubing with 6.4 mm (1/4 inch) i.d.
3	R412-2K	FMI 052	for tubing with 9.5 mm (3/8 inch) i.d.
4	R412-5K	FMI 053	for tubing with 1/4–28 ferrule fittings
5	H476K	FMI 054	for tubing with 3.2 mm (1/8 inch) o.d.

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
## Accessories

### Connecting cable

For connecting an MCP-CPF *Process* to a PC via the 9-pin RS232 connector.

Order No. AG0013

### Pressure control unit

consisting of a warning light and a pressure sensor (e.g. for pumping through a filter) for the pump-heads:

Pro-280, Pro-281, Pro-380, Pro-381

Order No. IS13825

### Foot switch

Order No. IS10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling bottles etc.

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## Removing the casing hood

⚠ In order to maintain the best possible seal of the casing, the pump should only be opened for replacing the fuses or changing the DIP switches. We recommend you proceed strictly as follows:

### Removing the casing hood

1. Remove the 4 screws with a Phillips screwdriver (size 2) on each side of the casing hood.
2. Spread **slightly** the casing hood at the bottom and remove it by carefully pulling it upwards.

### Mounting the casing hood

1. Before replacing the casing hood we recommend you to lubricate the sealing on the contact surface (e.g. with silicone oil spray - order no. SC0179).

In case that the two sealing ends have come apart, they must be stuck together again with Silicone sealing compound.

2. Spread the casing hood slightly at the bottom and push it carefully over the pump. Be careful that no cables are trapped.
3. Press the casing hood downwards and fasten the screws evenly from top to bottom. When tightening the screws, be careful that the profile sealing is not trapped. If the sealing should be damaged, it must be replaced, otherwise the seal may no longer be tight.

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## Removing the casing hood

### Sealing the casing hood

When the casing hood fits tight, each single screw must be removed again and stopped with silicone sealing compound. Then, the screws must be refastened again evenly from top to bottom.

➔ Remove the screw - fill the hole with silicone - fasten the screw again - remove the last bit of silicone with a rag. If necessary clean up with acetone.

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## Plug gaskets


If the interface ports are used under **IP-65 conditions**, the originally supplied gaskets must be inserted into the male plugs of the corresponding connecting cables.

➔ The interface »RS232 OUT« (male) is already equipped with a gasket.

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
- 1 Gasket for plug (15 pin, male) for analog interface (female)  
Spare part No. AD0069
- 2 Gasket for plug (9 pin, male) for RS232-IN interface (female)  
Spare part No. AD0071

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⚠ Choosing the correct interface cable (sealed plug, e.g. moulded plug) according to the IP-rating used, is entirely the user's responsibility.

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## Changing the fuses

The fuses are fixed to the control board (above) and the power supply print (below) as illustrated opposite.

### Fuse rating on control board

- a 4.0 A, slow-blow
- b 1.6 A, fast-blow

### Fuse rating power supply

- c 3.15 A, slow-blow

⚠ Switch the pump off, pull out the mains plug.

⚠ Remove the casing hood. For removing and remounting the casing hood, please refer to pages 55 – 56.

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

## Please note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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## Maintenance

We recommend the user to flush the pump-head thoroughly after pumping in order to prevent the piston from seizing by remaining particles (especially when media has been pumped that contained salt, proteins or particles).

Otherwise, the MCP-CPF *Process* requires no particular maintenance or mechanical adjustments.

Avoid pumping media containing abrasive solids.

## Service and repairs

Like all pumps, the MCP-CPF *Process* contains components that will wear over a period of time.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

## Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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## Repairs

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® distributor:

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial No., and model.

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### Mounting set

for mounting FMI pump-heads on the MCP-CPF *Process* drive, consisting of a mounting plate and assembly parts

For RH pump-heads	IS 1004
For Q pump-heads	IS 1003

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## Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Motor type DC motor

### Speed range

10–1800 rpm, digitally adjustable in steps of

0.1 rpm: 10.0 – 999.9 rpm

1 rpm: 1000 – 1800 rpm

### Flow rate

25 µl/min – 2.3 l/min (depending on pump-head)

### Differential pressure

max. 6.9 bar/100 psi (depending on pump-head)

### Remote control

via RS232 and analog interface (see pages 37 – 48)

### Mains connection

100 – 240 V<sub>AC</sub> / 50 – 60 Hz

### Fuse rating on control board

F1.6A/250V

T4A/250V

### Fuse rating power supply

T3.15A/250V

### Power consumption

Max. 100 W

### Protection rating

IP 65

### Operating conditions

Temperature +41 to 104°F

Rel. humidity max. 80%

at normal environmental conditions

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
**Assembly parts from ISMATEC®**

- 1 MCP-CPF *Process* drive
- 2 Coupling
- 3 Mounting plate
- 4 Countersunk head-screws  
(4 pcs.)
- 5 Allen key 2.5 mm
- 6 Coupling screw

**Assembly parts pump-head**


- a Pump-head  
RH00, RH0 or RH1
- b Screws (2 pcs)
- c Serrated lock washers (2 pcs.)

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### 1 Attaching the mounting plate

- Remove screws (b) and serrated lock washers (c) from the pump-head.
- Fasten the mounting-plate (3) on the pump-head with the 2 screws (b) and serrated lock washers (c).

⚠ The 4 countersinks must be concave to the pump-head.

### 2 Fixing the coupling

- Put the coupling (2) on the piston shaft. The hexagon socket screw (6) must face the flat part of the motor shaft.
- Push the coupling piece (2) towards the back of the mounting plate (3) until the inner face of the coupling meets the front of the head shaft. Tighten the coupling screw (6) with the Allen key 2.5 mm (5).

### 3 Mounting the pump-head

- The scale of the red flow control ring (f) must be pointing upwards
- Hold the mounting plate in such a way, that the screws (b) are on the top and at the bottom. Turn the coupling, so that the two coupling halves of pump-head and drive engage. Insert the pump-head assembly and push it on as far as it will go.
- Drive in the countersunk head-screws (4).

⚠ Entering the pump-head see pages 15 and 51.

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

## Type Q Mounting a pump-head

### Assembly parts from ISMATEC®

- 1 MCP-CPF *Process* drive
- 2 Coupling
- 3 Allen key 2.5 mm
- 4 Mounting plate
- 5 Socket head cap screws, 4 pcs. (for attaching the head to the mounting plate)
- 6 Socket head cap screws, 4 pcs. (for attaching the mounting plate to the drive)
- 7 Allen key 3.0 mm
- 8 Coupling screw


### Assembly parts pump-head

- a Pump-head Q0, Q1, Q2, Q3
- b Mounting base

#### 1 Detaching the mounting base

- Take off the four screws at the bottom of the mounting base (b).
- Remove the mounting base (b) from the pump-head (a). You must do this in order to reach the two lower screws (5) when you fix the pump-head to the mounting plate (see figure 4).

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



3

## 2 Fixing the pump-head

Fix the pump-head (a) with the 4 head cap screws and washers (5) to the mounting plate (4). For tightening use an Allen key 4.0 mm (7).

⚠ Please note that the slightly protruding cylinder bolts on one side of the plate serve to centre the pump-head. The upper contour of the mounting plate must be level with the top of the pump-head plate.

## 3 Mounting the coupling

– Put the coupling (2) on the piston shaft. The hexagon socket screw (8) must face the flat part of the motor shaft. Push the coupling piece (2) towards the back of the mounting plate (4) until the inner face of the coupling meets the front of the head shaft.

Tighten the coupling screw (8) with the Allen key 2.5 mm (3).

– Now tighten the complete unit with head cap screws using the 2.5 mm Allen key (6).

⚠ **Entering the pump-head see pages 15 and 51.**

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## **IDEX Health & Science**

Home of Upchurch Scientific &  
Ismatec Products

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Ask for the latest  
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH