**Pumpenköpfe für**  
**MCP *Standard + Process***  
**BVP *Standard + Process***

Pro 280 / Pro 281  
Pro 380 / Pro 381  
360 / 380 / 380 AD  
CA 4 / CA 8 / CA 12  
FMI Q0 / Q1 / Q2 / Q3  
Easy-Load  
Standard  
Quickload  
MS/CA 4-12 / MS/CA 8-6  
PTFE Tube 2 mm / 4 mm  
SB 2V / 3V  
WM 5

Deutsch

**Pump-heads for**  
**MCP *Standard + Process***  
**BVP *Standard + Process***

Pro 280 / Pro 281  
Pro 380 / Pro 381  
360 / 380 / 380 AD  
CA 4 / CA 8 / CA 12  
FMI Q0 / Q1 / Q2 / Q3  
Easy-Load  
Standard  
Quickload  
MS/CA 4-12 / MS/CA 8-6  
PTFE Tube 2 mm / 4 mm  
SB 2V / 3V  
WM 5

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## Please note

For detailed information on warranty and safety precautions please refer to the operating instructions for your pump drive.

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## Safety precautions

It is very important that the following safety precautions are observed:

### **Before starting up**

For changing or mounting tubing the pump must be switched off.

### **During operation**

Be particularly cautious that no parts of your body (e.g. fingers, hair, etc.) or loose objects (e.g. cables, tubing, etc.) can be trapped by any rotating parts of the pump (e.g. pump-rollers).

Cassettes and pump-head must not be opened while the pump is in operation.

### **General risks**

ISMATEC® do not admit responsibility for the handling of chemicals.

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
## Explanations of terms being used

- 1 Drive**  
MCP *Standard* or *Process*  
BVP *Standard* or *Process*
- 2 Release lever**
- 3 Fixing slits (2 pieces)**
- 4 Fixing pins (2 pieces)**
- 5 Mounting plate  
with (fixed) pump-head**
- 6 Only for MCP models**  
ID-Code (for identification of the pump-head)

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-----  
Caution: rotating  
pumphead creates a pinch  
and crush hazard.

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### Mounting with mains switch «OFF»

1 Push the mounting plate with the 2 fixing pins into the corresponding slits of the drive.

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
2 Turn the roller-head and press it gently against the pump-drive till the coupling engages.

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3 Press the mounting plate further against the drive and turn it now clockwise until the release lever snaps in.

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
«ON»

### 1 BVP alignment

- Run/Stop button on STOP
- for mounting press and maintain pressure on SLOW (Drive is turning at 10 rpm)

### 1b MCP alignment

- PUMP (set by MODE key)
- for mounting press and maintain pressure on MAX/MIN-key (Drive is turning at 10 rpm)

### 3 With free hand

- push the 2 pins into the fixing slits of the drive
- press the mounting plate gently against the drive till the head engages (rollers start to turn)

### 4 Release hold of SLOW or MAX/MIN button

- Press the mounting plate against the drive and turn it clockwise until the release lever snaps in.

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
## Decoupling the pump-head

Same procedure for MCP and BVP

- 1 Mains switch «OFF»
- 2 Push release lever to the left.
- 3 Turn mounting plate as far as possible to the left.  
➔ Now pull the pump-head straight out of the drive coupling (maybe by slightly rotating back and forth).

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


## ⚠ Pump-head identification

(➔ for each of the 4 programs)

In order to obtain correct pumping and dispensing values the individual identification code (ID-Code) of the mounted pump-head must be entered, **especially when a new pump is used for the first time** and each time the pump-head is changed.

- 1 Power switch »OFF«
- 2 Keep CAL-key pressed and switch power »ON«
- 3 Press the ▲▼ keys till the correct ID-Code of the mounted pump-head is set. For the ID-Code see the back of the pump-head or the table on pages 14 and 15.
- 4 Enter the correct ID-Code with the CAL-key  
(MODE display switches automatically to TUBE I.D.)

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
## Only for MCP Process from version 3.0

### Entering the pump-head

➔ In order to obtain correct pumping and dispensing values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

- 1 Select the program (1 to 4).
- 2 Switch the pump off.
- 3 Keep the settings key pressed and switch the pump on.  
➔ Pump is now in settings menu.
- 4 Select **HEAD** with the ▲▼ keys  
Head = ID code  
(Pump-head - Can be set separately for each program)

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8 Settings-Menu mit escape-Taste verlassen.

**Only for MCP Process  
from version 3.0 on**

5 Confirm with ok key.

6 Select the identification code (ID code) of the mounted pump-head with the ▲▼ keys. For the ID-Code see the back of the pump-head or the table on page 14–15.

7 Confirm with ok key  
➔ **In each of the 4 programs a different pump-head can be stored!**

8 Quit the settings menu with the escape key.

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

⚠ Before starting to pump it is essential  
to enter the following settings in order to  
obtain correct flow rates and dispensing  
volumes.

### Peristaltic pump-heads Entering the tubing i.d.

- 1 Change mode to  
TUBE I.D. mm (for setting the tubing  
inner diameter)
- 2 Press CAL-key  
(display blinks)
- 3a Compare blinking figure with  
required tube i.d.
- 4 Press ▲▼keys until correct i.d. is  
displayed
- 5 Confirm entered value with CAL-key  
➔ Select required operation mode with  
MODE-key

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## Entering the tubing i.d. or stroke angle

### FMI piston pump-heads Setting the stroke volume

In mode TUBE I.D. mm the piston angle can be entered as stated above from position 1 to 5:

- 3b Set piston angle (between 1 and 10)  
→ The same angle must also be manually set on the pump-head

4–5 (page 12)

- 6 Change to PROGRAM with the MODE key and select the next program with the CAL and ▲▼ keys (confirm with CAL).

→ Repeat procedure till the ID-Code(s) of the mounted pump-head(s) and the tubing i.d. are **entered in all 4 programs**.

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\*Only for programming via PC and RS232 interface.

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IDEX Health & Science  
Home of Upchurch Scientific &  
Ismatec Products  
619 Oak Street  
Oak Harbor, WA 98277  
United States  
Tel: +866 339 4653  
Fax: +360 679 3830  
fluidics.americas@idexcorp.com

IDEX Health & Science  
Ismatec Laboratoriumstechnik GmbH  
Futtererstrasse 16  
97877 Wertheim  
Germany  
Tel: +49 9377 9203-0  
Fax: +49 9377 1388  
fluidics.europe@idexcorp.com

www.idex-hs.com

Verlangen Sie den  
ISMATEC-Katalog.

Ask for the latest  
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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