MCP Process

Pumpen- Antrieb

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ISM 915

Pump drive

Microprocessor
controlled

ISM 915

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ISM 915

Schutzgrad IP 65

Für Pumpenköpfe:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS 3
MS/CA 4-12 / MS/CA 8-6
PTFE Membran HF/LF
PTFE Tube 2 mm / 4 mm
SB 2V / 3V / 5V
WM 5

Protection rating IP 65

For pump-heads:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS 3
MS/CA 4-12 / MS/CA 8-6
PTFE Diaphragm HF/LF
PTFE Tube 2 mm / 4 mm
SB 2V / 3V / 5V
WM 5

Classe de protection IP 65

Pour têtes de pompe:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS 3
MS/CA 4-12 / MS/CA 8-6
PTFE Diaphragme HF/LF
PTFE Tube 2 mm / 4 mm
SB 2V / 3V / 5V
WM 5

Deutsch

English

Français

ISMATEC®

IDEX
HEALTH & SCIENCE

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry.

As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.

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Safety precautions

- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend placing the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 8 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Pump drive: 2 years

Pump-head: 1 year

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- MCP *Process* drive Magnet
included Order No. ISM 915
incl. integrated power cord,
2 m long, with IEC 320 plug
(male)
- 1 power cord with an IEC 320
connector (female) and a country
specific mains plug
- 2 gaskets (in plastic bag) for
analog (AD0069) and
RS232-IN (AD0071) interface
- Operating manual

If ordered:

- pump-head with
mounting plate
- software
- other accessories

Please check the package
and its contents for transport
damage. If you find any signs of
damage, please contact your local
ISMATEC® representative
immediately.

➔ Complaints can only be accepted
within 8 days from receipt of the
goods.

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
⚠ Please note

For using the analog or RS232- interface under IP 65 conditions please refer to page 58.

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Rear panel

- 1 Mains switch (on/off)
- 2 Power cord
- 3 RS232 IN (female)
- 4 RS232 OUT (male)

5 Analog interface

■ input for:

- speed control
0 – 5 V or 0 – 10 V, and
0 – 20 mA or 4 – 20mA
- rotation direction
- Run/Stop
- Foot switch
- 2 digital inputs (TTL-level)

■ Output for:

- speed 0 – 10 V_{DC}
or 0 – 7.2 kHz
- 2 universal outputs

⚠ Mains voltage

100 - 240 V_{AC}
50 – 60 Hz
no adjustments necessary

Power consumption

Max. 100 W

Fuse rating on control board*

F1.6A/250V
T4A/250V

Fuse rating on power supply print*

T3.15A/250V

* See also page 59

⚠ Socket/Power cord

The socket must be earthed. (protective conductor contact)

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Starting the pump

- Mount the pump-head according to the mounting instruction manual supplied with the pump-head.
- Enter the ID Codes of the mounted pump-heads in the program memory currently used (see pages 15-16 and 53).
- Enter the tubing i.d. in each program used (see page 19).
- If using an FMI piston pump-head, adjust the angle for the piston stroke volume (see page 19).
- Insert the tubing into the peristaltic pump-head.
- Connect the pump tubing to the system.
- Connect the pump to the mains and switch it on with the power supply switch.

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5

6

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	●Time	◦Volume
PAUSE	◦Time	◦TOTAL

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Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Identification Code (ID-Code) of the pump-head entered in the program last used (e.g. pump-head 380AD = 380A)
- 4 Roller back-steps
➔ Are only displayed if at least 1 or more back-steps are entered
- 5 Program last used
- 6 Displays the operating mode last used


⚠ Before starting the pump the first time or after mounting a new pump-head, the ID-Code must be entered in each program used with this pump-head (see page 15–16 and 53).

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TUBE I.D.	mm	PROGRAM
PUMP	rpm	Flow rate
DISP	Time	Volume
PAUSE	Time	TOTAL

Anzeige der Betriebsart
LED's for operating modes
Affichage du mode d'exploitation

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Operating panel



- 1 Digital LED display
- 2 Reduce value
* settings
- 3 Increase value
- 4 **RUN/STOP**
– Starts and stops the pump.
– Interrupts or continues a dispensing cycle (pressing the RESET key aborts the interrupted dispensing process).
➔ The Start/Stop function can also be triggered by an external device such as a foot switch.
- 5 **CAL**
Calibrating key for the flow rate or dispensing volume (pages 23/26).
* ok
- 6 **MAX/MIN/RESET**
Multi-function key for:
– max. speed¹⁾ (maintain pressure on this key while pump is running) ideal for filling or emptying the tubing system
– min. speed (when pump is idle) turns slowly
– Stops an interrupted dispensing cycle for good
– Resets the accumulated volume in mode TOTAL to zero. Resets a flow rate or volume calibration to the default value (page 27)
* default

¹⁾ functions in the operating modes PUMP rpm, PUMP Flow rate and TOTAL

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TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	◦Volume
PAUSE	◦Time	◦TOTAL

Anzeige der Betriebsart
LED's for operating modes
Affichage du mode d'exploitation

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
Operating panel

- 7 Rotation direction**
Changes the rotation direction
➔ Each time the rotation direction is changed, TOTAL is reset to »0«.
* escape
- 8 MODE**
Changes between the operating modes
- 9 LEDs for active operating mode**
- **TUBE I.D. mm**
the tubing inner diameter
 - **PROGRAM**
Program selection 1 – 4
 - **PUMP rpm**
Pumping by speed
1.0 – 240.0 rpm
 - **PUMP Flow rate**
Pumping by flow rate in µl/min, ml/min or l/min
 - **DISP Time**
Dispensing by time
0.10 s to 999 h
 - **DISP Volume**
Dispensing by volume in µl, ml or litres
 - **PAUSE Time**
for intermittent dispensing from 0.10 s up to 999 h
 - **TOTAL**
Read-out of totally delivered volume in µl, ml or litres

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Control keys

- a RUN/STOP**
Starts and stops the pump.
- b MODE**
Changes between operating modes (see page 12).
- c Rotation direction**
Changes the rotation direction.
➔ This function is blocked when using the diaphragm pump-head.
➔ Each time the rotation direction is changed, TOTAL is reset to »0«.
- d CAL**
Key for saving a set value.
- e MAX/MIN/RESET**
Multi-function key for:
➔ max. speed (when pump is running)
➔ min. speed (when pump is idle) (turns slowly, e.g. for mounting a pump-head)
➔ Reset button resets the calibration to the standard values (default), page 27
➔ Reset button for accumulated volume in Mode TOTAL
- f Increment/Decrement keys**
▲ = increase value
▼ = decrease value
Maintaining pressure on ▲ or ▼ key changes display read-out into fast mode.

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e ▲▼ keys

- Select the basic settings
- Select value

➔ For selecting the basic settings refer to the following pages.

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Selecting the basic settings

- 1 Select the program (p. 20)
- 2 Switch the pump off
- 3 Keep the settings key pressed and switch the pump on

4 Select the basic settings by using the ▲▼ keys

Can be set separately for each program

- Cycles (Number of dispensing steps)
- Backstep (Number of roller back-steps)
- Head = ID code (pump-head)
- FS (foot switch)
- Digits (display)
- Init (resets all basic settings to the default values)

Valid for all 4 programs

- Init4 - in all 4 programs the above mentioned basic settings as well as
- Address (RS232, pump address) and
- Baud rate (RS232) are reset to their default values.

- 5 ok = confirm your selection
 - 6 Select the value by using the ▲▼ keys, or
 - press **escape** for cancelling wrong input
 - certain values can be reset to the default value by using the **default** key (for values, refer to the following pages)
 - 7 ok = confirm entered value
 - 8 escape = quit the menu basic settings
- ➔ For selecting the basic settings refer to the following pages.

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
Program specific basic settings

→ **Individually** adjustable for each of the 4 programs.

Cycles

Number of dispensing steps by time or volume 0..9999 (0 = Default value, infinity) (See intermittent dispensing on pages 29 and 30).

Backsteps

Number of roller back-steps (for drip-free dispensing in the mode DISP Volume) from 1 .. 100 (0 = Default value).

Head

Identification code (ID code) of the pump-head. For the ID-Code see the back of the pump-head or the table on page 53.

→ In order to obtain correct pumping and dispensing values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

→ In each of the 4 programs a different pump-head can be stored.

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Foot switch

- Toggle (default value)
 - Step on foot switch
 - = pump starts running
 - Press foot switch again
 - = pump stops
- Direct
 - the pump runs as long as the foot switch is pressed

Digits (display)

Number of significant digits, 3 or 4 (3 = default value)

Initializing

Pressing the ok key resets the following basic settings in the **currently** used program to the default values:

Mode:	PUMP rpm
Speed:	100.0 rpm
Dispensing time:	4.50 sec.
Dispensing volume and flow rate:	depends on the pump-head
Pause time:	2.00 sec.
Number of dispensing cycles:	0

→ The pump-head ID code entered in the currently used program remains stored.

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Entering the tubing i.d. or stroke angle

⚠ Before starting to pump it is essential to enter the following settings in order to obtain correct flow rates and dispensing volumes.

ONLY for Peristaltic pump-heads

- 1 Change mode to TUBE I.D. mm (for setting the tubing inner diameter)
- 2 Press CAL-key (display blinks)
- 3a Compare blinking figure with required tube i.d.
- 4 Press ▲▼ keys until correct i.d. is displayed
- 5 Confirm entered value with CAL-key
➔ Select required operation mode with MODE-key.

ONLY for FMI piston pump-heads

Setting the stroke volume

In mode TUBE I.D. mm the piston angle can be entered as stated above from position 1 to 5:

- 3b set piston angle (between 1 and 10)
➔ The same angle must also be manually set on the pump-head.

ONLY for Diaphragm pump-heads

Tube i.d. is not relevant. In mode TUBE I.D. mm the display shows »0.0«.

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Program selection

When switching the pump on, it always selects the previously used program.

1 Change mode to PROGRAM by using the MODE-key

2 Press CAL button, currently set program starts to blink
Change to the required program with the ▲▼ buttons

3 Confirm with the CAL-key

4 The pump automatically returns to the last entered operating parameters of the selected program

⚠ **From now on, any changes carried out in the operating modes are automatically stored in the currently selected program.**

TOTAL Volume

Change with MODE key to TOTAL. The totally delivered volume is displayed.

Three display readings are available:

Microlitre: e.g. 12.0 µl = **12.0µ**

Millilitre: e.g. 1200 ml = **1200**

Litre: e.g. 12.0 l = **12.0L**

(Display in litres from 9999 ml)

For zero-setting, either

- press the reset key twice, or
- switch the pump off.

➔ Each time the rotation direction is changed, TOTAL is reset to »0«.

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1

MODE

TUBE I.D.	◊mm	◊ PROGRAM
PUMP	◊rpm	◊ Flow Rate
DISP	◊Time	◊ Volume
PAUSE	◊Time	◊ TOTAL

2

CAL

Pro4

Pro2

3

CAL

4

TUBE I.D.	◊mm	◊ PROGRAM
PUMP	◊rpm	◊ Flow Rate
DISP	◊Time	◊ Volume
PAUSE	◊Time	◊ TOTAL

z.B. / e.g. / p.e. PUMP rpm

MODE

TUBE I.D.	◊mm	◊ PROGRAM
PUMP	◊rpm	◊ Flow Rate
DISP	◊Time	◊ Volume
PAUSE	◊Time	★ TOTAL

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
Pumping by drive speed

- 1 Change mode to PUMP rpm, 1.0 – 240.0 rpm, adjustable in steps of 0.1 rpm
- 2 – Enter required speed with ▲▼ keys
– Start pump with RUN/STOP
➔ The speed can also be adjusted while the pump is running.
- 3 Maintain pressure on MAX/MIN-key while the pump is running.
➔ Pump is running at maximum speed (ideal for fast filling or emptying the system).
- 4 Maintain pressure on MAX/MIN-key when pump is idle (power switched on)
➔ Pump is turning at 10 rpm (for mounting the pump head).

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1

TUBE I.D.	◊mm	◊PROGRAM
PUMP	◊rpm	★Flow Rate
DISP	◊Time	◊Volume
PAUSE	◊Time	◊TOTAL

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Pumping by flow rate

1 Change mode to PUMP Flow rate.

2 Enter the required flow rate with ▲▼ keys (is displayed in µl/min or ml/min).

Start pump with RUN/STOP.

➔ The flow rate can also be adjusted while the pump is running.

➔ For an accurate flow rate we recommend to calibrate the pump (see page 23).

➔ During the pumping Process the following values can be retrieved with th MODE-key:

- Rotation speed (PUMP rpm)
- Totally pumped volume (TOTAL)

3 Maintain pressure on MAX/MIN-key while the pump is running.

➔ Pump is running at maximum speed (ideal for fast filling or emptying the system).

4 Maintain pressure on MAX/MIN-key when pump is idle (power switched on).

➔ Pump is turning at 10 rpm (for mounting the pump-head).

⚠ Entering the ID-Code of the currently mounted pump-head (pages 16 and 53) and the inner diameter of the tubing used allows the user to set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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Calibrating the flow rate

1 Change mode to PUMP Flow rate.
Enter the required flow rate with the ▲▼ keys.

2 Change mode to DISP Time

- Enter 60 seconds by using the ▲▼ keys.
- Start pump with RUN/STOP.

➔ Pump stops automatically after 60 seconds.

➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:

3 Return to mode PUMP Flow rate.

4 Press the CAL-button (displayed value blinks)
Enter the weighed or measured value with the ▲▼ keys.

5 Confirm with the CAL-key (the flow rate setting returns automatically to the initially entered set point).

6 Start with the RUN/STOP-key.

➔The flow rate setting depends on the number of pump rollers and the tubing i.d. It is possible that the required set point cannot be entered accurately.
If necessary choose a tubing with a smaller i.d. or a pump-head with more rollers.

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Dispensing by time

The dispensing time can be entered from 0.10 s to 999 h.

- Change mode to DISP Time.
 - Enter the required dispensing time with the ▲▼ keys.
 ➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.
- | | | |
|----------|---------------|-------------------|
| Seconds: | 0.10" – 99.9" | (in 0.01 s steps) |
| | 0.1" – 899.9" | (in 0.1 s steps) |
| Minutes: | 15' – 899' | (in 1 min steps) |
| Hours: | 15 h – 999 h | (in 1 h steps) |
- Start pump with RUN/STOP


- ➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect from the subsequent dispensing cycle.
- ➔ During the dispensing process the following values can be retrieved by pressing the MODE-key:
 - speed (PUMP rpm)
 - flow rate (PUMP Flow rate)
 - totally dispensed volume (TOTAL)
- ➔ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.
- ➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible. We recommend to use a tubing with a smaller i.d. and to increase the dispensing time.

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Dispensing by volume


- 1 Change mode to DISP Volume
- 2 Use the ▲▼ keys for entering the required dispensing volume. Three display readings are available:
 Microlitre: e.g. 12.5 µl = **12.5µ**
 Millilitre: e.g. 230.5 ml = **230.5**
 Litre: e.g. 12.5 Litre = **12.5L**
 (Display in litres from 9999 ml)
- 3 Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

- ➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.
- ➔ With the ▲▼ keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears shortly in the display taking effect, however, only with the subsequent dispensing step.
- ➔ During the dispensing process the following values can be retrieved with the MODE-key:
 - Rotation speed (PUMP rpm)
 - Flow rate (PUMP Flow rate)
 - totally dispensed volume (TOTAL)
- ➔ With the ▲▼ keys the rotation speed and flow rate can be changed even during the dispensing process.

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Calibrating the volume

- 1 Change MODE to DISP Volume.
 - ➔ Use the ▲▼ keys for entering the required dispensing volume.
- 2 Start with RUN/STOP
 - ➔ Pump stops automatically.
 - ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.
 - ➔ Measure the dispensed liquid by volume or weight.
- 3 Stay in the same mode (DISP Volume) and press the CAL-key (display blinks).
- 4 Enter the effectively dispensed volume with the ▲▼ keys.
- 5 Save setting with CAL-key
 - ➔ Pump is calibrated and returns automatically to the range of the initially entered set point.
- 6 Start pump with RUN/STOP
 - ➔ Depending on the application conditions repeated calibration may be necessary.
 - ➔ The volume setting depends on the number of pump rollers and the tubing i.d. It is possible that the required set point cannot be entered accurately. If necessary choose a tubing with a smaller i.d. or a pump-head with more rollers.

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➔ At the same time, the volume is also reset to the default value.

Calibrant reset of volume

- 1 Change mode to DISP Volume.
 - 2 Press the CAL button (displayed value blinks).
 - 3 Press Reset/MAX/MIN button
The volume is reset to the default value.
 - 4 Confirm with the CAL button.
- k At the same time, the flow rate is also reset to the default value.

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






Dispensing by volume within a pre-set time

- Change to mode DISP Time. Enter the required time with the ▲▼ keys (0.10 s – 999 h) (see page 24).
- Change mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys.
- Push the MODE-button
➔ Both LEDs light up simultaneously.
- Start pump with RUN/STOP
➔ The pump dispenses the required volume within the pre-set time.

a Changes / Corrections

- Press the Reset key
➔ For changing the dispensing time or dispensing volume, MODE returns directly to DISP Time
➔ Repeat steps 1 to 3

b Error message

- LEDs are blinking alternately
9999 = volume too large
time too short
1111 = volume too small
time too long

Despite volume and time adjustments, the LEDs keep on blinking:
➔Change the tubing i.d. or the pump-head model in order to obtain another flow rate.
➔Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible. We recommend to use a tubing with a smaller i.d. and to increase the dispensing time.

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
Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

- 1 Go to mode PAUSE Time
 - Enter the pause time (between 0.10 s and 999 h) with ▲▼ keys.
- 2 Press MODE key once
 - ➔ Both LEDs for PAUSE and DISP Time light up simultaneously
 - Enter the dispensing time (0.10 s to 999 h) with ▲▼ keys
- 3 Start pump with RUN/STOP key
 - ➔ Pump dispenses at intervals until the RUN/STOP button is pushed again
 - ➔ Active mode blinks

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the subsequent dispensing cycle.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- rotation speed (PUMP rpm)*
- flow rate (PUMP Flow rate)*
- totally dispensed volume (TOTAL)

*The rotation speed and flow rate can be changed with the ▲▼ keys.

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

Please note

For entering the number of dispensing cycles see page 31.

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Please note

For entering the number of dispensing
cycles see page 31.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

1 Change mode to PAUSE Time

Enter pause time (between 0.10
s and 999 h) with
the ▲▼ keys.

2 Press MODE-key twice.

➔ Both LEDs for PAUSE and
DISP Volume light up
simultaneously. Enter the volume
with the ▲▼ keys.

3 Start pump with the
RUN/STOP key.

➔ Pump dispenses at intervals
until the RUN/STOP button is
pushed again.

➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing
volume can be changed even du-
ring the dispensing process. The
newly entered set point appears
shortly on the display taking ef-
fect, however, only with the next
dispensing cycle.

➔ During the dispensing process the
following values can be retrieved
with the MODE-key:

- rotation speed (PUMP rpm)*
- flow rate (PUMP Flow rate)*
- totally dispensed volume
(TOTAL)

*The rotation speed and flow rate
can be changed with the ▲▼ keys.

If necessary, calibrate the volume as
indicated on page 26.

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

settings

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	*Time	◦Volume
PAUSE	*Time	*TOTAL

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	*Volume
PAUSE	*Time	*TOTAL

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

After entering the dispensing cycles (basic settings page 16) the pump can be started with the Run/Stop key.

According to the pre-set pause and dispensing time or dispensing volume, the pump now dispenses at intervals and stops when the number of dispensing cycles has elapsed.

➔During the pauses the remaining number of dispensing cycles is indicated on the display.

- During long pauses the remaining pause time can be retrieved by pressing the CAL key.
- The dispensing or pause cycle can be interrupted by pressing RUN/STOP.
- Pressing RUN/STOP again continues the process.
- Pressing the RESET key aborts the interrupted dispensing process.

➔The number of dispensing cycles remains stored until changes are made in the basic settings (page 16).

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
Drip-free dispensing

Reverse roller steps (backsteps) allow drip-free dispensing (1 to 100 roller steps). We specially recommend this operating mode if corrosive media is pumped or in case that only few sample material is available.

- 1 Switch the pump off.
- 2 Maintain pressure on the settings key and switch the pump on again.
- 3 Select BACKSTEPS by using the ▲▼ keys (basic settings, page 16).
ok = confirm selection
- 4 Enter the number of reverse roller steps with the ▲▼ keys.
➔ In case of wrong input return to BACKSTEPS with the escape key.
ok = confirm selection
- 5 escape = quit menu basic settings

Reverse roller steps are active in the following modes:

- a Dispensing by volume (see page 25)
- b Dispensing a volume within a pre-set time (see page 28)
- c Intermittent dispensing by volume (see page 30)

⚠ Reject the first dispensing step. Afterwards, the number of reverse roller steps are compensated. Deleting the reverse roller steps

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

△ Please note

For our pressure control unit, see page 48.

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Pumping against pressure

The MCP *Process* can be used for continuous duty at a differential pressure of max. 1.5 bar. For short-time operation higher differential pressures may be managed depending on the mounted pump-head and tubing i.d. (smaller i.d. is preferable).

Depending on the application conditions some FMI pump-heads allow operations at differential pressures of up to 6.9 bar.

In case of any doubts please contact your ISMATEC® agent.

When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

- a Tubing used with cassettes do not need to be removed. Releasing the cassette on one side is sufficient.
- b The single-channel pump-heads 360, 380 and Easy Load allow opening the tube-bed in order to release the tubing pressure.

△ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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Overcurrent protector

The drive MCP *Process* features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL S = Speed set-point cannot be reached

OL I = Current Overload

Whenever this situation occurs, the pump must be switched off immediately. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure or tension lever tightened too much, etc.).

➔ Only when the operating conditions have been changed accordingly may the pump be started again.

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Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running in period.

If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

- In normal operation (pin 2 open) pin 3 operates the foot-switch (settings see page 17)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

Pin 5, speed IN

External speed control (0 – 5V, 0 – 10V, 0 – 20mA, 4 – 20mA)

Input impedance and input range can be selected via a dip-switch inside the pump(see page 38)

Pin 6, input 1

Pin 8, input 2

Digital inputs (TTL-level). Active when connected to pin 1. Are only retrievable via RS232 interface; e.g. with a software like the *ProgEdit*.

Pin 7, +26V_{DC}

About +26 V_{DC} are available (max. current 1A).

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Analog interface

Pin 9, speed OUT

The default setting is 0 – 10 V_{DC}, proportionally to the motor speed 0 – 240.0 rpm.

Alternatively a frequency range from 0 – 7.2 kHz is available. Alternatives with sliding switch S2 inside the pump (see page 38)

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.5 A)

Pin 13, speed intern

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated (Normal operation, i.e. pin 2 is open) Pin 13 serves as **auto-start function**.

If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated

(Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

Pin 14, output 1 (PWM)

Pin 15, output 2 (PWM)

Digital outputs in connection with pin 7 (+26 V_{DC}), implemented as **pulse-width-modulation** (open collector). Can only be activated via RS232 interface; e.g. with a software like the *ProgEdit*.

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
Please note

The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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This switch affects Pin 9, speed OUT

Position A: 0 – 10 V_{DC} (Standard)
Position B: 0 – 7.2 kHz

➔ In order to maintain the optimal casing seal, the pump should not be opened without a particular reason. This is, however, necessary in order to gain access to the switches. For this purpose we recommend that you proceed exactly as indicated on pages 57 – 58.

➔ Make sure that the pump is disconnected from the mains supply.

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
Please note

When switching on the MCP *Process* drive, the serial interface will only respond after a delay of 3 seconds.

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Kommunikationsparameter Communication parameters Paramètres de communication

Standardeinstellung
Default setting
Réglage standard

**9600 baud, 8 bit,
1 stopbit, no parity**
(DIP Switch 1 off)

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1 – 8).

It is factory-set at 1 and can be changed with @.

The address will be stored permanently (even after the pump has been switched off).

The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterisk *.

Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address (see also the basic settings on page 18).**

➔ Alternatively, the addressing can also be accomplished via the serial interface.

For that, each pump must firstly be connected individually to a PC and the address must be allocated with @. Only then may the pumps be daisy-chained.

Each pump requires a definite address in order to avoid interference of their responses.

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Each pump requires an individual address.
(see page 39).

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➔ Please ask for more detailed software-information.

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Programming Software

Software *ProgEdit*

- Create your individual application program.
- Download the program via RS232 interface to your MCP *Process* pump drive.
- Now run your application PC-independently wherever needed.

System requirements:

3 1/2" floppy drive, Win95 or higher, Pentium or higher

Order No. SOF104

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.
<http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website:
<http://www.ismatec.com>

Accessories

Pressure control unit

consisting of warning light and pressure sensor (e.g. for pumping against a filter)

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
Accessories

Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments

Order No. ISA IM0015C168
Order No. PCI IM0030

Adaptor 25 F / 9M

25 pins, female / 9 pins, male for IM0015C168

(1 adaptor for each controlled instrument)
Order No. XC0009

Extension cable 9F / 9M

9 pins, female / 9 pins, male for XC0009 (2 m long)

Order No. AG0013

➔ Plug gaskets for IP 65 rating, see page 56

Foot switch

Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

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
MS/CA Click'n'go IS 3510
MS/CA Anpresshebel ¹⁾ IS 0649

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CA Anpresshebel ¹⁾ IS 3820
Adapter **PVDF** IS 3861
(pro Kasette 2 Stk. bestellen)

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Spare-cassettes

Spare-cassettes in POM-C

CA Click'n'go IS 3710
CA pressure lever ¹⁾ IS 0122
Adaptor **POM-C** IS 0123
(order 2 adaptors per cassette)

MS/CA Click'n'go IS 3510
MS/CA pressure lever ¹⁾ IS 0649

Spare-cassettes in PVDF

CA pressure lever ¹⁾ IS 3820
Adaptor **PVDF** IS 3861
(order 2 adaptors per cassette)

MS/CA pressure lever ¹⁾ IS 3629

¹⁾ The cassettes with pressure lever are available on request. This type of cassette provide better results at elevated differential pressure conditions (>1 bar).

⚠ Click'n'go-Cassettes

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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MS/CA Click'n'go IS 3510
MS/CA levier de pression ¹⁾ IS 0649

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

Standard-Schläuche
Standard tubing / Tubes au mètre

Pump Tubing

Quality Tubing for Pumping
in Laboratories and Production

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Please refer to our detailed tubing
documentation or our website

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Accessories

Materials POM-C

- Polyoxymethylene-Copolymer
- good chemical resistance to many organic solvents and strong alkaline chemicals
 - is affected by strong acids and oxidizing substances
 - UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

Materials PVDF

- Polyvinylidene fluoride
- very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
 - not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
 - stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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Mounting

Each pump-head is supplied with an instruction manual.

Do not forget to enter the correct ID code (see page 16 and 53)

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Interchangeable Pump-heads

These pump-heads can be ordered separately and mounted on the MCP *Process* drive.

➔For the piston pump-heads the suffix „-W“ resp. „-WT“ means: with isolation (W) or temperature gland (WT) for crystallizing media.

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

⚠ Must be entered for each of the 4 programs using the particular pump-head.

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Please note

In the following tables the indicated values are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®

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Opening the casing hood

⚠ In order to maintain the best possible seal of the casing, the pump should only be opened for replacing the fuses or changing the DIP-switches. We recommend to proceed strictly as follows:

Removing the casing hood

1. Remove the 4 screws with a Phillips screw driver (size 2) on each side of the casing hood.
2. Spread **slightly** the casing hood at the bottom and remove it by carefully pulling it upwards.

Mounting the casing hood

1. Before replacing the casing hood we recommend you to lubricate the sealing on the contact surface (e.g. with silicone oil spray). In case that the two sealing ends have come apart, they must be stuck together again with Silicone sealing compound.
2. Spread the casing hood slightly at the bottom and push it carefully over the pump. Be careful that no cables are trapped.
3. Press the casing hood downwards and fasten the screws evenly from top to bottom. When tightening the screws, be careful that the profile sealing is not trapped. If the sealing should be damaged, it must be replaced, otherwise the seal may no longer be tight.

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Opening the casing hood

Sealing the casing hood

When the casing hood fits tight, each single screw must be removed again and stopped with silicone sealing compound. Then, the screws must be refastened again evenly from top to bottom.

➔ Remove the screw - fill the hole with silicone - fasten the screw again - remove the last bit of silicone with a rag. If necessary clean up with acetone.

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Plug gaskets

If the interface ports are used under **IP 65 conditions**, the originally supplied gaskets must be inserted into the male plugs of the corresponding connecting cables.

➔ The interface »RS232 OUT« (male) is already equipped with a gasket.

1. Gasket for plug (15 pin, male) for analog interface (female)
Spare part No. AD0069
2. Gasket for plug (9 pin, male) for RS232-IN interface (female)
Spare part No. AD0071

⚠ Choosing the correct interface cable (sealed plug, e.g. moulded plug) according to the IP-rating used, is entirely the user's responsibility.

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


Changing the fuses

The fuses are fixed to the control board (above) and the power supply print (below) as illustrated opposite.

Fuse rating on control board

- a T4A/250V
- b F1.6A/250V

Fuse rating power supply

- c T3.15A/250V

⚠ Switch the pump off, pull out the mains plug.

⚠ Remove the casing hood. For removing and remounting the casing hood, please refer to pages 57 and 58.

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Maintenance

Provided the MCP *Process* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective MCP *Process* drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., and model.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Drive

Motor type DC motor
Speed range 1.0–240.0 rpm
digitally adjustable
in steps of 0.1 rpm

Differential pressure

1.0–1.5 bar (14.5- 21.8 psi) (depending on pump-head) for short-time use and with small tubing i.d.'s higher pressures are possible.

Remote control

via RS232 and analog interface (see page 36–46)

Mains connection

100 - 240 V_{AC} / 50–60 Hz

Fuse rating on control board

F1.6A/250V
T4A/250V

Fuse rating power supply

T3.15A/250V

Power consumption

Max. 100 W

Protection rating

IP 65

Operating conditions

Temperature 41 to 104°F,
Rel. humidity max. 80%
- not condensing, at normal environmental conditions

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⚠ Please note

The Quick Start instruction on your right helps you control the drive instantly.

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* TOTAL MAX MIN RESET 2x, or/oder OFF

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH