MCP Standard

Pumpen- Antrieb

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ISM 404

Pump drive

Microprocessor
controlled

ISM 404

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ISM 404

Schutzgrad IP 30

Für Pumpenköpfe:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Protection rating IP 30

For pump-heads:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Classe de protection IP 30

Pour têtes de pompe:
Pro 280 / Pro 281
Pro 380 / Pro 381
360 / 380 / 380 AD
CA 4 / CA 8 / CA 12
FMI Q0 / Q1 / Q2 / Q3
Easy-Load
Standard
Quickload
MS/CA 4-12 / MS/CA 8-6
PTFE Tube 2 mm / 4 mm
SB 2V / 3V
WM 5

Deutsch

English

Français

ISMATEC®

IDEX
HEALTH & SCIENCE

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△ Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Caution: rotating pumphead creates a pinch and crush hazard

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.

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Safety precautions

- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tubebed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 9 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC®
from date of delivery: **2 years**

All other parts, excluding consumables: **1 year**

In case of any queries, please contact your ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

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Warranty terms

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- MCP *Standard* drive
Order No. ISM 404
- 1 power cord with an IEC 320
connector (female plug) and a
country specific mains plug
- Operating manual

If ordered:

- pump-head with
mounting plate
- other accessories

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➔ Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 Socket for footswitch and hand dispenser
- 2 Ventilator
- 3 Fuse-holder with voltage selector 115/230 V
- 4 Mains socket
- 5 RS232 In (female)
- 6 RS232 Out (male)
- 7 Valve connector
- 8 Analog interface
 - input for:
 - speed control
0 – 5 V or 0 – 10 V, and
0 – 20 mA or 4 – 20 mA
 - rotation direction
 - Run/Stop
 - Output for:
 - speed 0 – 10 V_{DC}
or 0 – 12 kHz

⚠ Mains voltage

Mains voltage	Voltage setting	Fuse rating
220 – 240 V _{AC}	230 V 50/60 Hz	2 x T1.25A/250V
110 – 120 V _{AC}	115 V 50/60 Hz	2 x T2.5A/250V

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

⚠ Socket/Power cord

Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

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Voltage setting 115/230 V Changing the fuses

➔ Switch the pump off,
pull out the mains plug.

- 1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).
- 2 Take out the voltage selector plate. Turn it and reinsert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.

- 3 Insert 2 new fuses
➔230 V: 2 x T1.25A/250V
➔115 V: 2 x T2.5A/250V

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

- 4 Shut the fuse-holder. The voltage rating is visible in the window.

Changing the fuses on the control board

The 2 fuses are fixed to the control board as illustrated opposite.

- a T4A/250V
- b F1.6A/250V

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

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
Starting the pump

- Mount the pump-head according to the mounting instruction manual supplied with the pump-head.
- Enter the ID Codes of the mounted pump-heads in the program memory currently used (see pages 16 and 50).
- Enter the tubing i.d. in each program used (see page 17).
- If using an FMI piston pump-head, adjust the angle for the piston stroke volume (see page 17).
- Insert the tubing into the peristaltic pump-head.
- Connect the pump tubing to the system.
- Connect the pump to the mains and switch it on with the power supply switch.

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TUBE I.D.	mm	PROGRAM
PUMP	rpm	Flow rate
DISP	Time	Volume
PAUSE	Time	TOTAL

Anzeige der Betriebsart
LEDs for operating modes
Affichage du mode d'exploitation

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Operating panel

- 1 Mains switch (on/off)
- 2 Digital LED display
- 3 Reduce value
- 4 Increase value
- 5 Run / Stop
- 6 Calibration / Data saving
- 7 MAX/MIN/RESET
(Multi-function, see page 15)
- 8 Rotation direction
- 9 Selecting the operating mode
➔MODE key

10 LEDs for active operating mode

- **TUBE I.D. mm**
Tubing i.d. selection
- **PROGRAM**
Program selection 1 – 4
- **PUMP rpm**
Pumping by speed
1.0 – 240.0 rpm
- **PUMP Flow rate**
Pumping by flow rate
- **DISP Time**
Dispensing by time
- **DISP Volume**
Dispensing by volume
- **PAUSE Time**
Pause for intermittent dispensing
- **TOTAL**
Read-out of totally delivered volume

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6 Pro4

7

TUBE I.D.	mm	PROGRAM
PUMP	rpm	Flow Rate
DISP	*Time	Volume
PAUSE	Time	TOTAL

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Start-up information

After switching on the power supply, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Baud rate
- 4 Identification Code (ID-Code) of the pump-head entered in the program last used (e.g. pump-head 380AD = 380A)
- 5 Roller back-steps
➔ are only displayed if at least 1 or more back-steps are entered
- 6 Program last used.
- 7 Displays the operating mode last used.

⚠ **Before starting the pump the first time or after mounting a new pump-head,** the ID-Code must be entered in each program used with this pump-head (see page 16 and 50).

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Resetting the parameters of the currently used program to the default values.

- 1 Power switch »OFF«
- 2 Push the **MODE key** while switching the drive on.

The following parameters are reset:

Mode:	PUMP rpm
Rotation direction:	clockwise
Speed:	100 rpm
Dispensing time:	4.5 sec.
Dispensing volume:	volume depending on pump-head
Pause time:	2.0 sec.

Number of dispensing cycles:	12
Roller back-steps:	0

➔ **Both pump-head and tubing i.d. entered in the currently set program remain stored.**

Resetting default values of all 4 programs at the same time

- 1 Power switch »OFF«
- 2 Push the **▼ key** while switching the drive on.

The following parameters are reset (in all 4 programs):

Pump-head:	Model MS-3
Tubing i.d.:	0.13 mm

➔ as well as all parameters mentioned above

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
Control keys

- a RUN/STOP**
Starts and stops the pump.
- b MODE**
Changes between operating modes (see page 12).
- c Rotation direction**
Changes the rotation direction.
- This function is blocked when using the diaphragm pump-head
 - Each time the rotation direction is changed, TOTAL is reset to »0«.
 - Functions only in the operating modes:
 - PUMP rpm
➔ minus sign = reverse direction
 - PUMP Flow rate
- d CAL**
Key for saving a set value.
- e MAX/MIN/RESET**
Multifunction key for:
- ➔ max. speed (when pump is running)
 - ➔ min. speed (when pump is idle) turns slowly, e.g. for mounting a pump-head
 - ➔ Reset button resets the calibration to the standard values (default), page 14
 - ➔ Reset button for accumulated volume in Mode TOTAL
- f Increment/Decrement keys**
- ▲ = increase value
 - ▼ = decrease value
- Maintaining pressure on ▲ or ▼ key changes display read-out into fast mode.

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Pump-head identification

In order to obtain correct pumping and dispensing values the individual identification code (ID-Code) of the mounted pump-head must be entered (**▲ for each of the 4 programs**), especially when a new pump is used for the first time and each time the pump-head is changed.

- 1 Power switch »OFF«
- 2 Keep CAL-key pressed and switch power »ON«
- 3 Press ▲ or ▼ key till the correct ID-Code of the mounted pump-head is set. For the ID-Code see the back of the pump-head or the table on page 50.
- 4 Enter the correct ID-Code with the CAL-key
(MODE display switches automatically to TUBE I.D.)
- 5 Enter the correct tubing i.d. with the ▲▼ keys (see also page 17).
- 6 Confirm with the CAL-key
- 7 Change to PROGRAM with the MODE key and select the next program with the CAL and ▲▼ keys (confirm with CAL).

➔ Repeat procedure 1 – 6 till the ID-Code(s) of the mounted pump-head(s) and the tubing i.d. are entered in all 4 programs.

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Entering the tubing i.d. or stroke angle

⚠ Before starting to pump it is essential to enter the following settings in order to obtain correct flow rates and dispensing volumes.

ONLY for Peristaltic pump-heads

- 1 Change mode to TUBE I.D. mm (for setting the tubing inner diameter)
 - 2 Press CAL-key (display blinks)
 - 3a Compare blinking figure with required tube i.d.
 - 4 Press ▲▼keys until correct i.d. is displayed
 - 5 Confirm entered value with CAL-key
- ➔ Select required operation mode with MODE-key.

ONLY for FMI piston pump-heads

Setting the stroke volume

In mode TUBE I.D. mm the piston angle can be entered as stated above from position 1 to 5:

- 3b set piston angle (between 1 and 10)
- ➔ The same angle must also be manually set on the pump-head.

ONLY for Diaphragm pump-heads

Tube i.d. is not relevant. In mode TUBE I.D. mm the display shows »0.0«.

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Program selection

When switching the pump on, it always selects the previously used program.

1 Change mode to PROGRAM by using the MODE-key.

2 Press CAL button, currently set program starts to blink.
Change to the required program with the ▲▼ buttons.

3 Confirm with the CAL-key.

4 The pump automatically returns to the last entered operating parameters of the selected program.

⚠ From now on, any changes carried out in the operating modes are automatically stored in the currently selected program.

TOTAL Volume

Change with MODE key to TOTAL. The totally delivered volume is displayed.

Three display readings are available:

Microlitre: e.g. 12.0 µl = **12.0µ**

Millilitre: e.g. 1200 ml = **1200**

Litre: e.g. 12.0 l = **12.0L**

(Display in litres from 9999 ml)

For zero-setting, either

- press the reset key twice, or
- switch the pump off.

➔ Each time the rotation direction is changed, TOTAL is reset to »0«.

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1 MODE

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	◦Volume
PAUSE	◦Time	◦TOTAL

2 CAL

Pro4

▲

Pro2

▼

3 CAL

4

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	◦Volume
PAUSE	◦Time	◦TOTAL

z.B. / e.g. / p.e. PUMP rpm

MODE


TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	◦Volume
PAUSE	◦Time	◦TOTAL

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

Pumping by drive speed


- 1 Change mode to PUMP rpm, 1.0 – 240.0 rpm, adjustable in steps of 0.1 rpm
- 2 • Enter required speed with ▲▼ keys
• Start pump with RUN/STOP
➔ The speed can also be adjusted while the pump is running.
- 3 Maintain pressure on MAX/MIN-key while the pump is running
➔ Pump is running at maximum speed (ideal for fast filling or emptying the system).
- 4 Maintain pressure on MAX/MIN-key when pump is idle (power switched on).
➔ Pump is turning at 10 rpm (for mounting the pump-head).

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Pumping by flow rate

- 1 Change mode to PUMP Flow rate
- 2 Enter the required flow rate with ▲▼ keys (is displayed in µl/min or ml/min)

Start pump with RUN/STOP

- ➔ The flow rate can also be adjusted while the pump is running.
- ➔ For an accurate flow rate we recommend to calibrate the pump (see page 21).
- ➔ During the pumping process the following values can be retrieved with the MODE-key:
 - Rotation speed (PUMP rpm)
 - Volume (TOTAL)

- 3 Maintain pressure on MAX/MIN-key while the pump is running
 - ➔ Pump is running at maximum speed (ideal for fast filling or emptying the system).

- 4 Maintain pressure on MAX/MIN-key when pump is idle (power switched on!)
 - ➔ Pump is turning at 10 rpm (for mounting the pump-head).

⚠ Entering the ID-Code of the currently mounted pump-head (pages 16 and 50) and the inner diameter of the tubing used allows the user to set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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Calibrating the flow rate

- 1 Change mode to PUMP Flow rate
Enter the required flow rate with the ▲▼ keys

- 2 Change mode to DISP Time
 - Enter 60 seconds by using the ▲▼ keys
 - Start pump with RUN/STOP

➔ Pump stops automatically after 60 seconds

➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:

- 3 Return to mode PUMP Flow rate
- 4 Press the CAL-button (displayed value blinks)
Enter the weighed or measured value with the ▲▼ keys
- 5 Confirm with the CAL-key (the flow rate setting returns automatically to the initially entered set point)
- 6 Start with the RUN/STOP-key

➔ The flow rate setting depends on the number of pump rollers and the tubing i.d. It is possible that the required set point cannot be entered accurately.
If necessary choose a tubing with a smaller i.d. or a pump-head with more rollers.

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

Please note

For optimum dispensing accuracy, we recommend that you keep the last drip as small as possible, e.g. by tapering the end of the tubing (e.g. with a tubing connector).

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Dispensing by time

The dispensing time can be entered from 0.1 s to 999 h.

- Change mode to DISP Time
- Enter the required dispensing time with the ▲▼ keys
 ➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds:	0.1" – 899.9"	(in 0.1 s steps)
Minutes:	15' – 899'	(in 1 min steps)
Hours:	15 h – 999 h	(in 1 h steps)

3 Start pump with RUN/STOP

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point only takes effect from the subsequent dispensing cycle.

➔ During the dispensing process the following values can be retrieved by pressing the MODE-key:

- speed (PUMP rpm)
- flow rate (PUMP Flow rate)
- totally pumped volume (TOTAL)

➔ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.

➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible. We recommend to use a tubing with a smaller i.d. and to increase the dispensing time.

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- Drehzahl (PUMP rpm)
- Fließrate (PUMP Flow rate)
- Volumen (TOTAL)

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Dispensing by volume

- 1 Change mode to DISP Volume
- 2 Use the ▲▼ keys for entering the required dispensing volume. Three display readings are available:
Microlitre: e.g. 12.5 µl = **12.5µ**
Millilitre: e.g. 230.5 ml = **230.5**
Litre: e.g. 12.5 Litre = **12.5L**
(Display in litres from 9999 ml)
➔ For dispensing accurately we recommend you calibrate the pump (see page 24).

- 3 Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the modes PUMP rpm and PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can also be changed even during the dispensing process. The newly entered set point appears shortly in the display taking effect, however, only with the subsequent dispensing step.

➔ During the dispensing process the following values can be retrieved with the MODE-key:

- Rotation speed (PUMP rpm)
- Flow rate (PUMP Flow rate)
- Volume (TOTAL)

➔ With the ▲▼ keys the rotation speed and flow rate can be changed even during the dispensing process.

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1 **MODE**

TUBE I.D.	◦mm	◦PROGRAM
PUMP	◦rpm	◦Flow Rate
DISP	◦Time	* Volume
PAUSE	◦Time	◦TOTAL

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
Calibrating the volume

- 1 Change **MODE** to **DISP Volume**
 - ➔ Use the ▲▼ keys for entering the required dispensing volume
- 2 Start with **RUN/STOP**
 - ➔ Pump stops automatically
 - ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes
 - ➔ Measure the dispensed liquid by volume or weight
- 3 Stay in the same mode (**DISP Volume**) and press the **CAL**-key (display blinks).
- 4 Enter the effectively dispensed volume with the ▲▼ keys.
- 5 Save setting with **CAL**-key.
 - ➔ Pump is calibrated and returns automatically to the range of the initially entered set point.
- 6 Start pump with **RUN/STOP**
 - ➔ Depending on the application conditions repeated calibration may be necessary.
 - ➔ The volume setting depends on the number of pump rollers and the tubing i.d. It is possible that the required set point cannot be entered accurately. If necessary choose a tubing with a smaller i.d. or a pump-head with more rollers.

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Dispensing by volume within a pre-set time

- 1 Change to mode DISP Time. Enter the required time with the ▲ keys (0.1 s – 999 h) (see page 22)
- 2 Change mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys.
- 3 Push the MODE-button.
➔ Both LEDs light up simultaneously.
- 4 Start pump with RUN/STOP
➔ The pump dispenses the required volume within the pre-set time.

a Changes / Corrections

- Press the Reset key
➔ For changing the dispensing time or dispensing volume, MODE returns directly to DISP Time
➔ Repeat steps 1 to 3

b Error message

- LEDs are blinking alternately
9999 = volume too large
time too short
1111 = volume too small
time too long

Despite volume and time adjustments, the LEDs keep on blinking:

- ➔ Change the tubing i.d. or the pump-head model in order to obtain another flow rate.

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

Please note

For entering the number of dispensing cycles see page 28.

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

1 Go to mode PAUSE Time
 • Enter the pause time (between 0.1 s and 999 h) with ▲▼ keys.

2 Press MODE key once
 ➔ Both LEDs for PAUSE and DISP Time light up simultaneously
 • Enter the dispensing time (0.1 s to 999 h) with ▲▼ keys.

3 Start pump with RUN/STOP key
 ➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.
 ➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process. The newly entered set point appears shortly on the display taking effect, however, only with the subsequent dispensing cycle.

➔ During the dispensing process the following values can be retrieved with the MODE-key:
 • drive speed (PUMP rpm)*
 • flow rate (PUMP Flow rate)*
 • totally dispensed volume (TOTAL)
 *The rotation speed and flow rate can be changed with the ▲▼ keys.

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
Please note

For entering the number of dispensing
cycles see page 28.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

- Change mode to PAUSE Time.
Enter pause time (between 0.1 s
and 999 h) with the ▲▼ keys.
 - Press MODE-key twice.
➔ Both LEDs for PAUSE and
DISP Volume light up simul-
taneously Enter the volume with
the ▲▼ keys.
 - Start pump with the
RUN/STOP key.
➔ Pump dispenses at intervals
until the RUN/STOP button is
pushed again.
➔ Active mode blinks.
- ➔ With the ▲▼ keys the dispensing
volume can be changed even du-
ring the dispensing process. The
newly entered set point appears
shortly on the display taking ef-
fect, however, only with the next
dispensing cycle.
- ➔ During the dispensing process the
following values can be retrieved
with the MODE-key:
- drive speed (PUMP rpm)*
 - flow rate (PUMP Flow rate)*
 - totally dispensed volume
(TOTAL)
- *The rotation speed and flow rate
can be changed with the ▲▼ keys.
- If necessary, calibrate the volume as
indicated on page 24.

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TUBE I.D.	◊mm	◊ PROGRAM
PUMP	◊rpm	◊ Flow Rate
DISP	◊Time	* Volume
PAUSE	*Time	◊ TOTAL

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

1 Change into one of the two modes for dispensing at intervals:

- DISP Time/Pause, or
- DISP Volume/Pause

2 Push the CAL-key

The number 0 for an infinite number of dispensing cycles blinks.

- ➔ With the RESET-key you can switch between 0 and the number of repetitions last entered.

3 With the ▲▼ keys the number of dispensing cycles can be adjusted.

4 Confirm with the CAL-key.

5 Start dispensing at intervals with the RUN/STOP-key.

- ➔ During the pauses the display shows the remaining number of dispensing cycles.

- ➔ During long pauses, the remaining interval time can be checked with the CAL-key.

- ➔ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key; when resuming the process, the pump continues with the remaining amount of dispensing cycles.

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a

TUBE I.D.	◦ mm	◦ PROGRAM
PUMP	◦ rpm	◦ Flow Rate
DISP	◦ Time	* Volume
PAUSE	◦ Time	◦ TOTAL

b

TUBE I.D.	◦ mm	◦ PROGRAM
PUMP	◦ rpm	◦ Flow Rate
DISP	* Time	* Volume
PAUSE	◦ Time	◦ TOTAL

c

TUBE I.D.	◦ mm	◦ PROGRAM
PUMP	◦ rpm	◦ Flow Rate
DISP	◦ Time	* Volume
PAUSE	* Time	◦ TOTAL

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Drip-free dispensing

With programmable roller back-steps or piston stroke back-steps (FMI pump-head) (1 – 100 steps or strokes)

- 1 Switch the pump off.
- 2 Maintain pressure on **◆** key and switch the pump on again.
- 3 Display blinks.
Enter the number of roller back-steps or reverse piston strokes with the **▲▼** keys.

- 4 Save input with the **◆** key.

Roller/Stroke back-steps are active in the modes:

- a Dispensing by volume (see page 23).
- b Dispensing a volume within a pre-set time (see page 25).
- c Intermittent dispensing by volume (see page 27).

⚠ Reject the first dispensing step. Afterwards, the roller/stroke back-steps are compensated.

Deleting roller/stroke back-steps

- Repeat step 1 and 2
- Enter figure 0 with the **▲▼** keys
- Save input by pressing the **◆** key

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

Pumping against pressure

The MCP *Standard* can be used for continuous duty at a differential pressure of max. 1.5 bar. For short-time operation higher differential pressures may be managed depending on the mounted pump-head and tubing i.d. (smaller i.d. is preferable).

Depending on the application conditions the FMI pump-heads allow operations at differential pressures of up to 7 bar.

In case of any doubts please contact your ISMATEC® agent.

When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

a Tubing used with cassettes do not need to be removed. Releasing the cassette on one side is sufficient.

b The single-channel pump-heads 360, 380 and Easy Load allow opening the tube-bed in order to release the tubing pressure.

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A digital display with a black background and white characters showing the text "OL 1".

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Overcurrent protector

The drive MCP *Standard* features an overload protector. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

OL I = Current Overload

Whenever this situation occurs, the pump must be switched off immediately. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure or tension lever tightened too much, etc.).

➔ Only when the operating conditions have been changed accordingly may the pump be started again.

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Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period.

If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used. Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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- Pin 2, local / remote
- Pin 3, stop / run
- Pin 4, direction
- Pin 13, speed IN

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs

Pin 2, local / remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, stop / run

- In normal operation (pin 2 open) pin 3 operates the foot-switch (see DIP switch settings, page 35)
- In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns clockwise; when connected to pin 1 (GND) it turns counter-clockwise.

Pin 5, speed IN

External speed control (0 – 5 V, 0 – 10 V, 0 – 20 mA, 4 – 20 mA). Input impedance and input range can be selected via a dip-switch inside the pump (see page 35).

Pin 7, +36V_{DC}

About +36 V_{DC} are available (max. current 1A). *)

*) The max. total current of the +36 V_{DC} supply on pin 7 and on pin 9 and 10 of the valve interface is 1 A. Fuse: F1.6A/250V, located on controller board, see diagram on page 36.

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Analog interface

Pin 9, motor speed

The default setting is 0 – 10 V_{DC}, proportionally to the motor speed 0 – 240.0 rpm.

Alternatively a frequency range from 0 – 12 kHz is available. Alternatives with sliding switch S2 inside the pump (see page 35).

Pin 10, +5 V_{DC}

About +5 V_{DC} are available. (max. current 0.5 A)

Pin 13, speed pin 5 internal

Depending on how the pump is operated, pin 13 has different functions:

Analog interface not activated

(Normal operation, i.e. pin 2 is open) Pin 13 serves as **auto-start function**.

If pin 13 is connected to pin 1 (GND), the pump can be started and stopped directly from the power supply (the power switch must be on).

Analog interface activated

(Pin 2 on GND)

- Pin 13 open:
The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND:
The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Please note

The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Settings of switch S1

Operation via foot-switch

With DIP switch 5 the user can switch between two possibilities:

- »FS toggle« (On/Off)
- »FS direct«

(On = as long as pressure is maintained on the foot-switch)

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5
Pin 3 (foot switch)	FS toggle						OFF*
	FS direct						ON
Pin 5 speed IN	0-5V	470kΩ		OFF*	OFF*	OFF*	
	0-10V	20kΩ		OFF	OFF	ON	
	0-20mA	240Ω		OFF	ON	OFF	
	4-20mA	240Ω		ON	ON	OFF	
Baud Rate RS232	9600 Bd		OFF*				
	1200 Bd		ON				

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Sliding switch S2

This switch affects Pin 9, motor speed

Position A: 0 – 10 V_{DC} (Standard)

Position B: 0 – 12 kHz

➔ **Make sure that the pump is disconnected from the mains supply.**

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➔Downloads

Valve interface

A detailed description of the valve interface is available for download on our web-site.

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Valve interface


Up to two valves can be operated directly from the VALVES interface with an adjustable effective voltage between 0 and ca. 26 V see page 37 'valve interface' command section. Commands X,Y,Z,& control the valves. Switching is done with an initial maximum voltage of 26 V_{DC} for 0.5 sec to make sure that the valve moves. The hold voltage can be adjusted with the 'Z _____' and 'X _____' commands through RS232 (but not at the pump's operating panel). A simple terminal program called 'IsmaTerm' is available for PC at http://www.ismatec.ch/int_e/downloads/downloads.htm.

The two 'sense' pins detect if a valve connector is plugged in. If just one valve is used between pin 1 (valve 2-) and pin 9 (valve 2+), the 'sense' pins 8 and 14 need to be linked by a piece of wire. If an additional second valve is present at pin 2 (valve 1-) and pin 10 (valve 1+), the 'sense' link should be a 47 ohm resistor instead.

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Valve interface

Pin 1, valve 2 –
PWM (000 .. 255) connection to GND for valve 2.

Pin 2, valve 1 –
PWM (000 .. 255) connection to GND for valve 1.

Pin 8, sense GND
GND for the detection of a valve plug via short-circuit with pin 14.

Pin 9, valve 2 +
About $+36 V_{DC}$ are available for valve 2 (max. current 1 A). *)

Pin 10, valve 1+
About $+36 V_{DC}$ are available for valve 1 (max. current 1 A). *)

Pin 14, sense
The valve plug must provide a short-circuit between pins 8 and 14. This enables the pump to detect the presence of a valve connector.

*) The max. total current of the $+36 V_{DC}$ supply on pin 7 and on pin 9 and 10 of the valve interface is 1 A. Fuse: F1.6A/250V, located on controller board, see diagram on page 36.

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⚠ Please note

When switching on the MCP *Standard* drive, the serial interface will only respond after a delay of 3 seconds.

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
Kommunikationsparameter Communication parameters Paramètres de communication

Standardeinstellung
Default setting
Réglage standard

**9600 baud, 8 bit,
1 stopbit, no parity**
(DIP Switch 1 off)

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1 – 8). It is factory-set at 1 and can be changed with @.

The address will be stored permanently (even after the pump has been switched off).

The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13).

The pump confirms most of the commands with an asterisk *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Each pump requires an individual address
(see page 38).

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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
Accessories

Software Labworldsoft®
is an innovative Windows software for controlling, measuring and evaluating operations in connection with a number of different laboratory instruments.

Hardware requirements
486 PC, at least 33 MHz and 8 MB RAM, mouse, VGA display, monochrome, with at least 16 grey scale values

Software requirement
Windows 3.1 or higher, Windows 95

Order No. SOF103

Interface card Moxa C168
including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments
Order No. ISA IM0015C168
Order No. PCI IM0030

Adaptor 25 F / 9M
25 pins, female / 9 pins, male for IM0015C168
(1 adapter for each controlled instrument)
Order No. XC0009

Extension cable 9F / 9M
9 pins, female / 9 pins, male for XC0009 (2 m long)
Order No. AG0013

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MS/CA Click'n'go IS 3510
 MS/CA Anpresshebel ¹⁾ IS 0649

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
CA Anpresshebel ¹⁾ IS 3820
 Adapter PVDF IS 3861
 (pro Kassette 2 Stk. bestellen)

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Accessories

Foot switch IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

Spare-cassettes in POM-C

CA Click'n'go IS 3710
 CA pressure lever ¹⁾ IS 0122
 Adaptor POM-C IS 0123
 (order 2 adaptors per cassette)

MS/CA Click'n'go IS 3510
 MS/CA pressure lever ¹⁾ IS 0649

Spare-cassettes in PVDF

CA pressure lever ¹⁾ IS 3820
 Adaptor PVDF IS 3861
 (order 2 adaptors per cassette)

MS/CA pressure lever ¹⁾ IS 3629

¹⁾ The cassettes with pressure lever are available on request. This type of cassette provides better results at elevated differential pressure conditions (>1 bar).

⚠ Click'n'go-Cassettes

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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MS/CA Click'n'go IS 3510
 MS/CA levier de pression ¹⁾ IS 0649

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


Standard-Schläuche
Standard tubing / Tubes au mètre

Pump Tubing

Quality tubing for tubing pumps
in laboratories and production

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Please refer to our detailed tubing
documentation or our website

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Accessories

Materials POM-C

- Polyoxymethylene-Copolymer
- good chemical resistance to many organic solvents and strong alkaline chemicals
 - is affected by strong acids and oxidizing substances
 - UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

Materials PVDF

- Polyvinylidene fluoride
- very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
 - not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
 - stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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Mounting

Each pump-head is supplied with an instruction manual.

Do not forget to enter the correct ID code (see page 16 and 50).

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
0.09 – 530 ml/min
SB 3V
ISM 734/011

Interchangeable Pump-heads

These pump-heads can be ordered separately and mounted on the MCP *Standard* drive.

➔ For the piston pump-heads the suffix „-W“ resp. „-WT“ means: with isolation (W) or temperature gland (WT) for crystallizing media.

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
⚠ Must be entered for each of the 4 programs using the particular pump-head.

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△ Please note

In the following tables the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®

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Maintenance

Provided the MCP *Standard* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective MCP *Standard* drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., and model.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our **Warranty Terms as well as our Terms and Conditions of Sale**.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Drive

Motor type DC motor
Speed range 1.0 – 240.0 rpm
digitally adjustable
in steps of 0.1 rpm

Differential pressure

1.0 – 1.5 bar (14.5 – 21.8 psi)
(depending on pump-head) for
short-time use and with small
tubing i.d.'s higher pressures are
possible.

Remote control

via RS232 and analog interface (see
page 33 – 43)

Mains connection

230 V_{AC} 50 Hz
115 V_{AC} 60 Hz

Fuse rating

2 x T1.25mA/250V(230 V) or
2 x T2.5mA/250V (115 V)

Fuse rating on control board

T4A/250V
F1.6A/250V

Protection rating

IP 30

Power consumption

Max. 100 W

Operating conditions

Temperature 41 to 104°F
Rel. humidity max. 80%
- not condensing, at normal laboratory conditions

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conformément à:
EN 61326-1, EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH