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CE 03.30.10 SP

# Operating Manual

Mode d'emploi

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Deutsch

**Tubing pump**  
with fix speed drive

**MS-CA**  
**Stand-mounted pump**  
with 2 or 4 channels  
6 or 8 rollers

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**Contents**

Safety precautions	4
Warranty terms	7
Product	8
Package contents	9
Rear panel	10
Mains voltage	10
Changing the fuses	11
Operating panel	
Starting the pump	12
Inserting the tubing	14
Pumping against back-pressure	15
Running in period for tubing	16
Tubing live	16

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**Contents**

Flow rates	17
Accessories	18
Maintenance	20
Repairs	20
Spare-parts	21
Disposal	21
Technical specifications	22

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## Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump **must not** be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes

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Caution: rotating pumphead creates a pinch and crush hazard

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## Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of pump tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pumphead). Therefore, the pump must not be operated before the cassettes are fully snapped-in.
- Do not manipulate the pumphead before the pump is switched off and disconnected from the mains.

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## Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 11 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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## Warranty

### 2 years from date of delivery

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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## Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC® the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

### Our warranty becomes invalid in case of:

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third-party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing

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## Product

The tubing pump is with 2 or 4 tube channels and 6 or 8 pump rollers available.

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## Product

### Package contents

- MS-CA pump (as ordered)  
including Click'n'go Cassettes
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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

### Socket/Power cord

Use exclusively the originally supplied power cord.

The socket must be earthed (protective conductor contact).

### Before starting-up

Check if the voltage setting rating on the type plate complies with your local mains voltage.

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
## Changing the fuses

➔ Switch the pump off,  
pull out the mains plug.

1 Pull out the fuse-holder by opening it  
with a small screw-driver (size 0).

2 Insert 2 new fuses  
230 V<sub>AC</sub>: 2 x T100mA/250V  
115 V<sub>AC</sub>: 2 x T200mA/250V

⚠ Use always 2 slow-blow fuses of  
the same type complying with the local  
mains voltage.

3 Shut the fuse-holder.

### ⚠ Before starting-up

Check if the voltage setting of the fuses  
complies with your local mains voltage.

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
## Operating panel

- 1 Mains socket
- 2 Fuse holder
- 3 Switch for Run/Stop and flow direction  
 Position A = clockwise direction  
 Position B = stop  
 Position C = counter-clockwise  
 direction

## Starting the pump

- 1 Check the voltage rating on the type plate (on rear panel). If necessary, change the rating as indicated on page 10–11.
- 2 Connect the pump tubing to the system.
- 3 Connect the pump to the mains.

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For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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## Starting the pump




4. Depending on the required flow direction, press the upper or lower part of the Run/Stop switch.  
If the pump is correctly connected to the power supply, the pump-head now starts to rotate in the required direction.
5. For stopping the pump, put the rocker switch back to its neutral middle position „B“.

 When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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Ersatz-Kassetten MS/CA Click'n'go  
Bestell-Nr. IS 3510

## Inserting the tubing

- 1 ➔ Switch the pump off
- 2 Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously
- 3 Insert the 3-stop color-coded tubing with one collar into the cassette
- 4 Let the tubing hang down (prevent it from being twisted)
- 5 Insert the tubing with the second collar at the other end of the cassette
- 6 Reinsert the cassette into the roller-head
- 7 ➔ Switch the pump on

➔ When the pump is idle,  
release all cassettes at the fixing-tongue  
(2).

Spare-cassettes MS/CA Click'n'go  
Order-No. IS 3510

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Cassettes suppl. MS/CA Click'n'go  
No de commande IS 3510

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### **Pumping against pressure**

The MS-CA pump can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures)

### **When the pump is idle**

When the pump is idle, we recommend you release the tubing from pressure. Releasing the cassette on the right side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

#### **⚠ Syphoning effect**

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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## Please note

Please refer to our detailed tubing documentation.

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## Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

## Tubing life

The service-life of the tubing depends on the application and the tubing material used.  
Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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## Flow rates

In the table listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® ST tubing, without differential pressure

We will be pleased to help the user select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®

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MS/CA Click'n'go Cassette  
MS/CA Click'n'go cassette

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MS/CA pressure lever\*

Order No. IS 3629

\* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

<sup>1)</sup> Materials POM-C / PVDF see page 19

### ⚠ Click'n'go-Cassettes

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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### Materials POM-C

Polyoxymethylene-Copolymer

- good chemical resistance to many organic solvents and strong alkaline chemicals
- is affected by strong acids and oxidizing substances
- UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

### Materials PVDF

Polyvinylidene fluoride

- very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
- not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
- stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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

#### Pump Tubing

Quality Tubing for Tubing Pumps  
in Laboratories and Production

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### Pump Tubing

3-stop tubing  
Please refer to our detailed tubing  
documentation or our website:  
[www.idex-hs.com](http://www.idex-hs.com)

### Maintenance

Provided the tubing pump MS-CA is  
operated properly and in compliance with  
this manual, the tubing is the only part  
that is subject to wear and tear.

### Repairs

For repairs please send the defektive  
MS-CA to your ISMATEC®  
representative.  
Please give information on defect, date of  
purchase, serial-no., and model.

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## Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask your ISMATEC® distributor for:

- parts lists
- replacement parts for pump drive
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

## Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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### Speed ranges:

2- and 4 channels	20 rpm
2- and 4-channels	40 rpm
only 2-channels	60 rpm

### Differential pressure

max. 1.0 bar (14.5 psi)  
Depending on tubing material.  
Tubing with small i.d. and/or pressure lever cassettes may enable higher pressures.

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## Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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## Technical Specifications

### Operating conditions

Temperature	+5 to +40°C,
Rel. humidity	max. 80%
– not condensing, at normal environmental conditions	

### Dimensions/Weight

	(DxWxH)
2 channels	4.9x3.5x5.3 inch
4 channels	5.7x3.5x5.3 inch

	Weight
2 channels	2.65 lb
4 channels	2.87 lb

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conformément à:  
EN 61326-1 / EN 61010-1

## **IDEX Health & Science**

Home of Upchurch Scientific &  
Ismatec Products  
619 Oak Street  
Oak Harbor, WA 98277  
United States  
Tel: +866 339 4653  
Fax: +360 679 3830  
fluidics.americas@idexcorp.com

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Ask for the latest  
ISMATEC-catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH