
## Mini-S

### Mini-S Pumps

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**620/640/660**

3 Kanäle  
6 Rollen

**820/840/880**

3 Kanäle  
8 Rollen

Deutsch

### Mini-S Pumps

Stand-mounted  
fixed-speed pumps

**620/640/660**

3 Channels  
6 Rollers

**820/840/880**

3 Channels  
8 Rollers

English

**ISMATEC®**

**IDEX**  
HEALTH & SCIENCE

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**Contents**

Safety Symbols	3
Safety precautions	4
Warranty terms	5
Product	6
Short description	6
Motor	7
Pump-head	7
Spring-loaded tube-bed	7
Mains voltage	8
Inserting the tubing	9
Starting your pump	10
Running the tubing in	11
Service life of tubing	11
Tubing quality	11
Type of tubing	12
Pumping against back-pressure	12
Pump tubing	12
When the pump is not in use	13
Tube connectors	13
Adhesives for tubing connections	13
Tubing life	13
Flow-rate chart	14
Maintenance	15
Changing a defective fuse	15
Spare parts	16
Disposal	16
Standards and Specifications	16
Technical Specifications	17

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## Safety Symbols

Below you find the following symbols used throughout this manual:

This "STOP sign" denotes a hazard. The user is instructed to refrain from a certain practice or procedure. Failure to comply with such warnings can result in personal injury and/or damage to the product.

This "caution sign" denotes a hazard. Certain practices or procedures may only be carried out by taking the necessary precaution.

**Caution:** rotating pumphead creates a pinch and crush hazard.

General information about something important.

➔ Wherever you see these signs, you are requested to follow the specific warnings and recommendations given.

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## ⚠ Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

-  **Caution:** rotating pumphead creates a pinch and crush hazard.
- The circuit between mains supply and pump has to be earthed (grounded).
- The pump must not be operated outside the designated operating and environmental conditions.
- Ensure that the pump is switched off and the power cable disconnected before changing pump tubing.
- The pump must not be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of rotating parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend placing the pump in a collecting basin.
- The permeability of pump-tubing depends on the material used and pressure conditions. Moreover, tubing can also become electrostatically charged.

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## ⚠ Safety precautions

- Only new fuses, according to the specifications stated on page 15 in this manual, may be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies, ISMATEC® denies any responsibility.

## Warranty Terms

### For all ISMATEC® panel-mounted tubing pumps excluding consumables:

**from date of delivery: 2 years**

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

### Our warranty becomes invalid in case of:

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third-party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications.

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## Product

➔ Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC<sup>®</sup> representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

Check if the following items are included in the delivery:

- Mini-S pump
- Power cord
- Operating instructions

## Short description

The **Mini-S** is a compact, stand-mounted peristaltic pump designed to easily meet the requirements of numerous pumping applications at low flow-rates.

This **3-channel fixed-speed pump system** allows you to deliver liquids or gases in the clockwise and counter-clockwise direction.

The following **6 models** are available:

All six pump models have the following features:

**Fixed-speed drive**  
**Spring-loaded tube-bed**

### Intended area of application:

The Mini-S is a standard tubing pump designed for suctioning, feeding and circulating liquids and gases in chemical and pharmaceutical laboratories and similar areas of applications. Depending on the tubing material chosen, it is suitable for almost all liquids (including liquids which contain particulates), corrosives and abrasive substances and gases. When using sterilisable tubing, the pump can also be used for pumping under sterile conditions.

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**This instrument must not be used for medical applications on human beings.**

#### **Motor**

The Mini-S pump has a synchronous motor with an excellent start/stop behaviour. The power consumption is 11 Watts.

#### **Pump-head**

This pump is equipped with a 3-channel pump-head with 6 or 8 self-lubricating pump rollers.

#### **Spring-loaded tube-bed**

The Mini-S pump series features a 3-channel tube-bed with a spring-loaded pressure lever. The pressure on the tubing can easily be adjusted with the knurled-head screw (see figure below).

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➡ The Mini-S pumps are equipped with a 12 mm rod and must be mounted to a stand (see picture above).

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
## Mains voltage

The Mini-S pump models are available in two voltage versions. Depending on your order it is supplied for a voltage of

**230 V (220-240 VAC) 50 Hz**  
**or**  
**115 V (110-120 VAC) 60 Hz**

### **⚠ Before operating your pump:**

Please check whether the voltage indicated on the identification plate (on the rear panel) agrees with your local mains voltage available.

**⚠ The mains socket used has to be earthed (grounded) (protective conductor contact),** as otherwise a safe operation is not ensured, or the instrument may be damaged. (For the power connection, we recommend you use the original power cord).

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

## Inserting the tubing

➔ The Mini-S pumps take 3-stop collared tubing, allowing the same tubing to be used twice, and hence extending service life by 100%.

For loading or changing tubing you proceed as follows:

1.  Switch off the pump and unplug it.
2. Open the tube-bed by turning the white lever (5) to the right or the left. The spring-loaded 3-channel tube-bed (4) is lowered and provides free access.
3. Insert one collar of the tubing into the tube-holding slit of the tube-bed. Feed the tubing over to the other end of the tube-bed and insert the second collar into that slit (see figure below).  
➔
  - Insert the first tubing in the back of the tube-bed.
  - The tube must lie correctly in the tube-bed (collars pointed upwards). It must not be twisted.
4. Do the same with the other 2 pump tubes.
5. Close the pump-head by turning the white lever (5) into the vertical position.
6. Start the pump as described on page 10.
7. Adjust the correct pressure to the tubing by turning the knurled-head screw (6) so that a smooth flow is achieved.
  - ➔ It is most important to adjust the tube-pressure under operating conditions, i.e. when all the tubes are connected and the pump system is completely installed. Apart from the tube size and the material, the correct pressure also depends on other parameters such as back-pressure, viscosity, etc.
  -  The knurled-head screw (6) must not be tightened more than necessary. Otherwise, tubing and pump are exposed to unnecessary strain, which reduces the service life considerably.
  - ➔ Depending on the application and tubing material, re-adjustment of the tube pressure is necessary from time to time.

## Offenes Schlauchbett Open tube-bed


Color Code Stopper  
Color-coded collar

Color Code Stopper  
Color-coded collar

Schlauchbett (4)  
Tube-bed (4)

Anpreßhebel (5)  
Lever (5)

## Geschlossenes Schlauchbett Closed tube-bed

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
### Starting your pump

⚠ When the pump is running, the pump-head must be closed. A revolving, open pump-head is dangerous as parts of your body (e.g. fingers, hair) or other objects may become trapped, which can cause serious personal injury or damage to the instrument.

1. Please ensure that the pressure lever (5) is locked and the tubes lie correctly in the tube-beds. Check if all tubing connections are tight.
2. ⚠ Check whether the location of the pump as well as the connection and safety measures agree with the "Warranty Terms" and "Safety Precautions" stated on pages 4 to 5.
3. ⚠ Connect the pump (socket 1) to the mains supply by using a protective earth cable. (We recommend you use the original power cord).
4. Press on the top or bottom of switch (3), depending on the flow-direction you require (see drawing below). If the pump is connected to the mains, the pump-head now starts to turn in the required direction.
5. To stop the pump, press switch (3) again so that it is back in the middle position (see drawing below).

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### Running the tubing in

➔ Every new tubing needs a certain time to run itself in. If **constant and reproducible flow-rates** are required, we recommend you **run new tubing in for at least 1 hour using water or the medium to be delivered.**

### Service life of tubing

The service life of tubing depends mainly on the application (chemicals being used), tubing material, pressure, temperature, viscosity, etc. (For more information about the service life, please refer to the current ISMATEC® catalog.

### Tubing quality

➔ The pump-heads of ISMATEC® peristaltic pumps are designed for ISMATEC® tubing. The use of ISMATEC® tubing guarantees reproducible results and long service life.

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## Type of tubing

These cassette pumps take tubing with 3 color-coded collars. Thanks to these collars the tubes are always inserted under the same conditions, which results in a high repeatability of the flow rate. What is more, the different color codes help to easily recognize the various tubing diameters.

Tubing with 3 collars can be used twice. Thanks to this advantage the service life of the tubes can be doubled.

**The selection of the right tubing material, especially with regard to the chemical compatibility, is entirely the customer's responsibility. For the chemical compatibility, please refer to the guidance chart in our current pump tubing selection guide.**

The pump-heads are designed for ISMATEC<sup>®</sup> tubing. The use of ISMATEC<sup>®</sup> tubing guarantees a long tubing life and reproducible results. For further information about pump tubing, please refer to the current pump tubing selection guide.

## Pumping against back-pressure

The Mini-S pumps **must not be used for applications with a differential pressure (back-pressure) exceeding 1.5 bar**. (When three tubes with the three largest inner diameters (**2.54, 2.79 and 3.17 mm**) are used simultaneously, the back-pressure **must not exceed 0.5 bar**).

## Pump tubing

These cassette pumps take tubing with 3 color-coded collars. Thanks to these collars the tubes are always inserted under the same conditions, which results in a high repeatability of the flow rate. What is more, the different color codes help to easily recognize the various tubing diameters.

Tubing with 3 collars can be used twice. Thanks to this advantage the service life of the tubes can be doubled.

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### When the pump is not in use

→ When the pump is idle, we recommend you release the tubing from pressure. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

#### ⚠ Syphoning effect

When the tubing is un-occluded, the fluid can flow back to the reservoir.

### Tube Connectors

A number of different tube connectors are available. Please refer to the latest ISMATEC® catalog.

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
#### Adhesives for connections with tube sleeves

Schlauchmaterial Tubing material	Klebstoff Adhesive	Inhalt Content	Bestell-Nr. Order-No.
Silikon Silicone	Silikon-Kautschuk-Kleber Silicone rubber adhesive	100 g	SC 1002
C-Flex C-Flex	Cyanacrylat Cyanacrylate	20 g	SC 1005
Tygon® R3603 Tygon® F-4040-A	Cyclohexanon oder/or Loctite-414	10 ml 20 g	SC 1000 SC 1001
Fluran F-5500 (Viton®)	Loctite-414	20 g	SC 1001
Ismaprene/PharMed®	Sicomet-750	20 g	SC 1003

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3-Stopper-Schläuche  
3-stop collared tubing

### Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

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The Mini-S pump models can be used with 26 different tube sizes which provide the following flow-rates per pump-channel:

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Please refer to our detailed tubing documentation or our website

**The selection of the right tubing material, especially with regard to the chemical compatibility, is entirely the customer's responsibility.**

**For chemical compatibility, please refer to the guidance chart in our current pump tubing selection guide.**

### △ Please note

The above mentioned flow-rates are approximate values, determined with water, at ambient temperature, without back-pressure and with Tygon tubing.

**Exact flow-rate values depend on various parameters such as quality and age of the tubing, back-pressure, temperature and viscosity of the liquid, etc.**

We will be pleased to help the user select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results, we recommend you use tubing from ISMATEC.

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
## Maintenance

Your ISMATEC® peristaltic pump does not require any special maintenance. Provided that you operate it properly and in accordance with the guidelines stated in this instruction manual, the tubing is the only part that needs replacing from time to time depending on the application and the operating conditions.

## Changing a defective fuse

The fuse-holder is included in the mains supply socket (see drawing below).

1. ⚠ **Make sure that the pump is disconnected from the mains.**
2. Press the two clips (3) on both sides against the fuse-holder (2) and, at the same time, pull the fuse-holder out.

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3. You will find two surge-proof fuses in this fuse-holder.

- T100mA/230 V  
**or**
- T200mA/115 V

Change **both** fuses.

⚠ **Ensure that the new fuses are also surge-proof and the rating is correct as given above.**

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## Spare-parts

After the warranty period, your ISMATEC<sup>®</sup> representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

➔ **In case of any queries or if anything is unclear regarding the handling of your ISMATEC<sup>®</sup> pump, please contact your local ISMATEC<sup>®</sup> representative.**

## Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

## Standards and Specifications

The ISMATEC<sup>®</sup> pump supplied with this instruction manual complies with the CE-specifications and has the **CE-mark**. Additionally, this instrument meets the following standards:

<b>Protection rating:</b>	<b>IP 30</b>
<b>Electromagnetic emission:</b>	<b>EN 50081-1 EN 55011, Class B</b>
<b>Electromagnetic immunity:</b>	<b>EN 50082-1</b>
<b>Safety standards:</b>	<b>EN 61010-1 IEC 61010-1</b>

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## Technical Specifications

### Motor type

Synchronous motor with fixed-speed

### Operating conditions

Temperature +4 up to +40°C

Rel. humidity max. 80 %

- not condensing, at normal environmental conditions

### Differential pressure

1.5 bar for tubing with an i.d. up to 2.29 mm

0.5 bar for tubing with an i.d. 2.54, 2.79 and 3.17 mm, when all three tube-channels are used simultaneously

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**Ismatec Laboratoriumstechnik GmbH**