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# Operating Manual

# Mode d'emploi

## Ventillose Taumelkolbenpumpe

### **REGLO-CPF Analog** Antrieb ISM 1014

#### **Kopf 2–25 µl/Hub**

FMI 009    Keramikkolben  
FMI 010    Stahlkolben

#### **Kopf 5–50 µl/Hub**

FMI 005    Keramikkolben  
FMI 013    LF = Low Flow

#### **Kopf 10–100 µl/Hub**

FMI 007    Keramikkolben

Deutsch

## Valveless Rotary Piston Pump

### **REGLO-CPF Analog** Drive ISM 1014

#### **Head 2–25 µl/stroke**

FMI 009    ceramic piston  
FMI 010    steel piston

#### **Head 5–50 µl/stroke**

FMI 005    ceramic piston  
FMI 013    LF = Low Flow

#### **Head 10–100 µl/stroke**

FMI 007    ceramic piston

English

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## Piston pump- heads require thorough flushing after use

We recommend you flush the pump-head **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

**Do not pump fluids containing suspended solids that are abrasive.** Please observe our note regarding filtration in the chapters Safety Precautions and Maintenance.

**Piston fractures are not covered by our warranty.**

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**Contents**

Safety precautions	6
Warranty terms	10
Product	11
Rear panel	13
Mains voltage	13
Voltage setting	14
Changing the fuses	14
Starting the pump	15
Zero point calibration	16
Setting the stroke volume	17
Overload protector	18
Pumping against pressure	19
When the pump is not in use	19

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**Contents**

Analog interface	20
Settings switch S1	23
Accessories	24
Flow-rates	25
Maintenance	26
Service and repairs	27
Disposal	27
Replacement parts	28
Technical specifications	29
<b>Appendix</b>	31
Mounting a ceramic piston pump-head on the REGLO-CPF Drive	

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## Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes.
- The circuit between mains supply and pump has to be earthed.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
  
- The pump must be switched off before taking off the pump-head.

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## Safety precautions

- Only new fuses, according to the specifications stated on page 14 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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## Safety precautions

### Media containing particles:

We recommend you rinse the pump-head **after every pumping process** to prevent the piston from seizing by remaining particles (especially after pumping media containing salt, protein or particles).

**Do not pump media containing particles that exceed a diameter of 0.8 mm.**

**Piston stroke or stroke volume must be set large enough with respect to the particle size!**

**Piston fractures are not covered by our warranty.**

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## Safety precautions

### **Corrosive liquids:**

Corrosive liquids may eventually produce leak paths around the sealing surfaces of the pump-head.

The pump should be inspected for leaks on a regular basis.

### **Dry running:**

Before starting to run the pump we advise you fill the pump-head with liquid.

Ensure that for priming the pump does not run dry for more than 15 sec.

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## Warranty

**For all parts manufactured by ISMATEC®**

**from date of delivery: 2 years**

**All other parts, excluding consumables: 1 year**

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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## Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

### Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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- FMI 007 Keramikkolben

## Product

The REGLO-CPF *Analog* is a valveless rotary piston pump.

### Package:

- Pump drive ISM 1014
- Power cord
- Operating instructions
- Pump-head as ordered

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## Product


Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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## Rear panel

- 1 Analog interface
- 2 Fuse-holder with voltage selector 115/230V
- 3 Mains socket

### Mains voltage

Mains voltage	Fuse rating
220-240 V <sub>AC</sub>	2 x T800mA/250V
110-120 V <sub>AC</sub>	2 x T800mA/250V

### Socket/Power cord

Use exclusively the originally supplied power cord.  
The socket must be earthed (protective conductor contact).

### Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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## Voltage setting 115V / 230V and changing the fuses

➔ Switch the pump off,  
pull out the mains plug.

1 Pull out the fuse-holder by opening  
it with a small screw-driver (size 0).

2 Take out the voltage selector plate.  
Turn it and re-insert it into the fuse-  
holder so that the required voltage  
rating is facing the window of the fuse-  
holder.

3 Insert 2 new fuses

230 V<sub>AC</sub>: 2 x T800mA/250V

115 V<sub>AC</sub>: 2 x T800mA/250V

⚠ Use always 2 slow-blow fuses of  
the same type complying with the local  
mains voltage.

4 Shut the fuse-holder. The voltage  
rating is visible in the window.

### ⚠ Before starting-up

Check if the voltage setting visible in the window  
of the fuse-holder complies with your local mains  
voltage.

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
## Starting the pump

- 1 Set the required speed on the 2-digit speed selector.  
99 = 1800 rpm  
(max. revolution)  
  
➔ The rotation speed can be adjusted while the pump is running.
- 2 Choose the rotation direction.
- 3 Start the pump with the main switch = run/stop function.
- 4 Overload indicator (red LED)

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## Zero point calibration

Before using the pump the first time, or whenever necessary, we recommend you calibrate the 0-point of the flow control ring.


- Start the pump and pump the liquid to be delivered.
- While the pump is in operation, turn the knurled adjustment ring to the right until the pump stops delivering the liquid (the motor is still running!).
- Loosen the fixing screw of the revolving flow control ring (black) and turn it to the value ,0'.
- Tighten the fixing screw again.

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## Setting the stroke volume

For setting the stroke volume, turn the knurled adjustment ring till the desired value (0...450) on the fixed and revolving flow control ring appears. (The fixing screw must not be loosened!)

Basically the following applies:

- **Small stroke volume/high speed:**
  - low pulsation
  - better resolution of the flow rate and dispensing volume
  - tendency to cavitation (insufficient flow-through of the liquid in the pump-head).
- **Large stroke volume/low speed:**
  - for elevated viscosity
  - for differential pressures
  - reduced resolution of the flow rate.

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## Overcurrent protector

The pump REGLO-CPF *Analog* features an overload protector. When an overload condition occurs, an activated overcurrent protector is indicated by the red LED and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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### **Pumping against pressure**

The REGLO-CPF *Analog* can be used for continuous duty at a differential pressure of max. 6.9 bar (100 psi).

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
- Pin 2, remote
- Pin 3, start
- Pin 4, direction
- Pin 13, internal speed

**Analog-Eingang**  
**Analog input/ Entrée analogique**  
 Pin 5, speed IN  
 0–5 V<sub>DC</sub> / 0–10 V<sub>DC</sub> / 0–20 mA / 4–20 mA

**Analog-Ausgang (TTL-Pegel)**  
**Analog output (TTL level)**  
**Sortie analogique (niveau TTL)**  
 Pin 9, speed OUT 0–9 kHz

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## Analog interface

**Pin 1, GND (ground)**  
 Reference potential for all other inputs.

**Pin 2, remote**  
 For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

**Pin 3, start**  
 In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND).

**Pin 4, direction**  
 In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

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
External speed control

(0–5V, 0–10V, 0–20mA, 4–20mA)

Input impedance and input range can be selected via a dip-switch inside the pump (see p. 24).

### Pin 7, +22V<sub>DC</sub>

About +22 V<sub>DC</sub> are available (max. current 0.2 A).

### Pin 9, speed OUT

Frequency proportional to the rotation speed: 0-9 kHz

(with each revolution of the piston, 300 impulses are generated).

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
## Analog interface

**Pin 10, +5V<sub>DC</sub>**  
 About +5 V<sub>DC</sub> are available.  
 (max. current 0.1 A)

**Pin 13, internal speed**  
 Analog interface activated  
 (Pin 2 on GND)

- Pin 13 open: The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
## Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

**⚠ The instrument should only be opened by a qualified technician!**  
Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	18k $\Omega$	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V	38k $\Omega$	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250 $\Omega$	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250 $\Omega$	OFF	OFF	OFF	ON	ON	OFF

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### Foot switch

Order No. ISM 891

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filing tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

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## Please note

In order to reach the max. flow rate, the stroke volume must be set to the max. value (25, 50 or 100 µl) depending on the pump-head model.

If, depending on the application, this stroke volume must be reduced we recommend you not to go below 10% of the maximum value.

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\*Range is outside the recommended min. stroke volume

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## Maintenance

We recommend the user to flush the pump-head thoroughly after pumping in order to prevent the piston from seizing by remaining particles (especially when media has been pumped that contained salt, proteins or particles).

Otherwise, the REGLO-CPF *Analog* requires no particular maintenance or mechanical adjustments.

Avoid pumping media containing abrasive solids.

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## Please note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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## Service and repairs

Like all pumps, the REGLO-CPF *Analog* contains components that will wear over a period of time.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

## Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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## Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask from your ISMATEC® distributor for:

### ■ Piston-/cylinder group

Type	Order No
H-00CKC-LF	FMI 311
H-00SKY-LF	FMI 310
H-0CKC	FMI 307
H-1CKC	FMI 308

### ■ Lip seals (2 pieces are needed)

Type	Order No
H408-00J (RH00)	FMI 325
H408-0A (RH0)	FMI 312
H408-1A (RH1)	FMI 319

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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### Speed

18–1800 rpm, digitally adjustable in steps of 1%.

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## Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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IP 30

### Operating conditions

Temperature +5°C to +40°C  
(+41°F to 104°F)  
Rel. humidity max. 80%  
not condensing, at normal laboratory conditions

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
## Mounting a ceramic piston pump-head on the REGLO-CPF drive

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### Assembly parts from ISMATEC

- 1 REGLO-CPF drive
- 2 Coupling completely assembled
- 3 Mounting plate
- 4 Countersunk head-screws (4 pcs.)
- 5 Allen key 2.5 mm

### from FMI

- a) CPF pump-head (RH00, RH0, RH1)
- b) Screws (2 pieces)
- c) Serrated lock washers (2 pieces)

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


## Mounting a ceramic piston pump- head on the REGLO-CPF drive

- 1 Attaching the mounting plate**  
Remove screws (b) and serrated lock washers (c) from the pump-head.

- 2 Fasten the mounting-plate (3) on the pump-head with the 2 screws (b) and serrated lock washers (c).**

### IMPORTANT

The 4 countersinks (5) must be concave to the pump-head.

- 3 Fixing the coupling**  
Place the completely assembled coupling (2) onto the pump-head shaft and push it on as far as it will go.  
  
The Allen screw (d) must be pointing towards the flat part of the pump-head shaft (e).

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
6

### Mounting a ceramic piston pump-head on the REGLO-CPF drive

- 4 Tighten the coupling (2) with the Allen key size 2.5 mm (5).

### 5 Mounting the pump-head

#### Important

- The scale of the red flow control ring (f) must be pointing upwards.
  - Hold the mounting plate in such a way, that the screws (b) are on the top and at the bottom.
  - Turn the coupling, so that afterwards the Allen screw (d) can be tightened through the aperture (g) with the Allen key (5).
  - Insert the pump-head assembly and push it on as far as it will go. Check the position of the Allen screw at the aperture (g).
- 6 Drive in the countersunk head-screws (4). Tighten the Allen screw (d) with an Allen key (5) size 2.5 mm.

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## **IDEX Health & Science**

Home of Upchurch Scientific &  
Ismatec Products

619 Oak Street  
Oak Harbor, WA 98277  
United States

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Ask for the latest  
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH