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Operating Manual

Mode d'emploi

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Deutsch

Gear pump with variable speed drive

REGLO-Z Analog
Drive without pump-head
Order No. ISM 895

REGLO-ZS Analog
Control unit with separate drive,
without pump-head
Order No. ISM 896

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Contents

Safety precautions	3–5
Warranty terms	6–7
Product	8
Operating panel	9
Rear panel	10
Mains voltage	10–11
Voltage setting	11
Changing the fuses	11
Starting the pump	12–13
Analog interface	14–16
Pumping against pressure	17
Before switching the pump off	17
Pump-head	18
Mounting the pump-head	19–20
Gear Pump Technology	21
Pump-head Specifications	21
Information on pump-heads	23
Pump-head with bypass	25
Accessories	23–27
Maintenance, Repairs	28
Spare-parts, Disposal	29
Technical specifications	30

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Safety precautions

ISMATEC® gear pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Prac- tice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion-proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when changing the pump-head.

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Safety precautions

- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.

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Safety precautions

- We advise you never pump liquid containing particles.
Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse-holder must not be shortcircuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

From date of delivery

REGLO-Z/-ZS drive: 2 years
Pump-head: 1 year
(without consumables)

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty Terms

We warrant the perfect functioning of the REGLO-Z/-ZS drive, provided it has been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

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Warranty Terms

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third-party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing.

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Product

Package contents

- Gear pump drive (REGLO-Z) or control unit with gear pump drive and connecting cable (REGLO-ZS) magnet included
- Power cord
- Operating manual

Accessories as ordered

- 1 Pump-head
- 2 Tubing adapters

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Operating panel

- 1 Mains switch and Start / Stop
- 2 Speed selector
99 = 5000 rpm (max. revolution)
- 3 Overload indicator (red LED)

Overcurrent protector

When an overload condition occurs, an activated overcurrent protector is indicated by the red LED (3) and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately**.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again) Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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⚠ Please note

The drive REGLO-Z runs only in clockwise direction.

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Rear panel

- 1 Analog interface
- 2 Mains socket
- 3 Fuse-holder with voltage selector 115/230V

⚠ Mains voltage

Mains voltage	Fuse rating
220-240 V _{AC}	2 x T800mA/250V
110-120 V _{AC}	2 x T800mA/250V

⚠ Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be earthed (protective conductor contact).

⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse-holder by opening it
with a small screw-driver (size 0).
 - 2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse-hol-
der so that the required voltage rating is
facing the window of the fuse-holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2 x T800mA/250V
115 V_{AC}: 2 x T800mA/250V
⚠ Use always 2 slow-blow fuses of
the same type complying with the
local mains voltage.
 - 4 Shut the fuse-holder. The voltage rating
is visible in the window.

⚠ **Before starting-up**
Check if the voltage setting visible in the window
of the fuse-holder complies with your local mains
voltage.

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
⚠ Please note

Before initial start-up, the pump-head must be filled with the medium or water. Never run the pump-head under dry conditions. **Never pump liquids that contain particles. If necessary, filter the medium before use.**

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Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on page 11.
- If necessary, mount the pump-head as stated on pages 19–21.
- Before screwing the tubing adapters into the pump-head, we recommend you **wrap PTFE tape around** the thread of the adapter.
- Connect the Tubing to your system. Fasten the tubing to the adapters by means of the tubing clamps.
- Connect the pump to the mains.
- Fill the pump-head with liquid before you start pumping.

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Starting the pump

- 1 Set the required speed on the 2-digit speed selector.
(99 = 5000 rpm)

➔ The rotation speed can be adjusted while the pump is running.

- 2 Press in the mains switch
= run/stop function

➔ The drive runs only in clock-wise direction.

⚠ Never run the pump-head under dry conditions.

⚠ Never pump liquids that contain particles. If necessary, filter the medium before use.

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
Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND).

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Input impedance and alternatives by means of DIP switches inside the pump (see page 16).

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
Analog interface

Pin 7, +22V_{DC}

About +22 V_{DC} are available (max.
current 0.2 A).

Pin 9, speed OUT

0–10 kHz

Frequency proportional to the rotation
speed: 0–5000 rpm.

Pin 10, +5V_{DC}

About +5 V_{DC} are available.
(max. current 0.1 A).

Pin 13, internal speed

Analog interface activated
(Pin 2 on GND).

- Pin 13 open: The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	18kΩ	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V	38kΩ	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250Ω	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250Ω	OFF	OFF	OFF	ON	ON	OFF

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Pumping against pressure

The max. differential pressure depends on the mounted pump-head model (see page 21).

Before switching the pump off

 If no non-return valve is mounted on the discharge side of the pump, a siphoning may occur at the moment the pump drive is stopped!
(see page 27)

If the pump will not be used for some time, it is strongly recommended to thoroughly rinse the pump-head with water.

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Pump-head

In general, the REGLO-Z gear pump is supplied as a complete unit, consisting of the REGLO-Z drive, and as ordered:

- 1 pump-head
- 2 driving magnet
- 3 tubing adapters (2 pcs.)

For the REGLO-Z, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on page 21-24.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

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Mounting the pump-head

Fig. 3

Screw the 2 tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, we recommend you **to wrap PTFE tape around** the thread of the adapter.

Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user and the arrow must point to the right.

➔ Please note that gear pumps can only pump in one direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ **Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.**

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Fig. 3

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Gear Pump Technology

»Cavity style« pump-head

- The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift.
- Excellent chemical resistance
- Smooth operation at a low noise
- Low inner friction
- Incl. B mount (flange)

»Suction shoe« pump-heads

- Enhanced pump performance at elevated differential pressures
- Fewer parts, easy to service
- Suitable for elevated temperature range
- Not suitable for applications requiring a suction lift
- Incl. B mount (flange)

Pump-head specifications

- Stainless steel pump-head housing SS 316/DIN 1.4571
- Temperature range
-46°C to +54°C / -51°F to +129°F (Cavity style)
-46°C to +177°C / -51°F to +350°F (Suction shoe)
(higher temperatures possible on request)
- Differential pressure
The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head.

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Pump-head specifications

- **Differential pressure**
The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head
- **System pressure**
The pressure within the pump-head max. 21 bar / 305 psi
(Exception: Z-1830, Z-150: 22 bar)
- **Gear material**
PTFE = Polytetrafluoroethylene
PPS = Polyphenylenesulphide
Graphite = Carbon graphite
PEEK = Polyetheretherketone
- **Ports**
(Inside thread) 1/8"-27 NPT
for tubing with 6.4 mm ID

The pump-heads are interchangeable. Different drives take different magnets. The REGLO-Z Analog drive is originally supplied with the correct magnet. When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC® representative. Some pump-heads are available with an internal **bypass**. This bypass is to be used for safety reasons only and must not be used as pressure control.

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
¹⁾ Flow rate without differential pressure

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Cavity Style

	Z-120	Z-140	Z-140 P	Z-130	Z-150	Für korrosive Medien For corrosive media Pour milieux corrosifs	
Pumpenkopf Nr. Pump-head No. Tête de pompe No.						Z-140 HC	Z-150 WI
Fließrate / Flow rate / Débit (ml/min) max. ¹⁾	3200	3200	3200	3200	3200	3200	3200
Fließrate / Flow rate / Débit (ml/min) min. ¹⁾	32	32	32	32	32	32	32
Differenzdruck max. (bar) Pression différentielle Differential pressure	2)	2)	2)	2)	2)	2)	2)
Zahnradmaterial Gear material Matériaux engrenage	PTFE	PTFE	PEEK	PPS	PPS	PTFE	NiC
Dichtungen Seals Étanchéités	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE	PTFE
Gehäusematerial Housing material Matériaux caisse	SS316	SS316	SS316	SS316	SS316	Hastelloy- C276	Surface hardened
Systemdruck max. (bar) System pressure Pression de système	21	21	21	21	21	21	21
Betriebstemperatur Operating temperature Temperature operation	-46...+54	-46...+54	-46...+54	-46...+54	-46...+54	-46...+54	-46...+54
Bypass Internal bypass Bypass interne	✓	--	--	✓	--	--	--
Anschlüsse (Innengewinde) Ports (internal thread) Ports (internal thread)	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT	1/8" - 27NPT
Bestell-Nr. Order No. No. de commande	MI0013	MI0016	MI0313	MI0019	MI0020	MI0284	MI0265

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Operating range for cavity style pump-heads
on the REGLO-Z Analog drive: see diagram
For higher pressures use the MCP-Z drive!

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
Pump-head with bypass

Some gear pump-head models have an internal fluid bypass.

- A fluid bypass valve helps protect against decoupling and damage from high pressure build-up.
- It allows adjustment of the max. differential pressure (from 0.7 bar up to the max. pressure valid for the individual pump-head).
- It should only be used for safety reasons and not for pressure control (bypass conditions create a temperature rise causing swelling in PTFE gears).

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
Bestell-Nr. AR 0001

Accessories

Tubing adapter 01 for all gear pump-heads

1/8" thread NPT
for tubing i.d. 6-7 mm (1/4")
(always order 2 pcs.)

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Foot-switch

Order No. ISM 891

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

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Please note

For maintenance, service, repairs and warranty of the pump-head we also refer to the enclosed instruction manual of the pump-head manufacturer. In case of doubt that manual shall prevail.

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Maintenance

We strongly recommend you thoroughly rinse the pump-head with water after use.

Provided the REGLO-Z/-ZS drive is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear by himself (see page 26).

Repairs

For repairs please send the defective REGLO-Z/-ZS drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., model.

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Spare-parts

If required your ISMATEC® representative will be pleased to send you on request:

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Drive

Motor type	DC-Motor
Speed range	50–5000 rpm
adjustable in 1 % steps	
Remote control	(see page 14–16)

Mains connection

Fuse rating

230 V 50 Hz	2 x T800mA/250V
115 V 60 Hz	2 x T800mA/250V

Protection rating	IP 30
-------------------	-------

Power consumption	max. 50 W
-------------------	-----------

Operating conditions

Temperature	41 to + 104°F
Rel. humidity	max. 80%
- not condensing, at normal laboratory conditions.	

Differential pressure

according to pump-head

Dimensions/Weight

	DxWxH (inches)
Pump drive	6.9x3.9x5.3
Control unit	6.4x2.6x3.2
	Weight (lb)
Pump drive	4.4
Control unit	3.5

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IDEX Health & Science
Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

IDEX Health & Science
Ismatec Laboratoriumstechnik
GmbH
Futtererstrasse 16
97877 Wertheim
Germany
Tel: +49 9377 9203-0
Fax: +49 9377 1388
fluidics.europe@idexcorp.com

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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