REGLO-Z /-ZS *Digital*

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Schutzgrad
IP30

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® drives are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions. In particular, it is in the user's responsibility not to exceed the specified maximum system pressure for the pump-head.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion-proof chambers or in the presence of flammable gases or fumes
- The pump must be switched off when changing the pump-head.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.

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Please note

Never pump media which contain particles!

If necessary, filter the media before starting your application.

Never run the pump-head under dry conditions!

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Safety precautions

- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- We advise you never to pump liquid containing particles. Even the smallest particles have an abrasive effect, which can ruin the gears in the course of a few seconds.
- Before initial start-up, the pump-head must be filled with the pumping liquid or water. Never run the pump-head under dry conditions.
- Only new fuses, according to the specifications stated on page 8 in this manual, must be used.
- The fuse-holders must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the potential hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Pump drive: 2 years

Pump-head: 1 year

The pump-heads are covered by the specific terms of Micropump Ltd. Gears and seals are regarded as consumables and as such they are not covered by the warranty.

In case of any queries, please contact your nearest ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Contents of the package

- **REGLO-Z *Digital*** drive magnet included
Order No. ISM 901
- 1 power cord with an IEC 320 connector (female) and a country specific mains plug
- 2 tubing clamps
- 1 non-return valve (PP) (For user instructions see page 21 to 22)
- Operating manual
- **REGLO-ZS *Digital*** control unit with separate drive, magnet and connecting cable included
Order No. ISM 1143
- 1 power cord with an IEC 320 connector (female) and a country specific mains plug
- 2 tubing clamps
- 1 non-return valve (PP) (For user instructions see page 21 to 22)
- Operating manual

If ordered:

- Pump-head with 2 tubing adapters
- Software
- other accessories

Please check the package

and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

➔ Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 RS232 IN (female) for operation controlling via PC
- 2 RS232 OUT (male) for cascading up to 8 pumps
Digital input (TTL level) Run/Stop, autostart
Digital output 0–10 kHz (speed output)
- 3 Mains socket
- 4 Fuse-holder

⚠ Mains voltage

Mains voltage	Fuse rating
100–240 V _{AC} 50–60 Hz	1 x T2A/250V

⚠ Socket/Power cord

Use exclusively the originally supplied power cord. **The socket must be earthed** (protective conductor contact).

Changing the fuses

- ➔ Switch the pump off, pull out the mains plug.
- 1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).
 - 2 Remove fuse (a) and replace it with the spare fuse (b): 1 x T2A/250V.

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
Starting the pump

- Mount the pump-head as stated on pages 44 and 45.
 - Wet the Pump inside, e.g. with a washing bottle
 - Before screwing the tubing adapters into the pump-head, we recommend you wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.
 - Connect the tubing to your system. Fasten the tubing to the adapters by means of the tubing clamps.
 - Connect the pump to the mains and switch it on with the power supply switch.
- ➔ Enter the ID code of the mounted pump-head in the settings-menu (see basic settings, page 16).

Socket/Power cord
Use exclusively the original power cord supplied with the pump. The socket must be earthed (protective conductor contact).

Before you start pumping, fill the pump-head with the medium or water. Never run the pump-head under dry conditions.

Never pump liquids that contain particles. If necessary, filter the medium before use.

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
Operating panel

(* = Settings page 12)

- 1 On/Off switch
- 2 Digital LED display
- 3 Increase value
- 4 Reduce value
* settings
- 5 RUN/STOP (double function)
 - a Starts and stops the pump
 - b Interrupts and continues a dispensing cycle (page 23)
➔ The function Start/Stop can also be triggered by an external device such as a foot switch. (Accessories on page 43).
- 6 RESET (double function)
 - a Rotation direction
(can only be selected for cavity style pump-heads)
 - b RESET
 - Resets the accumulated volume in mode TOTAL to zero (page 18)
 - Interrupts dispensing cycle for good (page 23)
- 7 CAL/MAX (double function)
 - a CALL = Calibrating key for the flow rate or dispensing volume
 - b MAX = Max. speed in the modes PUMP•Flow rate and TOTAL (keep the MAX key pressed down while the pump is running)

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Operating panel

(* = Settings page 12)

- 8 MODE key
Selecting operating mode
* default
- 9 LEDs for active operating mode
- **PUMP • Flow rate**
Pumping by flow rate (ml/min)
 - **DISP • Volume**
Dispensing by volume (ml)
 - **DISP • Time**
Dispensing by time
 - **PAUSE • Time**
Pause time for dispensing at intervals
 - **TOTAL**
Pumping with Read-out of totally delivered volume

Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Displays the operating mode last used

⚠ **Before starting the pump for the first time or after connecting a new pump-head,** the ID-code of the pump-head must be entered (see basic settings, page 16).

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
Control keys for basic settings

a settings
Allows access to the menu settings

b ok
Confirms a data entry

c escape

- Cancels an entry
- Quits the menu settings

d default
Sets a value to its default value:

- Cycles (dispensing cycles) 0
- Backsteps 0
- Address 1
- Foot switch tOGL
- Autostart AStb
- Digits (display) 3
- Display Speed FLOU
- Head (pump-head) N23

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⚠ Please note

The settings menu requires a break of **5 seconds** between power off and on.

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Selecting the basic settings

- 1 Switch the pump off and wait for 5 seconds
- 2 Keep the settings key pressed and switch the pump on
- 3 Select the basic settings by using the ▲▼ keys
 - Cycles (dispensing cycles)
 - Back-steps (number of back-steps for volume dispensing)
 - Address (RS232, pump address)
 - Foot switch (Toggle/Direct)
 - Autostart (Standby/RUN)
 - Digits (Number of digits displayed for flow rate and dispensing volume)
 - Display Speed (Type of speed display)
 - Head = ID-Code (pump-head)
 - Init (resets all basic settings to the default values)
- 4 ok = confirm your selection
- 5 Select the value by using the ▲▼ keys, or
 - press escape for cancelling wrong input
 - the values can be reset to the default value by using the default key

- 6 ok = confirm entered value
- 7 escape
Quit menu basic settings

➔ **Basic settings**
refer to the following pages

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Number of dispensing steps 0..9999
(0 = Default value, infinity)

Backsteps

Number of back-steps (for drip-
free dispensing in the mode DISP
Volume)
from 0 .. 100 (0 = Default value)

Address (RS232)

Pump address 1–8
(page 34 and 41)
(1=default value)

Foot switch

- Toggle (default value)
 - Step on foot switch
= pump starts running
 - Press foot switch again
= pump stops
- Direct
the pump runs as long as the
foot switch is pressed

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
Basic settings

Configuring Autostart

Autostart determines how the pump will react after being switched on, e.g. after a power failure, or if the pump is controlled via the electric power supply.

- Standby
Pump is stand-by when switched on (default value)
- RUN
Pumpe runs when switched on

Digits (display)

Number of significant digits, 3 or 4 (3 = default value) when displaying flow rate and dispensing volume

Display Speed

Configuring of the pump speed in mode Pump•Flow rate

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Basic settings

Head Identification code (ID) of the pump-head according to the spreadsheet on page 17 (N23 = default value).

➔ In order to obtain correct pumping and dispensing values, the individual ID code of the mounted pump-head must be entered; especially when a pump is used for the first time and each time the pump-head is changed.

Init (Initializing)

Pressing the ok key resets the following basic settings to the default values:

Mode	PUMP•Flow rate
Speed	400 rpm
Dispensing volume	30 revs
Dispensing time	4.50 sec
Pause time	2.00 sec
Rotation direction	clockwise
Back-steps	0
Number of dispensing cycles	0
Pump-head	N23
Digits (display)	3
Foot switch	Toggle
Autostart	Standby
Display	Flow rate

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Chart for basic settings "Head" (see also table on page 47)

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The pump-head can be configured according to either the Micropump or the ISMATEC convention. The resulting settings are equivalent.

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TOTAL volume

Change with the MODE key to TOTAL. The totally delivered volume is displayed.

Four display readings are available:

Microlitre:	e.g.	17.0 µl = 17.0µ
Millilitre:	e.g.	1700 ml = 17.00
Litre:	e.g.	17.0 l = 17.0L
Litre:	e.g.	1700 l = 1700.

For zero-setting, either

- press the reset key once, or
- switch the pump off.

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Pumping by flow rate

- 1 Change mode to PUMP • Flow rate
- 2 Enter the required flow rate with the ▲▼ keys (displayed in ml/min) Start pump with RUN/STOP
➔ The flow rate can also be adjusted while the pump is running.
➔ For an accurate flow rate we recommend to calibrate the pump. (see page 20)
- 3 For fast filling or emptying the system, keep pressing the CAL/MAX key while the pump is running.

➔ During the pumping process the following values can be retrieved with the MODE key:

- Flow rate (PUMP rpm)
- Totally pumped volume (TOTAL)

⚠ Entering the ID-Code of the currently mounted pump-head (page 17) allows you set the flow rate in ml/min. This, however, is an approximate value as not yet calibrated.

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
Calibrating the flow rate

- 1 Change the mode to PUMP • Flow rate. Enter the required flow rate with the ▲▼ keys.
- 2 Change the mode to DISP Time.
 - Enter 60 seconds by using the ▲▼ keys
 - Start the pump with RUN/STOP
 - ➔ The pump stops automatically after 60 sec.
 - ➔ Measure the dispensed liquid by volume or weight and calibrate the ascertained value as follows:
- 3 Return to the mode PUMP • Flow rate.
- 4 Press the CAL/MAX button (displayed value blinks) Enter the weighed or measured value with the ▲▼ keys.
- 5 Save with the CAL/MAX key (the flow rate setting returns automatically to the initially entered set point)
- 6 Start with the RUN/STOP key
 - ➔ Depending on the application conditions repeated calibration may be necessary.
 - ➔ During calibration, the following keys can be used:
 - ok
 - escape
 - default

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Medium: Water
Flow rate: 1000 mL/min
Pumphead: 140
Set Dispensing Volume 20 ml

Measuring result with	0 back-steps	15 back-steps
	19.950	20.042
	19.971	20.044
	19.971	20.042
	19.949	20.059
	19.932	20.046
	20.005	20.022
	19.959	20.007
	19.933	20.066
	19.975	20.036
	19.954	20.031
	19.981	20.003
	20.013	20.011
	19.981	20.063
	20.017	20.038
	19.994	20.079
	19.934	19.975
	20.003	20.019
	19.972	20.018
	19.929	20.041
	19.977	20.026
	19.948	19.982
	19.934	20.055
	19.924	20.007
	19.942	19.971
Mean value	19.965	20.029
Standard deviation (abs.)	0.026	0.026
Standard deviation (rel.)	0.14%	0.14%
Maximum	20.017	20.079
Minimum	19.924	19.971
Range abs. (Range) +/-	0.047	0.054
Range rel. (Range) +/-	0.12%	0.13%

Medium: Water
Flow rate: 1022 mL/min
Pumphead: 140
Set Dispensing Volume 100 ml

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Dosing with a valve

Precise and reproducible dosing with a gear pump requires a check valve on the discharge side. Such a valve is originally supplied as standard equipment.

Drip-free dispensing with back-steps

Back-steps can be programmed in the "Setting" menu (see page 14). At the end of each dispensing cycle the pump stops for a moment and afterwards runs back in the reverse direction by the pre-set number of back-steps. Consequently, the liquid is sucked back into the tube, which avoids the uncontrolled production of a drop at the tip of the tubing. Additionally, the reversing of the liquid presses the diaphragm of the check valve onto the valve seat which consequently avoids dripping.

Depending on the type of liquid and the volume size, the programming of back-steps can favour the repeatability of the dispensing volume. The correct number of back-steps depends on the volume and viscosity of the liquid.

Tests carried out with water have proved that about 15 back-steps are ideal for dosing volumes between 5 and 100 ml. With a set dosing volume of 20 ml results have been achieved with a standard deviation of as little as 0.14% and with 100 ml the standard deviation came even down to 0.05% (see fig. on the left).

➔ For accurate dosing the "Cavity Style" pump-heads are better suitable than the "Suction Shoe" models. The use of the back-steps mentioned above is only effective with the "Cavity Style" pump-heads. That is why for dosing applications, we only recommend using the pump-head models 120, 130, 140 and 150.

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
Mounting the tubing and check valve

In order to achieve optimum results, the following instructions for correct mounting are essential.


- The check valve must be mounted in a vertical position in the flow direction.
 ➔ On the check valve the flow direction is marked by two arrows. These arrows must point upwards.
- The liquid level of the reservoir on the inlet side of the pump must be in a lower position than the check valve. Ideally, the suction tube is positioned on a horizontal level to the pump or is slightly rising towards the pump-head.
- Movements on the tubing generate tensions which directly influence the volume of the liquid and thus have a negative effect on the repeatability. Therefore, we recommend fixing the tubing and check valve on the discharge side to a stand.
- In order to stop any draining we further recommend tapering the tip of the tubing with a tube adaptor. If necessary, we also recommend the programming of back-steps (see page 14).

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Dispensing by volume

- 1 Change the mode to DISP Volume
- 2 Use the ▲▼ keys for entering the required dispensing volume. Four display readings are available:

Mikroliter:	e.g.	17.0 µl	=	17.0µ
Milliliter:	e.g.	17.0 ml	=	17.0
Liter:	e.g.	17.0 l	=	17.0L
Liter:	e.g.	1700 l	=	1700.

➔ For dispensing accurately we recommend you calibrate the pump (page 24)

- 3 Start the pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the modes PUMP flow rate.

➔ With the ▲▼ keys the dispensing volume can be changed even during the dispensing process. The newly entered set point appears shortly in the display.

➔ During the dispensing process the following values can be retrieved with the MODE key:

- Flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the rotation speed or flow rate can be changed even during the dispensing process.

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⚠ Please note

➔ Use a check value for highly accurate dispensing.

With RUN/STOP a dispensing cycle can be interrupted and resumed.


With RESET an interrupted dispensing cycle can be definitely stopped.

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Calibrating the volume

- 1 Change the MODE to DISP Volume.
 ➔ Use the ▲▼ keys for entering the required dispensing volume
- 2 Start the pump with RUN/STOP.
 ➔ Pump stops automatically
 ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes
 ➔ Measure the dispensed liquid by volume or weight
- 3 Stay in the same mode (DISP Volume) and press the CAL/MAX-key (display blinks).
- 4 Enter the effectively dispensed volume with the ▲▼ keys.
- 5 Save with the CAL/MAX key.
 ➔ Pump is calibrated and returns automatically to the initially entered set point

- 6 Start the pump with RUN/STOP
 ➔ Depending on the application conditions repeated calibration may be necessary.
 ➔ During calibration, the following keys can be used:
 - ok
 - escape
 - default

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Dispensing by time

The dispensing time can be entered from 0.1s to 999h.

- Change the mode to DISP•Time
- Enter the required dispensing time with the ▲▼ keys
➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds: 0.1–899.9 (in 0.1s steps)
Minutes: 15'–899' (in 1min steps)
Hours: 15h–999h (in 1h steps)

- Start pump with RUN/STOP

➔ With the ▲▼ keys the dispensing time can be changed even during the dispensing process.

➔ During the dispensing process the following values can be retrieved by pressing the MODE-key:

- flow rate (PUMP Flow rate)
- totally dispensed volume (TOTAL)

➔ With the ▲▼ keys the speed and flow rate can be changed even during the dispensing process.

➔ Depending on the pump-head and the application, a very short dispensing time can result in dispensing volumes which are not reproducible.

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
➔ At the same time, the volume is also reset to the default value.

Default calibration of volume

- 1 Change the mode to DISP Volume.
- 2 Press the CAL/MAX button (displayed value blinks).
- 3 Press the default button (displayed value blinks).
- 4 Confirm with the CAL/MAX button.

➔ At the same time, the flow rate is also reset to the default value.

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
Dispensing by volume within a pre-set time

- 1 Change the mode to DISP Volume. Enter the required volume in ml with the ▲▼ keys
 - 2 Change the mode to DISP Time. Enter the required time with the ▲▼ keys (0.10s – 999h).
 - 3 Push the MODE button.
➔ Both LEDs light up simultaneously
 - 4 Start the pump with RUN/STOP.
➔ The pump dispenses the required volume within the pre-set time
- a Changes / Corrections**
Press the Reset key
➔ For changing the dispensing time or dispensing volume, the pump returns directly to the mode DISP Volume
➔ Repeat steps 1 to 3

b Error message

LEDs are blinking alternately

9999 = volume too large
time too short

1111 = volume too small
time too long

Despite volume and time adjustments, the LEDs keep on blinking:
➔ change the pump-head model in order to obtain another flow rate range.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time:

- Change mode to PAUSE Time. Enter the pause time (between 0.10s and 999h) with the ▲▼ keys.
- Press the MODE-key once ➔ Both LEDs for PAUSE and DISP Volume light up simultaneously. Enter the volume with the ▲▼ keys.
- Start the pump with the RUN/STOP key.
 - ➔ Pump dispenses at intervals until the RUN/STOP button is pushed again.
 - ➔ Active mode blinks.

➔ With the ▲▼ keys the dispensing volume or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

➔ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys. Furthermore, the total volume is displayed.

If necessary, calibrate the volume as indicated on page 24.
For entering the number of dispensing cycles see page 30.

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
Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time:

- Go to the mode PAUSE Time
Enter the pause time with the ▲▼ keys (between 0.10s and 999h).
- Press the MODE key twice
→ Both LEDs for PAUSE and DISP Time light up simultaneously.
Enter the dispensing time with the ▲▼ keys (0.10s to 999h).
- Start the pump with the RUN/STOP key
→ Pump dispenses at intervals until the RUN/STOP button is pushed again.

→ Active mode blinks.

→ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.

- During the dispensing process the following values can be retrieved with the MODE-key:
- flow rate (PUMP Flow rate)*
 - totally dispensed volume (TOTAL)
- * The flow rate can be changed with the ▲▼ keys.

For entering the number of dispensing cycles see page 30.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

See basic settings on page 14.

➔ During the pauses the display shows the remaining number of dispensing cycles.

➔ During long pauses, the remaining interval time can be checked with the CAL/MAX key and adjusted by using the ▲▼ keys.

➔ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key. When resuming the process, the pump continues with the remaining number of dispensing cycles.

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Pumping against pressure

The max. differential pressure depends on the mounted pump-head model (see page 46).

Pump-head with bypass

Some gear pump-head models have an internal fluid bypass.

- A fluid bypass valve helps protect against decoupling and damage from high pressure build-up.
- It allows adjustment of the max. differential pressure (from 0.7 bar up to the max. pressure valid for the individual pump-head).
- It should only be used for safety reasons and not for pressure control.

➔ Bypass conditions create a temperature rise causing swelling in PTFE gears.

Before switching the pump off

⚠ If no non-return valve is mounted on the discharge side of the pump, a siphoning may occur at the moment the pump drive is stopped! (See page 22)

If the pump will not be used for some time, it is strongly recommended to thoroughly rinse the pump-head.

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Overcurrent protector

The drive REGLO-Z *Digital* features an overload protector. It disables the pump drive in case of a motor blockage. When an overload condition occurs, it is indicated in the display by the letters 'OL' and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately**. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure, etc.). Wait for 5 seconds before switching on the pump again.

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
⚠ Please note

The analog interface is integrated in the port 'RS232 OUT' and can be used for remote control (e.g. by foot-switch).

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Analog interface

Pin 4: +5VDC

For connecting an external control unit (current <100mA).

Pin 5: GND

Point of reference for all signals described.

Pin 6: start

Connection for an external trigger device (foot-switch etc.). For operating modes, refer to Basic settings „Foot switch“ page 14.

Pin 8: busy

Informs about pump status.

Pumpe in operation

= output +5V

Pumpe in stand-by

= output on ground

(The output resistance is 4.7 kΩ)

Pin 9: speed OUT

A square-wave signal with a frequency proportional to the speed (rpm) is at your disposal:
120 pulses are generated per one pump gear revolution.

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
Please note

When switching on the REGLO-Z *Digital* pump, the serial interface will only respond after a delay of 5 seconds.

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1–8).

It is factory-set at 1 and can be changed with @.

The address will be stored permanently (even after the pump has been switched off).

The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface. (Page 41)

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed** by a carriage return (ASCII 13). The pump confirms most of the commands with an asterix *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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
RS232 OUT (male)

Pin 2: RS232 Rx

Is used for connecting additional pumps together with pin 3/5.

Pin 3: RS232 Tx

Is used for connecting additional pumps together with pin 2/5.

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address. (Page 34)**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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
Each pump requires an individual address. (see page 34)

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Bestell-Nr. XC0009

Programming Software

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.
<http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our website:
<http://www.ismatec.com>

Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments

Order No. IM0015C168

Connecting cable 9F / 9M

For connecting an REGLO-Z *Digital* to a PC via the 9-pin RS232 connector. (2 m)

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
Non-return valve (Polypropylene)

(included with the drive)

➔ should always be mounted in a vertical position

Order No. SD 0007

Tubing clamps

included with the drive (2 pcs.)

Order No. BN 1001

Tubing nipple 01

for all gear pump-heads

thread 1/8" NPT

for tubing i.d. 6–7 mm (1/4")

(always order 2 pcs.)

Order-No. AR 0001

Foot switch

Order No. ISM 894

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling bottles etc.

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Pump-head

In general, the REGLO-Z *Digital* gear pump is supplied as a complete unit, consisting of the REGLO-Z *Digital* drive, and as ordered:

- 1 pump-head
- 2 driving magnet
- 3 tubing adapters (2 pcs.)

For the REGLO-Z *Digital*, a wide range of different gear pump-heads with PTFE, PPS or graphite gears are available.

You find an overview chart of the pump-heads available on page 47.

⚠ Before mounting a new pump-head, make sure that the correct magnet is fitted to the drive.

Mounting the pump-head

Fig. 1

The following tools are required for mounting a pump-head:

- Allen key 3/32"
- Phillips screwdriver (Size 0)

Remove all screws from the fixing case (a) of the drive.

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Mounting the pump-head

Fig. 2

Fasten the driving magnet on the motor shaft (through the hole at the top).

➔ Ensure that there is a gap of about 3 mm between the magnet and the brim of the fixing case.

⚠ During operation the back of the magnet must not touch the housing of the pump drive.

Fig. 3

Screw the 2 tubing adapters into the pump-head.

➔ Before screwing the tubing adapters into the pump-head, we recommend you wrap **PTFE tape around the thread of the adapter**. This ensures optimum sealing and helps to protect the thread.

Fig. 4

Screw the pump-head to the fixing-case of the drive. The face of the pump-head with the engraved model number must be directed to the user. The engraved arrow indicates the preferred flow direction.

➔ Please note that gear pumps with suction shoe can only pump in one direction. Pump-heads “Cavity style” pump in both directions but with reduced performance in reverse direction. An arrow engraved on the front of the pump-head indicates the correct flow direction.

⚠ Fill the pump-head with liquid before you start pumping. Avoid dry-running of the pump-head.

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Fig. 2

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»Cavity style« pump-head

- The »Cavity style« pump-head can be used for viscous media and for applications with a certain suction lift
- Excellent chemical resistance
- Smooth operation at a low noise
- Low inner friction
- Incl. B mount (flange)

»Suction shoe« pump-heads

- Enhanced pump performance at elevated differential pressures
- Fewer parts, easy to service
- Suitable for elevated temperature range
- Not suitable for applications requiring a suction lift
- Incl. B mount (flange)

Pump-head specifications

- Stainless steel pump-head housing SS 316/DIN 1.4571
- **Temperature range**
-46°C to +54°C / -51°F to +129°F (Cavity style)
-46°C to +177°C / -51°F to +350°F (Suction shoe)
(higher temperatures possible on request)
- **Differential pressure**
The difference between the pressure on the inlet and the discharge port of the pump-head max. 5.2 bar (75 psi) depending on pump-head (page 47)
- **System pressure**
The pressure within the pump-head max. 21 bar / 305 psi
(Exception: Z-1830, Z-150: 22 bar)
- **Gear material**
PTFE = Polytetrafluoroethylene
PPS = Polyphenylenesulphide
Graphite = Carbon graphite
PEEK = Polyetheretherketone
- **Ports**
(Inside thread) 1/8"-27 NPT
for tubing with 6.4 mm ID

The pump-heads are interchangeable. Different drives take different magnets. The REGLO-Z/-ZS drive is originally supplied with the correct magnet. When ordering a new pump-head, please state the intended operation temperature. If you intend to use a Cavity style pump-head for temperatures exceeding 54°C (130°F), the correct sealing material must be ordered. Please contact your nearest ISMATEC® representative. Some pump-heads are available with an **internal bypass**. This bypass is to be used for safety reasons only and must not be used as pressure control (see page 31).

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
Operating range for cavity style pump-heads

on the REGLO-Z *Digital* drive: see diagram. For higher pressures use the MCP-Z drive!

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Contain the consummable parts (bushes, seals, gears)

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Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

After use, we strongly recommend you thoroughly flush the pump-head.

Provided the REGLO-Z *Digital* is operated properly and in compliance with this manual, the moving parts such as gears, etc. are the only parts that are subject to wear and tear.

The Service Kit corresponding to the individual pump-head allows the user to replace the parts that are subject to wear and tear by himself (see page 47).

Repairs

For repairs please send the REGLO-Z *Digital* drive to your ISMATEC® representative.

Please give information on defect, date of purchase, serial-no., and model.

Spare-parts

For repairs you intend to carry out by yourself (out of the warranty period) please ask your ISMATEC® distributor for:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical specifications

Drive

Motor type DC-Motor
Speed range 50–5000 rpm
digitally adjustable in steps
of 1 rpm

Differential pressure

max. 5.2 bar (75 psi) depending on
pump-head (page 47)

Remote control

via RS232 and analog interface
(pages 33 and 42)

Mains connection

100–240 V_{AC} / 50–60 Hz

Fuse rating power supply

1 x T2A/250V

Power consumption

Max. 140 W

Protection rating

IP 30

Operating conditions

Temperature +5 to +40°C
Rel. humidity max. 80%
at normal environmental conditions

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conformément à:
EN 61326-1 / EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

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Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH