**ISMATEC®**

**IDEX**  
HEALTH & SCIENCE

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### REGLO Analog

#### 2 Kanäle

ISM 830	6 Rollen
ISM 829	8 Rollen
ISM 795	12 Rollen

#### 4 Kanäle

ISM 828	6 Rollen
ISM 827	8 Rollen
ISM 796	12 Rollen

Deutsch

## Tubing pump

### REGLO Analog

#### 2 channels

ISM 830	6 rollers
ISM 829	8 rollers
ISM 795	12 rollers

#### 4 channels

ISM 828	6 rollers
ISM 827	8 rollers
ISM 796	12 rollers

English

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### REGLO Analog

#### 2 canaux

ISM 830	6 galets
ISM 829	8 galets
ISM 795	12 galets

#### 4 canaux

ISM 828	6 galets
ISM 827	8 galets
ISM 796	12 galets

Français

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**Contents**

Safety precautions	4–6
Warranty terms	7
Product	8
Rear panel	10
Mains voltage	10
Voltage setting	11
Changing the fuses	11
Inserting the tubing	12
Starting the pump	13
Overload protector	15
Pumping against pressure	16
When the pump is not in use	17
Running-in period for tubing	17
Tubing life	18

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**Contents**

Analog interface	19
Accessories	23
Flow-rates per channel	25
Maintenance	28
Spare-parts	28
Repairs	29
Disposal	29
Technical specifications	30

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## Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
  - for medical applications on human beings
  - in explosion proof chambers or in the presence of flammable gases or fumes
- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.

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## Safety precautions

- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.

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Caution: rotating pumphead creates a pinch and crush hazard.

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## Safety precautions

- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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## Warranty

**For all parts manufactured by ISMATEC®**

**from date of delivery: 2 years**

**All other parts,**

**excluding consumables: 1 year**

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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## Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packing or in a packet of equal quality. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded.

Shipping costs are charged to the customer.

### Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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## Product

The tubing pump REGLO *Analog* is with 2 or 4 tube channels and 6, 8 or 12 pump rollers available.

### REGLO *Analog* 2 channels

(Speed 3.2–160 rpm)

Order Number	Type	Flow rates ml/min
<b>ISM 830</b>	MS-2/06	0.005–68
<b>ISM 829</b>	MS-2/08	0.004–57
<b>ISM 795</b>	MS-2/12	0.003–38

### REGLO *Analog* 4 channels

(Speed 2.0–100 rpm)

Order Number	Type	Flow rates ml/min
<b>ISM 828</b>	MS-4/06	0.003–43
<b>ISM 827</b>	MS-4/08	0.003–35
<b>ISM 796</b>	MS-4/12	0.002–24

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## Product

### Package contents

- Pump with 2 or 4 cassettes (type as ordered)
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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## Rear panel

- 1 Analog interface
- 2 Mains socket
- 3 Fuse-holder with voltage selector 115/230V

### ⚠ Mains voltage

Mains voltage	Fuse rating
220-240 V <sub>AC</sub>	2x T500mA/250V
110-120 V <sub>AC</sub>	2x T500mA/250V

### ⚠ Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be earthed (protective conductor contact).

### ⚠ Before starting-up

Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
## Voltage setting 115V / 230V and changing the fuses




- ➔ Switch the pump off, pull out the mains plug.
- 1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).
- 2 Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
- 3 Insert 2 new fuses  
230 V<sub>AC</sub>: 2x T500mA/250V  
115 V<sub>AC</sub>: 2x T500mA/250V  
⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.
- 4 Shut the fuse-holder. The voltage rating is visible in the window.

⚠ **Before starting-up**  
Check if the voltage setting visible in the window of the fuse-holder complies with your local mains voltage.

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Ersatz-Kassetten MS/CA Click'n'go  
Bestell-Nr. IS 3510

## Inserting the tubing

- 1 ➔ Switch the pump off
- 2 Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously
- 3 Insert the 3-stop color-coded tubing with one collar into the cassette
- 4 Let the tubing hang down (prevent it from being twisted)
- 5 Insert the tubing with the second collar at the other end of the cassette
- 6 Reinsert the cassette into the roller-head
- 7 ➔ Switch the pump on

➔ When the pump is idle,  
release all cassettes at the fixing-  
tongue (2).

Spare-cassettes MS/CA Click'n'go  
Order-No. IS 3510

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Cassettes suppl. MS/CA Click'n'go  
No de commande IS 3510

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## Starting the pump

- Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on page 10–11.
- Connect the pump tubing to the system.
- Connect the pump to the mains.

⚠ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you to wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

➔ For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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## Starting the pump

- 1 Set the required speed on the 2-digit speed selector.  
99 = 100 rpm or 160 rpm  
(max. revolution)  
➔ The rotation speed can be adjusted while the pump is running.
- 2 Choose the rotation direction
- 3 Start the pump with the main switch  
= run/stop function
- 4 Overload indicator (red LED)

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## Overcurrent protector

The tubing pump REGLO *Analog* features an overload protector. When an overload condition occurs, an activated overcurrent protector is indicated by the red LED (4) and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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### Pumping against pressure


The REGLO *Analog* can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures).

In case of any doubts please contact your ISMATEC® agent.

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### When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure. Releasing the cassette on the right side is sufficient.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

#### **⚠ Syphoning effect**

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

### Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

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## ⚠ Please note

Please refer to our detailed tubing documentation.

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## Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our detailed tubing documentation (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
**Digitale Eingänge (TTL-Pegel)**

**Digital inputs (TTL-level)**

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- Pin 2, remote
- Pin 3, start
- Pin 4, direction
- Pin 13, internal speed

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**Analog interface**

**Pin 1, GND (ground)**

Reference potential for all other inputs.

**Pin 2, remote**

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

**Pin 3, start**

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND).

**Pin 4, direction**

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

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
External speed control

(0–5V, 0–10V, 0–20mA, 4–20mA)

Input impedance and input range can be selected via a dip-switch inside the pump (see p. 22).

### Input impedance

0...5 V 18 k $\Omega$

0...10 V 38 k $\Omega$

0...20 mA 270  $\Omega$

4...20 mA 270  $\Omega$

### Pin 7, +20V<sub>DC</sub>

About +20 V<sub>DC</sub> are available (max. current 0.2 A).

### Pin 9, speed OUT

Frequency proportional to the rotation speed:

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
## Analog interface

### Pin 10, +5V<sub>DC</sub>

About +5 V<sub>DC</sub> are available  
(max. current 0.1 A).

### Pin 13, internal speed

Analog interface activated  
(Pin 2 on GND)

- Pin 13 open: The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
## Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

**⚠ The instrument should only be opened by a qualified technician!**  
Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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
Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	18k $\Omega$	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V	38k $\Omega$	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250 $\Omega$	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250 $\Omega$	OFF	OFF	OFF	ON	ON	OFF

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### Foot switch

Order No. ISM 891

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

### 3-stop tubing

For selecting the correct tubing, please refer to our detailed tubing documentation.

➔ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you prime the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

For applications with a high differential pressure (> 1 bar) we recommend the use of cassettes with pressure lever.

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
MS/CA Anpresshebel  
MS/CA pressure lever  
MS/CA levier de pression

## Zubehör

### Ersatz-Kassetten

Kassette MS/CA Click'n'go Material  
**POM**

Bestell-Nr. **IS 3510**

Kassette MS/CA Anpresshebel  
Material **POM**

Bestell-Nr. **IS 0649**

Kassette MS/CA Anpresshebel  
Material **PVDF**

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## Accessories

### Spare-cassettes

Cassette MS/CA Click'n'go Material  
**POM**

Order No. **IS 3510**

Cassette MS/CA pressure lever  
Material **POM**

Order No. **IS 0649**

Cassette MS/CA pressure lever  
Material **PVDF**

Order No. **IS 3629**

➔ The cassettes with pressure lever are available on request. This type of cassette may provides better results both for the Tygon MH tubing and at elevated differential pressure conditions.

➔ The cassettes allow the user to insert tubes with different diameters and materials on the same roller-head.

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\*) Bestell-Nr. = Tygon® ST (R3603)

### △ Flow rates per channel

In the tables listed subsequently the values indicated are only approximate and determined as follows: ml/min, per channel, with water and Tygon® tubing, without differential pressure.

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### ⚠ Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®.

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## Maintenance

Provided the tubing pump REGLO *Analog* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

## Replacement parts

For repairs you intend to carry out by yourself (out of the warranty period), you can ask your ISMATEC® distributor for:

- parts lists
- replacement parts for pump drive
- wiring diagrams

Please give information on defect, date of purchase, serial-no., and model.

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## Service and repairs

Like all pumps, the REGLO *Analog* contains components that will wear over a period of time.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

## Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Motor type DC-Motor

Speed range

2 channels	3.2 – 160 rpm
4 channels	2.0 – 100 rpm adjustable in 1 % steps

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## Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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## Technical Specifications

### Operating conditions

Temperature +5 to +40°C,  
Rel. humidity max. 80%  
– not condensing, at normal environmental conditions

### Dimensions/Weight

(DxWxH)  
2 channels 178x100x143 mm  
(7"x4"x5<sup>2</sup>/<sub>3</sub>"")  
4 channels 190x100x143 mm  
(7<sup>1</sup>/<sub>2</sub>"x4"x5<sup>2</sup>/<sub>3</sub>"")

Weight  
2 channels 2.0 kg (5 lb)  
4 channels 2.1 kg (5 lb)

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conformément à:  
EN 61326-1, EN 61010-1

## IDEX Health & Science

Home of Upchurch Scientific &  
Ismatec Products

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Ask for the latest  
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH