
REGLO *Digital*

Schlauch- pumpe

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Schutzgrad
IP 30

2 Kanäle

ISM 831 6 Rollen
ISM 832 8 Rollen
ISM 596 12 Rollen

4 Kanäle

ISM 833 6 Rollen
ISM 834 8 Rollen
ISM 597 12 Rollen

Tubing pump

Microprocessor
controlled

Protection rating
IP 30

2 channels

ISM 831 6 rollers
ISM 832 8 rollers
ISM 596 12 rollers

4 channels

ISM 833 6 rollers
ISM 834 8 rollers
ISM 597 12 rollers

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Classe de protection
IP 30

2 canaux

ISM 831 6 galets
ISM 832 8 galets
ISM 596 12 galets

4 canaux

ISM 833 6 galets
ISM 834 8 galets
ISM 597 12 galets

mit Settings-Menü

with menu settings

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Deutsch

English

Français

ISMATEC®

IDEX
HEALTH & SCIENCE

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Please note (cassettes)

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes. For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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Contents

Dispensing by volume within a pre-set time	24
Intermittent dispensing by volume	25
by time	26
Number of dispensing cycles	27
Drip-free dispensing	28
Pumping against pressure	29
When the pump is idle	29
Overcurrent protector	29
Analog interface RS232 OUT	30
Serial interface RS232 IN	31
RS232 OUT	32
Commands	33
Cascading several pumps	38
Flow rates	39
Programming software	40
Accessories Foot switch	40
Spare cassettes	41
Cassette material	41
Preparing the cassettes with pressure lever	42
Service and repairs	43
Replacement parts	43
Disposal	43
Technical specifications	44

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⚠ Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Caution: rotating pumphead creates a pinch and crush hazard.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping and dispensing applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The pump must not be operated outside the destined operating and environmental conditions. In particular, it is the user's responsibility not to exceed the specified maximum differential pressure for the pump-head.
- The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when pump-heads, cassettes or tubing are inserted or changed.
- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in.

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Safety precautions

- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.

 The pump should not come in contact with the liquid.

- Only new fuses, according to the specifications stated on page 8 in this manual, must be used.
- The fuse holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC® 2 years

All other parts, excluding consumables: 1 year

In case of any queries, please contact your ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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
Contents of the package

- Pump REGLO *Digital* (type as ordered)
- 2 or 4 cassettes
- 1 power cord with country specific mains plug
- Operating manual

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 RS232 IN (female) for operation controlling via PC
- 2 RS232 OUT (male) for cascading up to 8 pumps
Digital input (TTL level) Run/Stop
Digital output 0 – 10 kHz (speed output), Busy (TTL Level)
- 3 Mains socket
- 4 Fuse holder

⚠ Mains voltage

Mains voltage	Fuse rating
100–240 V _{AC} 50–60 Hz	1 x T2A/250V

⚠ Socket/Power cord

Use exclusively the originally applied power cord. The socket must be earthed (protective conductor contact).

Changing the fuses

- ➔ Switch the pump off, pull out the mains plug.
- 1 Pull out the fuse holder by opening it with a small screwdriver (size 0).
 - 2 Remove fuse (a) and replace it with the spare fuse (b) 1 x T2A/250V.

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

Operating panel

(* = Settings page 11)

- 1 **On/Off switch**
- 2 **Digital LED display**
- 3 **Increase value**
- 4 **Reduce value**
* settings
- 5 **RUN/STOP (double function)**
 - a Starts and stops the pump
 - b Interrupts and continues a dispensing cycle (page 20)
➔ The function Start/Stop can also be triggered by an external device such as a foot switch. (Accessories on page 41).
- 6 **RESET (double function)**
 - a **Rotation direction**
Changes the rotation direction. By changing the rotation direction or starting the pump, the rotation direction is indicated in the display with a running hyphen.
 - b **RESET**
 - Resets the accumulated volume in mode TOTAL to zero (page 18)
 - Interrupts dispensing cycle for good (page 20)
- 7 **CAL/MAX (double function)**
 - a **CAL** = Calibrating key for the flow rate or dispensing volume
 - b **MAX** = Max. speed in the modes PUMP•Flow rate and TOTAL (keep the MAX key pressed down while the pump is running)

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

Operating panel

(* = Settings page 11)

- 8 **MODE key**
Selecting operating mode
* default
- 9 **LEDs for active operating mode**
- **PUMP • Flow rate**
Pumping by flow rate (ml/min)
 - **DISP • Volume**
Dispensing by volume (ml)
 - **DISP • Time**
Dispensing by time
 - **PAUSE • Time**
Pause time for dispensing at intervals
 - **TOTAL**
Puming with Read-out of totally delivered volume

Start-up information

After switching on the power supply switch, the following values are displayed:

- 1 LED test »8.8.8.8.«
- 2 Firmware version
- 3 Displays the operating mode last used

⚠ **Before starting the pump for the first time or after connecting a new pump-head,** the ID-code of the pump-head must be entered (see basic settings, page 15).

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Selecting the basic settings

- 1 Switch the pump off and wait for 5 seconds
- 2 Keep the settings key pressed and switch the pump on
- 3 Select the basic settings by using the ▲▼ keys
 - Cycles (dispensing cycles)
 - Back-steps (number of back-steps for volume dispensing)
 - Address (RS232, pump address)
 - Foot switch (Toggle/Direct)
 - Autostart (Standby/RUN)
 - Digits (Number of digits displayed for flow rate and dispensing volume)
 - Display Speed (Type of speed display)
 - Tube (Inner diameter)
 - Head = ID-Code (pump-head)
 - Init (resets basic settings to the default values except pump-head and tube diameter)
- 4 ok = confirm your selection
- 5 Select the value by using the ▲▼ keys, or
 - press escape for cancelling wrong input
 - the values can be reset to the default value by using the default key

- 6 ok = confirm entered value
- 7 escape
Quit menu basic settings

➔ **Basic settings**
refer to the following pages

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⚠ Please note

The settings menu requires a break of **5 seconds** between power off and on.

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Number of dispensing steps
0..9999 (0 = Default value,
infinity)

Backsteps

Number of roller back-steps (for
drip-free dispensing in the mode
DISP Volume)
from 0 .. 100 (0 = Default value)

Address (RS232)

Pump address 1 – 8 (page 32)
(1=default value)

Foot switch

- Toggle (default value)
 - Step on foot switch
= pump starts running
 - Press foot switch again
= pump stops
- Direct
the pump runs as long as the
foot switch is pressed.

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
Basic settings

Configuring Autostart

Autostart determines the pump action when the mains supply is on again, i.e. after power black-out or when the pump is controlled via mains power.

Standby

Pump is stand-by when switched on (default value)

RUN

Pumpe runs when switched on.

Digits (display)

Number of significant digits, 3 or 4 (3 = default value) when displaying flow rate and dispensing volume.

Display Speed

Display of the pump speed in mode Pump•Flow rate as

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Tubing inner diameter

Configuring the tubing inner diameter
(0.13 mm = default)

Head Identification code (ID) of the pump-head, see also table on page 39
(2-06 = default value)

➔ The pump-head ID is composed of:

- number of channels (2 or 4)
- number of rollers (6, 8 or 12)

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Basic settings

Init (Initializing)

Pressing the ok key resets the following basic settings to the default values (except pump-head and tube diameter):

Mode	PUMP•Flow rate
Speed	40 rpm
Dispensing volume	30 roller steps
Dispensing time	4.50 sec
Pause time	2.00 sec
Rotation direction	clockwise
Back-steps	0
Number of dispensing cycles	0
Pump-head digits display	3
Foot switch	Toggle
Autostart	Standby
Display	Flow rate

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Inserting the tubing

⚠ First enter the tube i.d. as described on page 15.

- 1 ➔ Switch the pump off
- 2 Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously
- 3 Insert the 3-stop color-coded tubing with one stopper into the cassette
- 4 Let the tubing hang down (prevent it from being twisted)
- 5 Insert the tubing with the second stopper at the other end of the cassette
- 6 Reinsert the cassette into the roller-head
- 7 ➔ Switch the pump on

➔ When the pump is idle, release all cassettes at the fixing-tongue (2).

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


3-stop color-coded tubing

Please refer to our detailed tubing documentation.

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Pumping by flow rate

- Change mode to PUMP • Flow rate
- Enter the required flow rate with the ▲▼ keys (displayed in ml/min) Start pump with RUN/STOP
 ➔ The flow rate can also be adjusted while the pump is running.
 ➔ For an accurate flow rate we recommend to calibrate the pump (see page 19).
- For fast filling or emptying the system, keep pressing the CAL/MAX key while the pump is running.

➔ During the pumping process the following values can be retrieved with the MODE key:

- Flow rate (PUMP rpm)
- Totally pumped volume (TOTAL)

TOTAL volume

Change with the MODE key to TOTAL. The totally delivered volume is displayed.

Four display readings are available:

Microlitre:	e.g.	17.0 µl = 17.0µ
Millilitre:	e.g.	17.0 ml = 17.0
Litre:	e.g.	17.0 l = 17.0L
Litre:	e.g.	1700 l = 1700.

For zero-setting, either

- press the RESET key once, or
- switch the pump off.

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
Calibrating the flow rate

- Change mode to PUMP • Flow rate
 Enter the required flow rate with the ▲▼ keys
 ➔ Check whether the correct tubing i.d. is selected (see page 15).
- Change mode to DISP • Time
 – Enter 60 seconds by using the ▲▼ keys
 – Start pump with RUN/STOP
 ➔ Pump stops automatically after 60 seconds
 ➔ Measure the dispensed liquid by volume or weight.

Calibrate the ascertained value as follows:

- Return to mode PUMP • Flow rate.
- Press the CAL/MAX-button (displayed value blinks)
 Enter the measured value (balance) with the ▲▼ keys
- Confirm with the CAL/MAX-key (the flow rate setting is readjusted to the initially entered set point as close as possible). The rotation speed has a resolution of $1/100 \text{ rpm}$. After calibration, the initially set flow rate can only be approached.

➔ Setting the default calibration: see page 22

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Please note

A dispensing cycle can be interrupted by pressing the RUN/STOP key

- the volume dispensed up to now is blinking
- by pressing the RUN/STOP key again, the dispensing cycle is continued
- by pressing the RESET-key the dispensing cycle is stopped for good

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➔ For dispensing accurately we recommend you calibrate the pump (see page 21).

- 3 Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➔ The dispensing speed can be adjusted in the mode PUMP • Flow rate.

➔ With the ▲▼ keys the dispensing volume can also be changed **even during the dispensing process**. The newly entered set point appears briefly in the display.

➔ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be retrieved by pressing the MODE key.

➔ The flow rate can be changed even during the dispensing process by using the ▲▼ keys .

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
Calibrating the volume

- 1 Change the MODE to DISP • Volume
 ➔ Use the ▲▼ keys for entering the required dispensing volume
- 2 Start with RUN/STOP
 ➔ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes
 ➔ Measure the dispensed liquid by volume or weight
- 3 Stay in the same mode (DISP • Volume) and press the CAL/MAX-key (display blinks)
- 4 Enter the effectively dispensed volume with the ▲▼ keys
- 5 Save with the CAL/MAX key.
 ➔ The calibrated volume is adjusted to the initially entered set point as close as possible. (The volume per roller step determines the smallest calibration unit. A dispensing volume is always the multiple of this unit. The initially set volume can often only be approached).
- 6 Start pump with RUN/STOP
 ➔ Depending on the application conditions repeated calibration may be necessary.
 ➔ The smallest possible dispensing volume corresponds to the volume of one roller step.

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Default calibration of flow rate

- 1 Change mode to PUMP • Flow rate.
 - 2 Press the CAL/MAX key (displayed value blinks)
 - 3 Press the default key. (Default value blinks).
 - 4 Confirm with the CAL/MAX key.
- ➔ At the same time, the volume is also reset to the default value.

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Default calibration of volume

- 1 Change mode to DISP • Volume
 - 2 Press the CAL/MAX key (displayed value blinks)
 - 3 Press the default key. (Default value blinks)
 - 4 Confirm with the CAL/MAX key.
- ➔ At the same time, the flow rate is also reset to the default value.

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Dispensing by time

The dispensing time can be entered from 0.1s – 999h.

- 1 Change the mode to DISP•Time
- 2 Enter the required dispensing time with the ▲▼ keys
 ➔ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds: 0.1" – 899.9" (in 0.1 s steps)
 Minutes: 15' – 899' (in 1 min steps)
 Hours: 15 h – 999 h (in 1 h steps)

- 3 Start pump with RUN/STOP

➔ With the ▲▼ keys the dispensing time can be changed while the dispensing cycle is in process.

➔ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be retrieved by pressing the MODE-key. The flow rate can be changed during the dispensing process by using the ▲▼ keys.

➔ Very short dispensing times can result in non-repeatable dispensing volumes.

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LED's blinken wechselsei-
tig

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1111 = Volumen zu klein
Zeit zu lang

Dispensing by volume within a pre-set time

1 Change to mode
DISP • Volume. Enter the
required volume (in µl or ml)
with the ▲▼ keys.

2 Change to the mode
DISP•Time.
Enter the required time with the
▲▼ keys (0.1 s – 999 h)

3 Push the MODE-key once
➔ Both LEDs light up simulta-
neously
➔ 1111 and 9999 indicate an er-
ror message (see below).

4 Start pump with RUN/STOP
➔ The pump dispenses the
required volume within the
pre-set time.

a Changes / Corrections

For changing the dispensing
time or dispensing volume
press the RESET-key when the
pump has stopped.

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Intermittent dispensing (by volume)

Intermittent dispensing by volume with a pre-set pause time

- Go to the mode PAUSE • Time. Enter the pause time (between 0.1 s – 999 h) with ▲▼ keys.
 - Press MODE key until both LEDs for PAUSE and DISP • Volume light up simultaneously.
 - ➔ Enter the dispensing volume with the ▲▼ keys.
 - Start with the RUN/STOP key.
 - ➔ The pump dispenses now at intervals.
 - ➔ Active mode blinks.
 - Pushing the RUN/STOP key again interrupts the dispensing process or pause period.
 - ➔ The display blinks.
 The user can now proceed as follows:
 - ➔ By pressing the RUN/STOP key again, the dispensing or pause period is continued.
 - ➔ By pressing the RESET-key the pump is stopped for good.
- ➔ With the ▲▼ keys the dispensing volume or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.
- ➔ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys . Furthermore, the total volume can be displayed this way.

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
Please note

For entering the number of dispensing cycles see page 13.

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Please note

For entering the number of dosing cycles see page 13.

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Intermittent dispensing (by time)

Intermittent dispensing by time with a pre-set pause time

- Go to mode PAUSE • Time Enter the pause time (between 0.1 s – 999 h) with ▲▼ keys.
- Press MODE key until both LEDs for PAUSE and DISP•Time light up simultaneously. ➔ Enter the dispensing time (0.1 s – 999 h) with ▲▼ keys.
- Start with the RUN/STOP key. ➔ The pump dispenses now at intervals. ➔ Active mode blinks.
- Pushing the RUN/STOP key again interrupts the dispensing process or pause period. ➔ The display blinks.

The user can now proceed as follows:
 ➔ By pressing the RUN/STOP key again, the dispensing or pause period is continued.
 ➔ By pressing the RESET-key the pump is stopped for good.

➔ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.
 ➔ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys.

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

Please note

With RUN/STOP a dispensing cycle can be interrupted and resumed.

With RESET an interrupted dispensing cycle can be definitely stopped.

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Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).

See basic settings on page 13.

➔ During the pauses the display shows the remaining number of dispensing cycles.

➔ During long pauses, the remaining interval time can be checked with the CAL/MAX key and adjusted by using the ▲▼ keys.

➔ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key. When resuming the process, the pump continues with the remaining number of dispensing cycles.

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Drip-free dispensing

With programmable roller back-steps (1 – 100 steps), see basic settings bStP (backsteps) page 13

Roller back-steps are active in the following modes:

- a Dispensing by volume (see page 20)
- b Dispensing a volume within a pre-set time (see page 24)
- c Intermittent dispensing by volume (see page 25)

⚠ Reject the first dispensing step. Afterwards, the roller back-steps are compensated.

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Pumping against pressure

The REGLO *Digital* can be used for continuous duty at a differential pressure of max. 1.0 bar (smaller tube sizes and/or lever type cassettes allow higher pressures).

When the pump is idle

When the pump is idle, we recommend you release the tubing from pressure. Releasing the cassette on the right side is sufficient. This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

Overcurrent protector

The drive REGLO *Digital* features an overload protector. It disables the pump drive in case of a motor blockage. When an overload condition occurs, it is indicated in the display by the letters ‚OL‘ and the pump is stopped.

Whenever this situation occurs, **the pump must be switched off immediately**. Wait for 5 sec. before switching on the pump again. Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure etc.).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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
Analog interface

Pin 4: +5VDC

For connecting an external control unit (current <100mA).

Pin 5: GND

Point of reference for all signals described.

Pin 6: start

Connection for an external trigger device (foot-switch, etc.) For operating modes, refer to Basic settings "Foot switch" page 13

Pin 8: busy

Informs about pump status:

Pumpe in operation

= output +5V

Pumpe in stand-by

= output on ground

(The output resistance is 4.7 kΩ)

Pin 9: speed OUT

A square-wave signal with a frequency proportional to the speed (rpm) is at your disposal:

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⚠ Please note


The analog interface is integrated in the port 'RS232 OUT' and can be used for remote control (e.g. by foot-switch).

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⚠ Please note

When switching on the REGLO *Digital* pump, the serial interface will only respond after a delay of 5 seconds.

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Serial interface

RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

Setting the pump address

Each command string begins with the pump address (1 – 8).

It is factory-set at 1 and can be changed with @.

The address will be stored permanently (even after the pump has been switched off).

The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed by** a carriage return (ASCII 13). The pump confirms most of the commands with an asterisk *. Yes/No inquiries are answered by + (yes) or – (no).

Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

Incorrect command strings are answered by #.

Numerical values are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➔ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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
Each pump requires an individual address.
(see page 31)

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Please note

The values indicated in this table are only approximate and determined as follows: ml/min, per channel, with water and Tygon® tubing, without differential pressure

Order No. = Tygon® ST (R3603)

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you use tubing from ISMATEC®.

Please ask for the detailed tubing documentation.

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Programming Software

LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.
<http://www.ni.com/>

LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our web-site:
<http://www.ismatec.com>

Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments.

Order No. IM0015C168

Adaptor 25 F / 9M

25 pins, female / 9 pins, male for IM0015C168 (1 adapter for each controlled instrument)
Order No. XC0009

Extension cable 9F / 9M

9 pins, female / 9 pins, male for XC0009 (2 m long)
Order No. AG0013

Accessories

Foot switch

Order No. ISM 894

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling bottles etc.

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


MS/CA Anpresshebel

MS/CA pressure lever
MS/CA levier de pression

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Spare cassettes

Cassette	MS/CA Click'n'go
Material	POM-C
Order No.	IS 3510
Cassette	MS/CA pressure lever
Material	POM-C
Order No.	IS 0649
Cassette	MS/CA pressure lever
Material	PVDF
Order No.	IS 3629

➔ The cassettes with pressure lever are available on request. This type of cassette may provides better results both for the Tygon MH tubing and at elevated differential pressure conditions.

➔ The cassettes allow the user to insert tubes with different diameters and materials on the same roller-head.

Cassette material

Materials POM-C

Polyoxymethylene-Copolymer
– good chemical resistance to many organic solvents and strong alkaline chemicals
– is affected by strong acids and oxidizing substances
– UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

Materials PVDF

Polyvinylidene fluoride
– very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
– not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
– stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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Preparing the cassettes with pressure lever

- 1 ➔ Switch the pump off
- 2 Push the pressure lever down (it must be relaxed)
Insert the tubing with 2 colour-coded collars
- 3 Reinsert the cassette(s) onto the rollers. Ensure that it snaps in on both sides.
- 4 ➔ Switch the pump on and start it
- 5 Lift the pressure lever step by step until the medium is drawn in and a continuous flow is achieved (no further step than necessary!)

Note

- Under certain circumstances (e.g. new tubing, elevated differential pressure, etc.) we recommend the user to lift the lever by just one more step
- Never over-tighten the pressure lever
- The setting of the pressure lever determines the occlusion
- The optimum occlusion depends on tubing material, inner diameter, wall thickness, age of tubing, and application conditions
- As a side-effect the occlusion setting also affects the flow rate. This, however, should not encourage the user to use the pressure lever for flow rate adjustments because long-term impact is not predictable
- If, however, the occlusion setting is used to adjust the flow rate, the channel-to-channel flow rate conformity may lead to good results in the beginning of the operation, yet the individual flow rates may increasingly diverge at a later stage
- For an optimum long-term flow rate stability and channel-to-channel conformity we suggest using the Click'n'go cassettes the flow rate

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Service and repairs

Provided the REGLO *Digital* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

Replacement parts

For repairs you intend to carry out by yourself (after the warranty period), you can ask from your ISMATEC® distributor:

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial No., and model.

Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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2 channels 1.6–160 rpm

4 channels 1.0–100 rpm

digitally adjustable in steps of 0.01 rpm

Differential pressure

max. 1.0 bar (14.5 psi)

Remote control

via RS232 and analog interface (see pages 30–38)

Mains connection

100 – 240 V_{AC} / 50 – 60 Hz

Fuse rating power supply

1 x T2A/250V

Power consumption

Max. 75 W

Protection rating

IP 30

Operating conditions

Temperature +41 to +104°F

Rel. humidity max. 80%

– not condensing, at normal laboratory conditions

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conformément à:
EN 61326-1, EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products
619 Oak Street
Oak Harbor, WA 98277
United States
Tel: +866 339 4653
Fax: +360 679 3830
fluidics.americas@idexcorp.com

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Ask for the latest
ISMATEC® catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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GmbH