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Operating Manual

Mode d'emploi

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REGLO *Quick*

1 Kanal, 4 Rollen
ISM 897

Deutsch

Tubing pump

REGLO *Quick*

1 channel, 4 rollers
ISM 897

English

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REGLO *Quick*

1 canal, 4 galets
ISM 897

Français

ISMATEC®

IDEX
HEALTH & SCIENCE

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Contents

Safety precautions	6
Warranty terms	7
Product	8
Rear panel	9
Mains voltage	9
Voltage setting	10
Changing the fuses	10
Inserting the tubing	12
Starting the pump	13
Overload protector	14
Pumping against pressure	15
When the pump is not in use	15
Flow rates / tubing	16
Running-in period for tubing	17
Tubing life	17

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Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- The pump must not be used:
 - for medical applications on human beings,
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- The pump must be switched off when tubing is inserted or changed.

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Caution: rotating pumphead creates a pinch and crush hazard

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Safety precautions

- The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Pump-heads consist of revolving parts. Therefore, the pump must not be operated before the protective cover is closed.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.

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Safety precautions

- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used.
- The fuse holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC®

from date of delivery: 2 years

All other parts,

excluding consumables: 1 year

In case of any queries, please contact your ISMATEC® representative.

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Warranty terms

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packing or in a packet of equal quality.

The duration of the warranty is not affected by making a claim for warranty service.

Further claims are excluded.

Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- improper operation by the user, or if the pump is diverted from its proper use
- unauthorized modification or misuse by the user or by a third-party
- improper site preparation and maintenance
- operation outside of the environmental and electrical specifications for the product
- use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

Package contents

- Tubing pump REGLO Quick
1 channel, 4 pump rollers
- Power cord
- Operating instruction

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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
Rear panel

- 1 Analog interface
- 2 Mains socket
- 3 Fuse holder with voltage selector 115/230V

Mains voltage

Mains voltage	Fuse rating
220-240 V _{AC}	2x T500mA/250V
110-120 V _{AC}	2x T500mA/250V

Socket/Power cord

Use exclusively the originally supplied power cord. The socket must be earthed (protective conductor contact).

Before starting up

Check if the voltage setting visible in the window of the fuse holder complies with your local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

- ➔ Switch the pump off,
pull out the mains plug.
- 1 Pull out the fuse holder by opening
it with a small screw-driver (size 0).
 - 2 Take out the voltage selector plate.
Turn it and re-insert it into the fuse
holder so that the required voltage
rating is facing the window of the
fuse holder.
 - 3 Insert 2 new fuses
230 V_{AC}: 2x T500mA/250V
115 V_{AC}: 2x T500mA/250V
⚠ Use always 2 slow-blow fuses
of the same type complying with the
local mains voltage.
 - 4 Shut the fuse holder. The voltage
rating is visible in the window.

⚠ Before starting-up


Check if the voltage setting visible in the window of the fuse holder complies with your local mains voltage.

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
- 6 Schutzblende (h) **schliessen**.

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Inserting the tubing

- 1 ➔ Switch the pump off.
- 2 Open the tube bed by snapping open the protective cover (h).
- 3 Insert the tube (standard tubing) (there are 4 different ways - see figures 1–4).
- 4 Place the tube in the centre of the tube bed between the two prism grooves (i).
- 5 The tube must be well stretched if it is inserted as shown in figure 3.

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Starting the pump

- Check the voltage rating in the window of the fuse holder (on rear panel).
If necessary, change the rating as indicated on page 10–11.
- Connect the pump tubing to the system.
- Connect the pump to the mains.

⚠ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case we recommend you wet the tubing and to run the pump with the tubing inserted for about 15 to 30 minutes.

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
Starting the pump

- 1 Set the required speed on the 2-digit speed selector.
99 = 160 rpm (max. revolution)
➔ The rotation speed can be adjusted while the pump is running.
- 2 Choose the rotation direction.
- 3 Start the pump with the main switch
= run/stop function
- 4 Overload indicator (red LED)

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Overcurrent protector

The tubing pump REGLO *Quick* features an overload protector.

When an overload condition occurs, an activated overcurrent protector is indicated by the red LED (4) and the pump is stopped.

Whenever this situation occurs, the pump must be switched off immediately.

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (e.g. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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Pumping against pressure

The REGLO *Quick* can be used for continuous duty at a differential pressure of max. 1.5 bar (smaller tube sizes allow higher pressures).

In case of any doubts please contact your ISMATEC® agent.

When the pump is not in use

When the pump is idle, we recommend you release the tubing from pressure. This helps to protect the tubing from unnecessary strain and prolongs its service life.

Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

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⚠ Flow rates per channel

In the following tables the indicated values are only approximate and determined as follows:
ml/min, per channel, with water (22°C) and Tygon® tubing, without differential pressure.

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➔ Only the following tubing material must be used:

- Tygon LFL
 - Tygon ST (R3603/R3607)
 - Tygon HC (F-4040-A) yellow
 - Ismaprene
 - Silicone peroxide
- ➔ Other tubing material must not be used

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Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you run in new tubing with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service life of the tubing depends on the application and the tubing material used. Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service life is stated in our detailed tubing documentation (without obligation!). In order to improve the lubrication and service life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, pin 2 must be connected with pin 1 (GND).

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND).

Pin 4, direction

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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- Pin 2, remote
- Pin 3, start
- Pin 4, direction
- Pin 13, internal speed

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Analog interface

Pin 5, speed IN

For external speed control
0–5V, 0–10V or 0–20mA, 4–20mA,
alternatives by means of DIP switches
inside the pump (see p. 21).

Pin 7, +20V_{DC}

About +20 V_{DC} are available (max.
current 0.2 A).

Pin 9, speed OUT

Frequency proportional to the rotation
speed: 0–8 kHz,
3.2–160 rpm

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
Analog interface

Pin 10, +5V_{DC}

About +5 V_{DC} are available.
(max. current 0.1 A)

Pin 13, internal speed

Analog interface activated
(Pin 2 on GND)

- Pin 13 open: The rotation speed is adjusted via pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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Settings of switch S1


→ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!
Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pins		Imp.	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0–5V	18k Ω	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0–10V	38k Ω	OFF	ON	OFF	OFF	OFF	ON
	0–20mA	250 Ω	OFF	OFF	ON	OFF	ON	ON
	4–20mA	250 Ω	OFF	OFF	OFF	ON	ON	OFF

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Foot switch

Order No. ISM 891

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

Maintenance

Provided the tubing pump REGLO *Quick* is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective REGLO *Quick* to your ISMATEC® representative. Please give information on defect, date of purchase, serial No., and model.

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Spare parts

After the warranty period, your ISMA-TEC® representative will be pleased to send you on request:

- spare parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial No., and model.

Disposal

Please retain packing materials until the product warranty ends. Afterwards please discard packing materials in an environment-friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Technical Specifications

Operating conditions

Temperature +41 to +104°F
Rel. humidity max. 80 %
– not condensing, at normal laboratory conditions

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conformément à:
EN 61326-1 / EN 61010-1

IDEX Health & Science

Home of Upchurch Scientific &
Ismatec Products

619 Oak Street
Oak Harbor, WA 98277
United States

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Ask for the
latest ISMATEC®
catalogue.

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ISMATEC® - Your competent partner for demanding metering and dispensing applications.

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Ismatec Laboratoriumstechnik GmbH