ProgEdit

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SOF 104

Software
for
MCP-Process drives

SOF 104

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SOF 104

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Introduction

The *ProgEdit* software enables the user to write individual command sequences (programs) and to download them to the drive.

The drive MCP-Process offers 4 program memories, each one enabling the user to store an individual program.

- ➔ The drive is capable of running these programs independently of the PC (as a stand-alone unit).
- ➔ In the PC the programs can be stored, reopened and printed («File» menu).

Installing *ProgEdit*

1. Start »Windows«
2. Insert the disk into the disk drive
3. From the Windows **Start** menu, choose **Run**
4. Type a:SETUP and click OK
5. The program will guide you through the installation process.

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Please note

For using the RS232- interface under IP-65 conditions please refer to page 51 in the Operating Manual of the MCP-Process.

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Connecting the drive to the PC

Connect a serial interface to the RS232 input on the back of the MCP-Process. Use the enclosed connecting cable AG0013. If your PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

Setting the serial interface

- Start *ProgEdit*.
- Open the menu Option / Settings. The dialog »Settings« appears.
- Click the button »Properties of communication port«.
- Select the interface your drive is connected to and confirm by pressing the OK button twice.

Example: Writing the first program

- In *ProgEdit* choose the command Drive / Run/Stop. In the command window appears Action{Run}.
- Select the command Times / Delay Seconds. DelaySec{Seconds=3} appears in the command window.
- Choose the command Drive / Run/Stop. Click »Stop« in the parameter area. Action{Stop} appears in the command window.
- Save the program by using the File menu.

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Downloading the program to the drive

Make sure the MCP-Process is connected and switched on.

1. Select the required program memory (a).
2. Press the button »Download« (b). Now your program is downloaded to the selected program memory of the drive.
 - ➔ A previously stored program is automatically overwritten.

Starting the program

If required the MCP-Process can now be disconnected from the PC.

- ➔ Push the RUN/STOP key on the drive. The pump now starts to carry out the program.
- ➔ Example of page 6: The drive starts at the pre-set speed and runs whilst the programmed time (3 seconds) is counted down in the display. Then the drive stops again. The program is now terminated and »Pro1« appears in the display.

Stopping, resuming and resetting

A running program can be interrupted any time by using the RUN/STOP key. An interrupted program can be resumed by pushing the RUN/STOP key again.

Push the Reset button if you want **to terminate** an interrupted program.

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
**Application example:
Pressure control**

Pumping against a filter
using a pressure control

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Pressure control unit

- 1 Warning light with sensor cable
- 2 Fastening rod
- 3 Double socket
- 4 Guiding rod
- 5 Adapter for pressure detector
- 6 Reducing fitting
- 7 Tube nipples (3 sizes) incl. tube clamp
- 8 Pressure sensor with cable

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⚠ Please note

Before mounting the pressure control unit, **wrap PTFE tape** around all threads of the pressure sensor, tube nipples and reducing fittings.

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Mounting the pressure control unit

We recommend you to proceed as follows:

1. Mount the fastening rod (2) to the pump-head (instead of the mounted guiding axle with hexagon head cap screw)
2. Screw the guiding rod (4) into the adapter
3. Screw a suitable tube nipple (7) into the left and right thread of the adapter (5) (with reducing fittings if necessary)
4. Mount the pressure sensor (8) onto the adapter (5)
➔ When mounted, the seal face of the pressure sensor must sit close on the adapter.
5. Mount the complete pressure control unit by means of the double socket on the fastening rod of the pump-head.
6. Connect the cable of the pressure sensor to the analog interface (on the back of the drive).

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Pressure control

There are two »pressure control« commands:

The actual pressure control functions **according to a set-point** by changing the pump speed. The set speed determines the maximum speed which is not exceeded even if the pressure is lower than the set-point.

With the second version, the drive **delivers at the set speed** as long as the measured pressure is smaller than the set pressure. As soon as the measured pressure exceeds the threshold value the subsequent command is carried out.

Pumping against a filter

A typical application for using pressure control is pumping against a filter.

The drive works at a constant speed as long as the filter is new.

As soon as the filter reaches a certain degree of clogging a warning light starts to blink and the speed is automatically reduced. This prevents the pressure from exceeding a certain level even if the filter cannot be changed immediately.

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Pumping against a filter

In principle this application can be programmed with four commands (see program example pressure control, page 8):

1. Define detector (I/O's / Define sensor)
2. The command with the pressure as upper threshold value (I/O's / Control)
3. Setting the output the warning light is connected to (Outputs / Output 1 value).
4. The command pressure control (I/O's / Control)

MODE key function for pressure commands

During both ways of pressure control the set pressure is indicated in the display. During operation it can be adjusted with the arrow keys. The user can retrieve more information by using the MODE key. Next to the »PROGRAM« LED a second LED lights up.

The following information appears in turn:

LED	Display
PUMP rpm	Speed
PUMP flow rate	Flow rate
TOTAL	Measured pressure

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LED	Affichage
PUMP rpm	Speed ¹⁾
PUMP flow rate	Flow rate ²⁾
TOTAL	Measured ³⁾ pressure

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Using the Online Help

In the Online Help all commands are explained in details and illustrated with examples.

Demonstration programs

Load the enclosed demonstration programs and use them as starting-point for your own programs.

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A Unit of IDEX Corporation
Laboratoriumstechnik GmbH